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DAY 1

TIME	TOPICS
09:00 – 10:00	INTRODUCTION OF OPTOTEK COMPANY, PRODUCTION AND SERVICE DEPARTMENT
10:00 – 11:00	INTRODUCTION TO USERS AND SERVICE MANUALS; INTRODUCTION OF TOOLS, MEASURING AND TESTING EQUIPMENT
11:00 – 12:00	PACKAGING INTRODUCTION TO THE SYSTEM, OVERVIEW OF BASIC CONFIGURATION, OPTIONS, ACCESSORIES OPTICAL FIBER PREPARATION HANDPIECE PREPARATION OVERVIEW OF USER MENUS AND OPERATION MODES
12:00 – 13:00	lunch
13:00 – 14:00	ERROR CODES AND THEIR TROUBLE SHOOTING OVERVIEW
15:00 – 16:00	CALIBRATION PROCEDURE ANNUAL CHECK AND MAINTENANCE PROCEDURE
16:00 – 17:00	SUB-ASSEMBLIES REPLACEMENT

DAY 2

TIME	TOPICS
9:00 – 10:00	SOFTWARE INSTALLATION PROCEDURE
10:00 – 12:00	FAULTS SIMULATIONS – REPAIR PROCEDURE
12:00 – 13:00	lunch
13:00 – 16:00	SOLVING REAL SERVICE PROBLEMS - EXAM

16:00 – 17:00	SERVICE TRAINING OVERVIEW, DISCUSSION
---------------	---------------------------------------

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2.2 依照原廠以下 Installation and Annual Check Test Report 填寫安裝記錄。

Installation and Annual Check Test Report – LacriMax	
SERIAL NUMBER:	_____
SW version:	_____
<div style="border: 1px solid black; padding: 5px; margin-bottom: 10px;">1. INTRODUCTION</div> <p>SERVICE MANUAL WAS READ AND UNDERSTOOD: YES</p> <p>PERSON RESPONSIBLE FOR INSTALLATION / ANNUAL CHECK (IN CAPITAL LETTERS): _____</p> <p>DATE: _____ SIGNATURE: _____</p> <div style="border: 1px solid black; padding: 5px; margin-bottom: 10px;">2. UNPACKING THE SYSTEM (do not fill in this part in the case of Annual check)</div> <p>TRANSPORT BOX WAS NOT DAMAGED: <input type="checkbox"/>YES <input type="checkbox"/>NO (IF NO, PHOTO OF THE TRANSPORT BOX DAMAGE MUST BE ATTACHED)</p> <p>SHOCK WATCH INDICATOR WAS NOT RED: <input type="checkbox"/>YES <input type="checkbox"/>NO (IF NO, PHOTO OF THE SHOCK WATCH INDICATOR MUST BE ATTACHED):</p> <p>TRANSPORT BOX WAS STORED FOR THE FUTURE TRANSPORTATION: <input type="checkbox"/>YES</p> <p>PERSON RESPONSIBLE FOR INSTALLATION (IN CAPITAL LETTERS): _____</p> <p>DATE: _____ SIGNATURE: _____</p> <p>All procedures mentioned afterwards refer to the Service Manual if not defined otherwise.</p> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;">3. ELECTRICAL TEST</div>	

No.	ITEM	TEST	CONFORMITY
1	EMERGENCY SWITCH	Press Emergency switch. Unit must terminate	
2	INTERLOCK SWITCH	Disconnect Interlock switch. "Door Interlock" warning is displayed	
3	FIBER INTERLOCK	Disconnect Fiber. "Fiber Interlock" warning is displayed	
4	FIBER TYPE DETECTOR	Reconnect Fiber. "Fiber 400/600µm detected" warning is displayed	
5	FOOTSWITCH	Disconnect footswitch or hold footswitch in standby mode. "Footswitch Interlock" warning is displayed	
6	COOLING FAN	Run unit full power for 2 minutes. Rotation of cooling fan must be increased	

PASS: ☐YES ☐NO

4. LASER POWER TEST (use new, optical fiber defined in User's Manual)

Laser output measurement (Default Repeat pulse mode ON 90ms / OFF 50ms)

AVERAGE POWER	Measured laser output power	Required $\pm 10\%$	CONFORMITY
6.42W @ 400um FIBER		5.78 W – 7.06 W	
6.42W @ 600um FIBER		5.78 W – 7.06 W	

PASS: ☐YES ☐NO

5. AIMING BEAM QUALITY / POWER (use new, optical fiber defined in User's Manual)

PARAMETER @ 400um FIBER	REQUIREMENT	Meas. Result [µW]	PASSED
QUALITY	Round & Homogeneous	Visual inspection	
MAX. POWER	$P = 300 \pm 100 \mu W$		
PARAMETER @ 600um FIBER			

QUALITY	Round & Homogeneous	Visual Inspection	
MAX. POWER	P = 300 +/- 100 μ W		

PASS: ☐YES ☐NO

OVERALL CONFORMITY CONFIRMATION

POINTS 1 TO 5 FROM INSTALLATION / ANNUAL CHECK REPORT ARE CONFORMED:

☐YES ☐NO

IF NO, CONTACT OPTOTEK DISTRIBUTOR FOR SUPPORT

LASER INSTALLATION RECORD (from USER'S MANUAL page 62) FILLED IN CASE OF THE NEW INSTALLATION: ☐YES ☐NO

LASER MAINTENANCE / SERVICE RECORD TABLE (from USER'S MANUAL page 57) FILLED IN CASE OF THE ANNUAL CHECK: ☐YES ☐NO

PERSON RESPONSIBLE FOR INSTALLATION / ANNUAL CHECK

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Workshop - Why use DCR ?

topics :

- The Method developed by prof. dr. Brigita Drnovšek-Olup has been successfully used since 2006 at the University Medical Centre Ljubljana - Eye Hospital.
- For the effective treatment of watery eyes (up to a high success rate)
- To avoid unwanted charring during the DCR treatment.

- Minimal invasive technique with less bleeding and charring.
- No scarring and fast patient recovery.
- The operation can be performed under local anesthesia on an outpatient basis.

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