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1. Kusumawati, I. and Indrayanto, G. 2013. Chapter 15 Natural antioxidants in cosmetics. Stud. Nat. Prod. Chem. 40: 485-505.
2. Deniz, F.S.S., Orhan, I.E. and Duman, H. 2021. Profiling cosmeceutical effects of various herbal extracts through elastase, collagenase, tyrosinase inhibitory and antioxidant assays. Phytochem. Lett. 45: 171-183.
3. Aburjai, T. and Natsheh, F.M. 2003. Plants used in cosmetics. Phytother. Res. 17: 987-100.
4. Russo, E.B. 2013. Cannabis and cannabinoids: pharmacology, toxicology, and therapeutic potential. Routledge. NY, USA.
5. Newton, D.E. 2013. Marijuana: a reference handbook. ABC-CLIO. Santa Barbara, CA, USA.
6. Roth, B.L., Baner, K., Westkaemper, R., Siebert, D. and *et al.* 2002. Salvinorin A: a potent naturally occurring nonnitrogenous κ opioid selective agonist. Proc. Natl. Acad. Sci. U. S. A. 99: 11934-11939.
7. Chavkin, C., Sud, S., Jin, W., Stewart, J. and *et al.* 2004. Salvinorin A, an active component of the hallucinogenic sage *Salvia divinorum* is a highly efficacious κ -opioid receptor agonist: structural and functional considerations. J. Pharmacol. Exp. Ther. 308: 1197-1203.
8. Wang, Y., Tang, K., Inan, S., Siebert, D. and *et al.* 2005. Comparison of pharmacological activities of three distinct κ ligands (Salvinorin A, TRK-820 and 3FLB) on κ opioid receptors in vitro and their antipruritic and antinociceptive activities in vivo. J. Pharmacol. Exp. Ther. 312: 220-230.
9. 林柏翔。2012。台灣市面上販售之*Salvia divinorum*及台灣原生種鼠尾草中迷幻成分salvinorin A之含量檢測及管理探討。高雄醫學大學藥學研究所碩士論文。
10. Medana, C., Massolino, C., Pazzi, M. and Baiocchi, C. 2006. Determination of salvinorins and divinorins in *Salvia divinorum* leaves by liquid chromatography/multistage mass spectrometry. Rapid Commun. Mass Spectrom. 20: 131-136.
11. Cspers, M.J., Williams, T.D., Lovell, K.M. and *et al.* 2013. LC-MS/MS quantification of salvinorin A from biological fluids. Anal. Methods 5: 7042.
12. Aubin, A.J. and Waite, M. 2008. UPLC-MS/MS analysis of salvinorin

A from *Salvia divinorum*. (Application note 720002571-zh_tw). [https://www.waters.com/content/dam/waters/en/appnotes/2008/720002571/720002571-zh_tw.pdf.]

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Development of the Test Method for Salvinorin A in Cosmetics

YUN-HUA HSU, SHU-CHING CHENG, SHU-HAN CHANG,
SHOU-CHIEH HUANG, SU-HSIANG TSENG AND DER-YUAN WANG

Division of Research and Analysis, TFDA

ABSTRACT

Salvinorin A is the main component of the powerful hallucinogenic plant *Salvia divinorum*. The hallucinatory effect of salvinorin A in *Salvia divinorum* is similar to marijuana, and its drug abuse has become a new trend and needs more concern. According to the Cosmetic Hygiene and Safety Act and relevant regulations, the ingredients must be safe for consumers under customary conditions of use, companies and individuals who market cosmetics have responsibility to ensure the safety of their products. The aim of this study was to develop a method for the determination of salvinorin A in cosmetics. Nine products were used for evaluating robustness of the method. The liquid-liquid extraction was used for sample preparation. Gradient elution of the mobile phase was carried out through an Accucore RP-MS column with total run time of 16 minutes. The recovery tests were carried out by spiking 20 µg/kg and 100 µg/kg standards into blank samples (lotion, soap and base oil). The recoveries of salvinorin A were between 83.1-103.1% and 87.4-108.8%, respectively, and their intraday and interday coefficients of variation were 3.6 and 7.1% for lotion, 5.4 and 6.7% for base oil, and 3.0 and 6.5% for soap at the spiked level of 20 µg/kg. Both accuracy and precision met the guideline for validation of chemical methods. The limit of quantitation was 20 µg/kg for salvinorin A.

Key words: cosmetics, *Salvia divinorum*, salvinorin A