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1. European Commission Directorate General for Health and Food Safety. 2019. Analytical quality control and method validation procedures for pesticide residues analysis in food and feed. SANTE/12682/2019. [<https://www.eurl->

pesticides.eu/userfiles/file/EurlALL/AqcGuidance_SANTE_2019_1268_2.pdf

2. U.S. Food and Drug Administration. 2019. Guidelines for the validation of chemical methods in food, feed, cosmetics, and veterinary products. [<https://www.fda.gov/media/81810/download>]
3. U.S. Food and Drug Administration. 2020. Inductively coupled plasma-mass spectrometric determination of arsenic, cadmium, chromium, lead, mercury, and other elements in food using microwave assisted digestion. [<https://www.fda.gov/media/87509/download>]
4. Codex Alimentarius. 2018. Codex Alimentarius - Joint FAO/WHO Food Standards Programme. [<http://www.fao.org/documents/card/en/c/I8608EN/>]
5. National Association of Testing Authorities. 2018. General guidelines for the validation and verification of chemical test methods. [<https://www.nata.com.au/phocadownload/gen-accreditation-guidance/Validation-and-Verification-of-Quantitative-and-Qualitative-Test-Methods.pdf>]
6. Royal Society of Chemistry. 2008. The estimation and use of recovery factors. [https://www.rsc.org/images/recovery-factors-technical-brief-21A_tcm18-214880.pdf]
7. 衛生福利部食品藥物管理署。2013。食品化學檢驗方法之確效規範。102 年 9 月 9 日公布。
8. European Union Reference Laboratories. 2016. Analysis of carbofuran (sum) by applying hydrolysis on QuEChERS extracts. EURL-SRM-Analytical Method Report. [http://www.eurl-pesticides.eu/userfiles/file/EurlSRM/EPRW2016_Lemke_PD_003_Carbofuran.pdf]

9. European Union Reference Laboratories (EURL). 2017. Analysis of bifenazate (sum) by the QuEChERS method using LC-MS/MS. EURL-SRM-Analytical Method Report. [http://www.eurl-pesticides.eu/userfiles/file/EurlSRM/meth_Bifenazate_EurlSRM.pdf]
10. Sandín-España, P., Sevilla-Moran, B., Lopez-Goti, C., Mateo-Miranda, M. M. and Alonso-Prados, J. L. 2016. Rapid photodegradation of clethodim and sethoxydim herbicides in soil and plant surface model systems. Arab. J. Chem. 9(5): 694-703.
11. Sevilla-Morán, B., López-Goti, C., Alonso-Prados, J. L. and Sandín-España, P. 2013. Degradation of cyclohexanedione oxime herbicides. In Herbicides—Advances in Research. pp. 101-130. IntechOpen, London, UK.
12. Australian Pesticides and Veterinary Medicines Authority. 2008. Dichlorvos environmental assessment. [<https://apvma.gov.au/sites/default/files/publication/15081-dichlorvos-enviro.pdf>]
13. European Union Reference Laboratories. 2013. Analysis of dicofol via QuEChERS -use of isotope labeled dicofol to improve precision. EURL-SRM-Analytical Method Report. [http://www.eurl-pesticides.eu/library/docs/srm/EurlSrm_Observations_dicofol.pdf]
14. European Union Reference Laboratories. 2013. Analysis of organotin-pesticides by the QuEChERS method- impact of acidifying on the recoveries. EURL-SRM-Analytical Method Report. [http://www.eurl-pesticides.eu/library/docs/srm/EURL_observations_Organotins.pdf]
15. Divito, C. B., Davies, S., Masoudi, S. and Muhoro, C. N. 2007. Relative stability of formamidine and carbamate groups in the bifunctional pesticide formetanate hydrochloride. J. Agric. Food Chem. 55(14): 5377-5382.

16. European Union Reference Laboratories. 2015. Analysis of pyridate and its metabolite pyridafol by QuEChERS and LC-MS/MS. EURL-SRM-Analytical Method Report. [http://www.eurl-pesticides.eu/userfiles/file/EurlSRM/EURL_Observations_Pyridate.pdf]
17. CVUA Stuttgart. 2004. QuEChERS-a mini-multiresidue method for the analysis of-pesticide residues in low-fat products. [<https://www.quechers.eu/pdf/reality.pdf>]
18. European Union Reference Laboratories. 2016. Analysis of pymetrozine by the QuEChERS method- impact of pH on recovery rate. EURL-SRM-Analytical Method Report. [http://www.eurl-pesticides.eu/userfiles/file/EurlSRM/EurlSrm_Observations_Pymetrozine.pdf]
19. U.S. Environmental Protection Agency. 1998. Fact sheet for thiodicarb. [https://www3.epa.gov/pesticides/chem_search/reg_actions/reregistration/fs_PC-114501_1-Dec-98.pdf]
20. Rahman, M. M., Abd El-Aty, A. M., Kim, S. W., Lee, Y. J., Na, T. W., Park, J. S., Shin, H. C, and Shim, J. H. (2017). Simultaneous determination and identity confirmation of thiodicarb and its degradation product methomyl in animal-derived foodstuffs using high-performance liquid chromatography with fluorescence detection and tandem mass spectrometry. Journal of Chromatography B, 1040, 97-104.
21. European Committee for Standardization. 2018. Foods of plant origin– Multimethod for the determination of pesticide residues using GC- and LC- based analysis following acetonitrile extraction/partitioning and clean-up by dispersive SPE – Modular QuEChERS-method. NF EN 15662: 2018 (English version).