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102.10.29

製造廠名稱 : STORZ OPHTHALMICS, INC.

21 Park Place Boulevard North,
Clearwater, FL 34619
U.S.A.

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PRODUCT IDENTITY AND INDICATIONS

Storz ultraviolet absorbing posterior and anterior chamber intraocular lenses manufactured by Storz Ophthalmics Inc., are optical implants for the replacement of the human crystalline lens. Storz intraocular lenses are indicated for primary implantation for the visual correction of aphakia in patients 60 years of age or older where a cataractous lens has been removed. Posterior chamber lenses are indicated after cataract removal by extracapsular extraction methods. Anterior chamber lenses are indicated subsequent to either a primary intra- or extracapsular cataract extraction or after prior surgery such that a secondary implant procedure is required. Storz intraocular lenses are made from polymethylmethacrylate (PMMA) with a UV absorbing compound added. Lenses are available in powers (in situ) from +4.0 to +34.0 diopters in 0.5 diopter increments.

CONTRAINdications

Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the risk to the patient's eyesight: congenital bilateral cataracts; patients with recurrent anterior or posterior segment microphthalmos or certain forms of chronic angle-closure glaucoma; for anterior chamber implantation, previous history of, or predisposition to, retinal detachment; patients in whom the intraocular lens may interfere with the ability to observe, diagnose or treat posterior segment disease; surgical difficulties at the time of cataract extraction which might increase the potential for complications, e.g., persistent bleeding, significant iris damage, inability to clean the anterior chamber of vitreous, uncontrollable positive pressure, significant vitreous loss, or significant anterior chamber bleeding; only one eye with potentially good vision; medically uncontrollable glaucoma; corneal endothelial dystrophy; proliferative diabetic retinopathy.

WARNINGS

As with any surgical procedure there is risk involved. Potential complications accompanying cataract or implant surgery may include, but are not limited to the following: lens dislocation, nonpigment precipitates, corneal macular edema, corneal edema, pupillary membranes, flat anterior chamber, iris prolapse, hypopyon, and secondary glaucoma.

The safety of intraocular lens implantation has not been substantiated in patients with pre-existing ocular conditions (chronic drug miosis, glaucoma, amblyopia, diabetic retinopathy, previous corneal transplant, history of retinal detachment or iritis, etc.). Physicians considering lens implantation in such patients should explore the use of unsatisfactory to meet the needs of the patient. Some adverse reactions which have been associated with the implantation of anterior chamber intraocular lenses are as follows: hypopyon, acute corneal decompensation, lens repositioning, lens removal due to corneal touch, lens removal due to inflammation, corneal transplant, intraocular infection, lens replacement, vitreous aspiration for pupillary block, iridectomy for pupillary block.

The long-term effects of intraocular lens implantation have not been determined. Therefore, physicians should continue to monitor implant patients post-operatively on a regular basis.

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88.

8. 11

Storz Ophthalmics' si
mation. Li... to yo...
for in

Physicians considering lens implantation in patients with senile macular or retinal degeneration should be aware that the patient may not achieve any improvement in central visual acuity following intraocular lens implantation.

Patients with ocular pathology, e.g. glaucoma, or corneal diseases, may not achieve the visual acuity of patients without such problems. Secondary glaucoma has been reported occasionally in patients with controlled glaucoma who received posterior chamber lens implants. The intraocular pressure of implant patients with ocular pathology should be carefully monitored postoperatively.

Hypema, macular edema and vitritis have been reported in patients who experience surgical complications associated with the cataract extraction procedure (posterior capsule rupture, detached Descemet's membrane, anterior chamber bleeding, iris damage, and vitreous bulge or loss). Patients who experience operative complications should be carefully monitored postoperatively for the occurrence of these complications.

The safety and effectiveness of a posterior chamber lens, if placed in the anterior chamber, have not been established. Implantation of posterior chamber lenses in the anterior chamber has been shown in some cases to be unsafe.

The need for a secondary iridectomy for pupillary block may be prevented by one or more iridectomies at the time of intraocular lens implantation.

The effectiveness of UV-absorbing intraocular lenses in reducing the incidence of retinal disorders has not been established.

Some Storz intraocular lenses are manufactured with optics that are 5.00mm. Small amounts of lens decentration occurring with an intraocular lens having a narrow or small optic may result in the patient experiencing glare or other visual disturbances under certain lighting conditions. Surgeons should consider this potential before implanting an intraocular lens with a small or narrow optic.

HOW SUPPLIED

Storz intraocular lenses are supplied STERILE in a dry, heat-sealed package. The inner package is terminally sterilized and should be opened only under sterile conditions.

DIRECTIONS FOR USE

To remove the lens, peel open the outer pouch and remove the tray in a sterile environment.

Carefully lift the lid (clear tray) or slide the lid backward (gray tray) to expose the lens.

Remove the lens by grasping the optic, using smooth-edged forceps, and lifting. Locking forceps or needle holders should never be used.

Do not soak or rinse the lens in solutions other than balanced salt solutions or their equivalent.

The lens should be carefully examined to insure that particles have not adhered to it.

Storz one-piece PMMA lenses are designed to compress in the plane of the lens only. To avoid breakage, do not attempt to expand the haptics, flex the haptics out of the plane of the lens or twist or torque the lens.

Prior to implantation, examine the lens for type, power, proper configuration and optical surfaces.

PRECAUTIONS

1. Do not attempt to resterilize this lens.
2. Do not store the lens in direct sunlight or at temperatures greater than 40°C. Keep from freezing.
3. Do not soak or rinse lens in solutions other than balanced salt solution or equivalent.
4. A high level of surgical skill is required for intraocular lens implantation. A surgeon should have observed and/or assisted in numerous surgical implantations and should have completed one or more courses on intraocular lens implantation before attempting to implant intraocular lenses.
5. Do not use if sterile pouch is opened or damaged.
6. The patient must be advised that the doctor or the medical center should be informed of any side effects not referred to in the information the patient received prior to surgery.

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

3. 8. 11
88. 8. 11

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Storz Ophthalmics, Inc.
21 Park Place Boulevard North
Clearwater, Florida 33759
USA

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Storz Ophthalmics' sincere commitment to ophthalmology brings to you the finest quality intraocular for implantation in the human eye.

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Printed on Recycled Paper

BUSINESS REPLY MAIL

FIRST-CLASS MAIL PERMIT NO. 331 CLEARWATER FL
POSTAGE WILL BE PAID BY ADDRESSEE

REGULATORY AFFAIRS
STORZ OPHTHALMICS INC
21 PARK PLACE BLVD N
CLEARWATER FL 33759-9881

NO POSTAGE
NECESSARY
IF MAILED
IN THE
UNITED STATES

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88. 8. 11

Storz Ophthalmics' sincere commitment to ophthalmology brings to you the finest quality intraocular lenses available for implantation in the human eye.

21 PARK PLACE BLVD. N.
CLEARWATER, FL 33759
(813) 724-6600

stotz

STORZ OPHTHALMICS, INC.

STORZ OPHTHALMICS, INC.

According to Federal (USA) regulations, this form must be completed and returned to Storz Ophthalmics, Inc. for all lenses implanted, lost, or destroyed.

Patient Name															
Last															
First															
Patient S.S. No.			-				-				M.I.				
Date of Birth															
Month			Day							Year					
Hospital Name															
Hospital Address															
Surgeon Name															
Last													State		
Date of Surgery															
Month			Day							Year					
OD	<input type="checkbox"/>			OS	<input type="checkbox"/>										
Anterior Chamber	<input type="checkbox"/>			Posterior Chamber	<input type="checkbox"/>										
Primary Implantation	<input type="checkbox"/>			Secondary Implantation	<input type="checkbox"/>										
Principal Investigator															
Last													First		

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
INTRAOKULARLINSE
LENTILLE INTRA-OCTUAIRE
LENTE INTRAOCULARE
LENTE INTRAOCULAR

OPHTHALMICS

Do not store the lens in direct
sunlight or at temperatures
higher than 43°C/110°F

OPHTHALMIC


OPHTHALMIC

EZE-FIT[®]

Intraocular Lens

Equiconvex Optic
One-Piece PMMA
Flexible Haptics

INTRAOKULARLINSE
LENTILLE INTRA-OCTULAIRE
LENTE INTRAOCULARE
LENTE INTRAOCULAR

Do not store the lens in direct
sunlight or at temperatures
higher than 43°C/110°F

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

BLUEVISTA®

Intraocular Lens

Equiconvex Optic
One-Piece PMMA
With Blue Haptics

Do not store the lens in direct
sunlight or at temperatures
higher than 43°C/110°F

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0013645

STORZ Ophthalmics, Inc.
21 Park Place Boulevard North
Clearwater, Florida 33759
USA

Contents: One intraocular lens. Sterile unless enclosed package has
been opened or damaged. See insert for professional information.

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nent GmbH
ergewann 4
lberg, Germany
ackage has
rmation.

ent GmbH
ergewann 4
lberg, Germany
ackage has
rmation.

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BLUE VISTA®

Blue Vista® lenses combine the benefits of STORZ[®] conventional one-piece lens design with the enhanced haptic visibility from traditional three-piece lens models.

Blue Vista:

Feature:

- High Molecular Weight PMMA
- Blue Haptic
- PMMA Haptic
- Equiconvex Optic

Benefit:




- Yag Compatible
- Easier Visibility
- More Durable than Polypropylene
- Added Stability
- Full Aperture Refractive Surface
- Stable Optical A-Constant Plane
- 1:1 Anterior-to-Posterior Curvature Ratio
- Thinnest Possible Optic Design.

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH 4.96mm	UV ABSORBING	OTHER FEATURES
	95BUV	EQUICONVEX	7.0	0	13.75	10°	MODIFIED C	118.0 4.96mm	YES	BLUE HAPTICS
	BV359	EQUICONVEX	5.5	0	12.25	DOUBLE FLAT STEP-VAULTED 2.4°	MODIFIED C	118.0 4.96mm	YES	BLUE HAPTICS
	P325BUV	EQUICONVEX	6.5	0	13.75	10°	MODIFIED C	118.0 4.96mm	YES	BLUE HAPTICS








ASSEMBLED POLYPROPYLENE - PMMA

Three-Piece Polypropylene

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH	UV ABSORBING	OTHER FEATURES
	71UV	BICONVEX	6.0	4	13.75	10°	MODIFIED J POLYPROPYLENE	119.0 5.55mm	YES	
	P013UV	EQUICONVEX	6.5	1	14.00	10°	MODIFIED J POLYPROPYLENE	118.5 5.26mm	YES	
	127UV	EQUICONVEX	7.0	2	13.75	10°	MODIFIED C POLYPROPYLENE	118.0 4.96mm	YES	

ONE-PIECE DESIGNS

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH	UV ABSORBING	OTHER FEATURES
	P359UV	EQUICONVEX	5.5	0	12.25	DOUBLE FLAT STEP-VAULTED 2.4°	MODIFIED C	118.0 4.96mm	YES	MULTI-DIRECTIONAL FLEXIBILITY HAPTIC AVAILABLE IN EXTENDED HIGH DIOPTER POWERS
	P379UV	EQUICONVEX	6.0	0	13.00	STEP-VAULTED 6.5°	MODIFIED C	118.4 5.20mm	YES	MULTI-DIRECTIONAL FLEXIBILITY
	68UV	BICONVEX	6.0	2	13.75	STEP-VAULTED 5.9°	MODIFIED C	118.5 5.26mm	YES	AVAILABLE IN EXTENDED LOW DIOPTER POWERS
	P010UV	EQUICONVEX	6.5	1	13.75	10°	MODIFIED C	118.0 4.96mm	YES	MULTI-DIRECTIONAL FLEXIBILITY
	P366UV	EQUICONVEX	6.5	0	13.4 12.0 EYELET TO EYELET	STEP-VAULTED 6.6°	MODIFIED C	118.5 5.26mm	YES	MULTI-DIRECTIONAL FLEXIBILITY EYELET ON HAPTICS (43MM)
	P525UV	EQUICONVEX	6.5	2	13.75	10°	MODIFIED C	118.0 4.96mm	YES	MULTI-DIRECTIONAL FLEXIBILITY
	95UV	EQUICONVEX	7.0	0	13.75	10°	MODIFIED C	118.0 4.96mm	YES	MULTI-DIRECTIONAL FLEXIBILITY
	P389UV	EQUICONVEX	5.5	0	13.00	DOUBLE FLAT STEP-VAULTED 2.1°	MODIFIED C	118.0 4.96mm	YES	MULTI-DIRECTIONAL FLEXIBILITY

Additional Lens Models on back of page

ONE-PIECE DESIGNS

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	CONSTANT AC DEPTH	UV ABSORBING	OTHER FEATURES
	P351UV	EQUICONVEX	6.5	0	12.50	STEP-VAULTED 7.6°	MODIFIED C	118.5 5.20mm	YES	MULTI-DIRECTIONAL FLEXIBILITY
	P399UV	EQUICONVEX	5.5	0	12.25	DOUBLE FLAT STEP-VAULTED 2.4°	MODIFIED C	118.0 4.96mm	YES	MULTI-DIRECTIONAL FLEXIBILITY NOTCH ON HAPTIC
	P484UV	EQUICONVEX	6.5	1	13.00	10°	MODIFIED C	117.9 4.91mm	YES	MULTI-DIRECTIONAL FLEXIBILITY
	P563UV	EQUICONVEX	6.5	2	13.50	STEP-VAULTED 6.0°	MODIFIED C	118.4 5.20mm	YES	MULTI-DIRECTIONAL FLEXIBILITY
	P329UV	EQUICONVEX	5.0 X 6.0	0	12.00	STEP-VAULTED 7.6°	MODIFIED C	118.4 5.20mm	YES	MULTI-DIRECTIONAL FLEXIBILITY
	P337UV	EQUICONVEX	5.0 X 6.0	0	13.00	STEP-VAULTED 6.5°	MODIFIED C	118.4 5.20mm	YES	MULTI-DIRECTIONAL FLEXIBILITY

Laser Ridge

	68RUUV	BICONVEX	6.5 X 7.0	2	13.75	STEP-VAULTED 6.3°	MODIFIED C	117.3 4.56mm	YES	RIDGED OPTIC EYELET ON HAPTIC WITH CONTROL TIP
--	--------	----------	-----------	---	-------	-------------------	------------	-----------------	-----	--

®

OPHTHALMICS

INTRA- CULAR LENSES

C A T A L O G

FEATURING

THE EZE-FITTM IOL

EXTENDED
RANGE DIOPTER
POWERED INTRACULAR LENSES

&

BIFOCAL
INTRACULAR LENSES

OPHTHALMICS

DIOPTRIC POWER RANGES

Many IOL designs from STORZ Ophthalmics are available from +10.0 to +30.0 D in .5 D steps. Most designs are available from +14.0 to +27.0 D in .5 D steps.

Extended Range Powers available from -18.0 D to +45.0 D on selected models. See page 15 for more information.

"A" CONSTANTS AND AC DEPTHS

The suggested "A" Constants and estimated AC Depths listed in this catalog are presented as guidelines only. It is recommended that personalized "A" Constants and AC Depths be developed based on individual experience with an implant design, surgical technique, measuring equipment and postoperative results.

RETURN POLICY

STORZ Ophthalmics will accept all currently manufactured lenses for exchange only. Discontinued or special purchase lenses may not qualify for exchange.

Qualified lenses in the original unopened package may be returned and exchanged for another lens. Please call 1-800-325-9500 for a return authorization number.

CAUTIONS

All IOLs are limited by Federal (USA) law to sale by or on the order of a physician.

TRADEMARKS/PATENTS

EZE-FIT™ SYSTEM, BLUE VISTA®, CAPSULORBLUE™, STORZ®, and PHACO PROFILE®, are trademarks of STORZ Ophthalmics.

Products are protected under one or more of the following patents: 4,863,465; 4,476,591; 4,481,431; RE 31,640; RE 31,626; 4,504,981; 4,298,994; 4,687,485; 4,813,955, and/or pending patents.

REGULATION POLICY

In order to comply with Medical Device Reporting (MDR) regulations, STORZ Ophthalmics has instituted a system to facilitate reporting of patient injury and malfunctions involving IOLs. All Events including particulate contamination, haptic breakage, and optic scratches due to damage, defect, or handling should be reported to STORZ Customer Service at 1-800-252-7890.

A large, stylized graphic of an intraocular lens (IOL) is shown. It features a central circular optic with a textured surface, surrounded by several thick, curved haptics that fan outwards. The entire assembly is set against a white background with a black rectangular border around the main graphic area.

THE EZE-FIT™ IOL SYSTEM

SERIES 0 ...Most flexible (comparable to assembled PMMA haptics)

SERIES 1 ...Intermediate flexibility

SERIES 2 ...Flexible (comparable to current one-piece PMMA haptics)

Each series maintains comparable flexibility of the haptics throughout the different optic sizes.

GENERATION OF ONE-PIECE PMMA LENSES

Features on all EZE-FIT™ Lenses

- Equiconvex Full Aperture Refractive Surface
- Double flat, step-vaulted haptics
- High molecular weight material










- The EZE-FIT™ IOL System is designed to accommodate any surgical technique.
- The comparable flexibility of haptics on different size optics within a series allows surgeons to change optic sizes without changing insertion technique.
- The model number describes the lens optic size and flexibility of the haptic.
- “A” Constant remains the same throughout the EZE-FIT™ IOL System allowing for no diopter change with different lenses.

OPTIC SIZE	SERIES 0	SERIES 1	SERIES 2
5.0mm	EZE-50 EZE-50S	EZE-150 EZE150A EZE-150S	P434UV P526UV P534UV P408UV
5.5mm	EZE-55 EZE-55N EZE-55S	EZE-155 EZE-155N EZE155A EZE-155S EZJ-155A	P359UV P399UV P499UV P389UV
5.0mm X 6.0mm	EZE-56		P328UV P329UV P336UV P337UV
6.0mm	EZE-60 EZE-60N	EZE-160 EZE-160N EZE160A EZJ-160A	P504UV P508UV P574UV
6.5mm	EZE-65	EZE-165 EZJ-165A	P517UV
7.0mm		EZJ-170A	

THE EZE-FIT™ IOL SYSTEM

SERIES 0 - Most flexible haptic (comparable to assembled PMMA haptics)






One-Piece PMMA Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	EZE-50S	EQUICONVEX	5.0	0	12.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-50	EQUICONVEX	5.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-55S	EQUICONVEX	5.5	0	12.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-55	EQUICONVEX	5.5	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-55N	EQUICONVEX	5.5	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH
	EZE-56	EQUICONVEX	5.0 X 6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-60	EQUICONVEX	6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-60N	EQUICONVEX	6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH
	EZE-65	EQUICONVEX	6.5	0	13.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	VERY FLEXIBLE MULTI-DIRECTIONAL HAPTIC

THE EZE-FIT™ IOL SYSTEM

SERIES 1 - Intermediate haptic flexibility





One-Piece PMMA Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	*A* CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	EZE-150S	EQUICONVEX	5.0	0	12.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-150	EQUICONVEX	5.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE150A	EQUICONVEX	5.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-155S	EQUICONVEX	5.5	0	12.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-155	EQUICONVEX	5.5	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE 155A	EQUICONVEX	5.5	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC

THE EZE-FIT™ IOL SYSTEM

SERIES 1 - Intermediate haptic flexibility continued





One-Piece PMMA Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	EZE-155N	EQUICONVEX	5.5	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH
	EZE-160	EQUICONVEX	6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-160N	EQUICONVEX	6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH
	EZE-160A	EQUICONVEX	6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZE-165	EQUICONVEX	6.5	0	13.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC

THE EZE-FIT™ IOL SYSTEM

SERIES 1 - Intermediate haptic flexibility continued



One-Piece PMMA Modified J Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC/DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	EZJ-155A	EQUICONVEX	5.5	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED J	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZJ-160A	EQUICONVEX	6.0	0	13.00	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED J	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZJ-165A	EQUICONVEX	6.5	0	13.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED J	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	EZJ-170A	EQUICONVEX	7.0	0	13.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED J	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC

THE EZE-FIT™ IOL SYSTEM

SERIES 2 - Flexibility (comparable to current one-piece PMMA haptics)







One-Piece PMMA Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC/DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	P434UV	EQUICONVEX	5.0	0	12.25	DOUBLE FLAT STEP-VAULTED 2.2°	MODIFIED C	118.0 4.96mm 1.22	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	P526UV	EQUICONVEX	5.0	0	12.75	STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC

THE EZE-FIT™ IOL SYSTEM

SERIES 2 - Flexibility (comparable to current one-piece PMMA haptics) continued

One-Piece PMMA Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	P408UV	EQUICONVEX	5.0	0	13.00	STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	P359UV	EQUICONVEX	5.5	0	12.25	DOUBLE FLAT STEP-VAULTED 2.4°	MODIFIED C	118.0 4.96mm 1.22	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC EXTENDED DIOPTER POWERS AVAILABLE PLANO TO +45D
	P399UV	EQUICONVEX	5.5	0	12.25	DOUBLE FLAT STEP-VAULTED 2.4°	MODIFIED C	118.0 4.96mm 1.22	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH
	P499UV	EQUICONVEX	5.5	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH
	P389UV	EQUICONVEX	5.5	0	13.00	DOUBLE FLAT STEP-VAULTED 2.1°	MODIFIED C	118.0 4.96mm 1.22	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC EXTENDED DIOPTER POWERS AVAILABLE AS CUSTOM ORDER PLANO TO +45D
	P512UV	EQUICONVEX	5.5	0	13.25	DOUBLE FLAT STEP-VAULTED 2.5°	MODIFIED J	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC

THE EZE-FIT™ IOL SYSTEM

SERIES 2 - Flexibility (comparable to current one-piece PMMA haptics) continued




One-Piece PMMA Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	P328UV	EQUICONVEX	5.0 X 6.0	0	12.00	STEP-VAULTED 7.6°	MODIFIED C	118.4 5.20mm 1.45	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH
	P329UV	EQUICONVEX	5.0 X 6.0	0	12.00	STEP-VAULTED 7.6°	MODIFIED C	118.4 5.20mm 1.45	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	P336UV	EQUICONVEX	5.0 X 6.0	0	13.00	STEP-VAULTED 6.5°	MODIFIED C	118.4 5.20mm 1.45	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH
	P337UV	EQUICONVEX	5.0 X 6.0	0	13.00	STEP-VAULTED 6.5°	MODIFIED C	118.4 5.20mm 1.45	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	P504UV	EQUICONVEX	6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	P508UV	EQUICONVEX	6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC WITH NOTCH

THE EZE-FIT™ IOL SYSTEM

SERIES 2 - Flexibility (comparable to current one-piece PMMA haptics) continued

One-Piece PMMA Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	P574UV	EQUICONVEX	6.0	0	13.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC HIGH MINUS DIOPTER POWERS AVAILABLE -1 TO -18 WHOLE DIOPTER STEPS
	P517UV	EQUICONVEX	6.5	0	13.25	STEP-VAULTED 3°	MODIFIED C	118.0 4.96mm 1.22	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC HIGH MINUS DIOPTER POWERS AVAILABLE -1 TO -18 WHOLE DIOPTER STEPS CUSTOM ORDER
	P366UV	EQUICONVEX	6.5	0	13.4 12.0 EYELET TO EYELET	STEP-VAULTED 6°	MODIFIED C	118.5 5.26mm 1.51	YES	MULTI-DIRECTIONAL FLEXIBILITY EYELET ON HAPTICS

BLUE VISTA® IOLs

One-Piece Blue PMMA Haptic Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	BVR-150S	EQUICONVEX	5.0	0	12.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	BVR-150M	EQUICONVEX	5.0	0	12.75	DOUBLE FLAT STEP-VAULTED 2.9°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	BVR-155S	EQUICONVEX	5.5	0	12.25	DOUBLE FLAT STEP-VAULTED 3.1°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	BVR-155M	EQUICONVEX	5.5	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	BVR-160S	EQUICONVEX	6.0	0	12.25	DOUBLE FLAT STEP-VAULTED 2.9°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	BVR-160M	EQUICONVEX	6.0	0	12.75	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	BVR-165L	EQUICONVEX	6.5	0	13.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC
	BVR-170L	EQUICONVEX	7.0	0	13.25	DOUBLE FLAT STEP-VAULTED 3°	MODIFIED C	118.1 5.02mm 1.28	YES	FLEXIBLE MULTI-DIRECTIONAL HAPTIC





CAPSULORBLUE™ IOLs

Assembled Blue PMMA Haptic Posterior Chamber Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	P519UV	EQUICONVEX	5.5	0	12.50	3°	MODIFIED J ASSEMBLED PMMA	117.9 4.91mm 1.17	YES	
	P454UV	EQUICONVEX	5.5	0	12.75	4°	MODIFIED C ASSEMBLED PMMA	118.0 4.90mm 1.22	YES	
	P541UV	EQUICONVEX	5.0 X 6.0	0	13.50	7°	MODIFIED J ASSEMBLED PMMA	118.4 5.20mm 1.45	YES	
	P391UV	EQUICONVEX	5.0 X 6.0	0	13.25	3°	MODIFIED C ASSEMBLED PMMA	118.0 4.90mm 1.22	YES	
	560CUV	EQUICONVEX	5.0 X 6.0	0	13.75	10°	MODIFIED C ASSEMBLED PMMA	118.5 5.26mm 1.51	YES	
	P453UV	EQUICONVEX	6.0	0	12.75	4°	MODIFIED C ASSEMBLED PMMA	118.0 4.96mm 1.22	YES	
	P524UV	EQUICONVEX	6.0	0	13.50	3°	MODIFIED J ASSEMBLED PMMA	117.9 4.91mm 1.17	YES	



CAPSULORBLUE™ IOLs continued

Assembled Blue PMMA Haptic Posterior Chamber Lenses


LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	P506UV	EQUICONVEX	6.0	0	13.50	10°	MODIFIED C ASSEMBLED PMMA	118.5 5.26mm 1.51	YES	
	P518UV	EQUICONVEX	6.5	0	13.25	4°	MODIFIED C ASSEMBLED PMMA	118.0 4.96mm 1.22	YES	
	P538UV	EQUICONVEX	6.5	0	13.50	3°	MODIFIED J ASSEMBLED PMMA	117.9 4.91mm 1.17	YES	
	P507UV	EQUICONVEX	7.0	0	13.50	10°	MODIFIED C ASSEMBLED PMMA	118.0 4.96mm 1.22	YES	

ANTERIOR CHAMBER IOLs

One-Piece PMMA AC Lenses

LENS DESIGN	MODEL NUMBER	OPTIC DESIGN	OPTIC SIZE (mm)	POSITIONING HOLES 0.34mm STANDARD	LENGTH (mm)	EFFECTIVE HAPTIC ANGLE	HAPTIC DESIGN	"A" CONSTANT AC DEPTH "S" FACTOR	UV ABSORBING	OTHER FEATURES
	S122UV	EQUICONVEX	6.0	0	12.50	STEP-VAULTED 4.4°	MODIFIED S	115.8 3.68mm -0.022	YES	4PT FIXATION
	L122UV	EQUICONVEX	6.0	0	13.75	STEP-VAULTED 3.7°	MODIFIED S	115.8 3.68mm -0.022	YES	4PT FIXATION

STORZ Equiconvex IOLs – Shaped for Success

STORZ equiconvex optics vs. "standard" biconvex optic. Both 1:1 and 3:1 ratios provide good optical performance; however, only a 1:1 ratio permits a perfectly balanced, thin lens design.

1:1 Anterior-to-Posterior Curvature Ratio

- Outstanding optical performance—less spherical aberration than conventional IOLs.¹
- Minimal image degradation due to decentration or tilt; minimal internal reflection.¹

Stable Optical Planes

- Principal plane designed for minimal detectable shift—regardless of power.


Equiconvex Profile Permits Thinnest Possible Designs

- Thin edges designed to minimize edge glare.
- Low lens mass facilitates implantation.

Maintains Posterior Capsule Contact

- Enhances YAG compatibility, fewer shots required to open posterior capsule.²

STORZ Step-Vaulted IOLs – Secure Placement by Design

- (1) Planar (0°) angulation does not position optic away from iris or anterior capsule;
- (2) 10° angulation holds lens against posterior capsule, but may sacrifice stability;
- (3) Step-vaulted haptics are designed to achieve maximum stability and facilitate implantation, but increase pressure on posterior capsule;
- (4) Double flat, step-vaulted haptics offer a unique solution.

Step-Vaulted Haptic Balances Flexibility with Stability

- High quality one-piece PMMA lens permits easy insertion into the bag.

Facilitates In-The-Bag Implantation

- Flat end of haptic prevents skating over the viscoelastic.²

Stress is Distributed Evenly Along Flat Portion of Haptic

- Step vaulting is designed to assure stability in the bag.

Step-Vaulted Design Conforms to the Anatomy of the Eye

- Haptic is initially angled upward, then flattened.
- Controlled flexibility.

Designed for Controlled Vault to Minimize Tilt and Decentration

- Flattened portion of haptic rests in fornix of capsule to stabilize lens.

Holds IOL Gently Against Posterior Capsule

- Stretches posterior capsule smoothly across optic surface.

Double Flat, Step Vaulted Design Offers These Features:


- Ease of insertion of a planar design.
- The benefits of an angled design.
- The stability of a step-vaulted design.

References:

1. Holladay JT. Evaluating the intraocular lens optic. *Surv Ophthalmol*. 1986;30:385-390.
2. Data on file, Storz Ophthalmics, St. Louis, Mo.

Extended Diopter Range Lenses From STORZ


P359UV

Diopter range available
Plano to +45D. half
diopter steps

Optic Size: 5.5mm
Overall Length: 12.25mm
A-Constant: 118.0

P574UV


Diopter range available
-1D. to -18D. in whole
diopter steps

Optic Size: 6.0mm
Overall Length: 13.25mm
A-Constant: 118.1

Extended Diopter Range Lenses are not available for consignment. Storz will stock the above IOL models in the extended range for immediate delivery. It is recommended that IOLs are ordered before scheduling surgery.

The Following Two IOL Models can be Custom Ordered in Extended Diopter Range


P389UV

Diopter range available
Plano to +45D. half
diopter steps

Optic Size: 5.5mm
Overall Length: 13.00mm
A-Constant: 118.0

P517UV

Diopter range available
-1D. to -18D. in whole
diopter steps

Optic Size: 6.5mm
Overall Length: 13.25mm
A-Constant: 118.0

The above IOL models can be custom manufactured with an expected delivery date of 45 days from the date of order.

*For more information on extended diopter range IOLs from Storz contact your
Storz Sales Representative or call Customer Service at 1-800-325-9500.*

storz

OPTICAL MICROS

10000 N.W. 2nd Avenue • Suite 100 • Miami, FL 33169 • Telephone: 305/667-0000 • Telex: 855-9560

68UV storz +04.0
Posterior Chamber IOL
SN: x00101 Length:13.75

Optic Size: 6.00mm
 Optic Design: Biconvex
 Haptic Design:
 Short C
 Placement: Ciliary sulcus or capsular bag
 Material:
 PMMA with UV absorber

storz Model 68UV
 Length 13.75 mm Power +04.0
 SN: x00101

**MODEL
68UV
diopter
+04.0**

1-x mm-l LENGTH
13.75mm
 OPTIC
6.00mm

EXP. DATE
2003-04
 SN
x00101

REF # **69040** LOT **OTEST**

**SN
x00101**

"A" Factor : P.C. 118.5
 Optic Design : Biconvex
 Haptic Design : Short C
 Material : PMMA with UV absorber
 Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.

REF # 69040

[01] 10757770133212 [21] x00101

50168 Rev A

This certifies that _____
 received the following intraocular lens implant in the
 Left eye
 Right eye
**68UV Posterior Chamber IOL
Length: 13.75 mm
SN: x00101**

Haptic angulation: Step-vaulted
 Optic size: 6.00mm
 Optic design: Biconvex
 Haptic design: Short C
 Placement: Ciliary sulcus or capsular bag
 Material: PMMA with UV absorber

Physician: _____

127UV storz +04.0
Posterior Chamber IOL
SN: x00201 Length:13.75

Optic Size: 7.00mm
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Ciliary sulcus or capsular bag
 Material:
 PMMA optic with UV absorber, polypropylene haptics with blue dye

storz Model 127UV
 Length 13.75 mm Power +04.0
 SN: x00201

**MODEL
127UV
diopter
+04.0**

1-x mm-l LENGTH
13.75mm
 OPTIC
7.00mm

EXP. DATE
2003-04
 SN
x00201

REF # **48040** LOT **OTEST**

**SN
x00201**

"A" Factor : P.C. 118.0
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA optic with UV absorber, polypropylene haptics with blue dye
 Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.

REF # 48040

[01] 10757770111104 [21] x00201

50159 Rev A

This certifies that _____
 received the following intraocular lens implant in the
 Left eye
 Right eye
**127UV Posterior Chamber IOL
Length: 13.75 mm
SN: x00201**

Haptic angulation: 10 degrees
 Optic size: 7.00mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Ciliary sulcus or capsular bag
 Material: PMMA optic with UV absorber

Physician: _____


95UV storz +04.0
Posterior Chamber IOL
SN: x00301 Length:13.75

Optic Size: 7.00mm Haptic angulation: 10 degrees
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Ciliary sulcus or capsular bag
 Material:
 PMMA with UV absorber

storz Model 95UV
Length 13.75 mm Power +04.0
SN: x00301

**MODEL
95UV
DIOPTER
+04.0**

LENTH
13.75mm
OPTIC
7.00mm

EXP. DATE
2003-04
SN
x00301

REF # **59040** LOT **0TEST**

SN
x00301

"A" Factor : P.C. 118.0
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA with UV absorber

Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.

REF # 59040

[01] 10757770121592 [21] x00301

50177 Rev A

This certifies that
 received the following intraocular lens implant in the
 Left eye Right eye Left eye
95UV Posterior Chamber IOL
Diopter: +04.0
 Length: 13.75 mm
 SN: x00301

Haptic angulation: 10 degrees

Optic size: 7.00mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Ciliary sulcus or capsular bag
 Material: PMMA with UV absorber

Physician:


P010UV storz +04.0
Posterior Chamber IOL
SN: x00401 Length:13.75

Optic Size: 6.50mm Haptic angulation: 10 degrees
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Ciliary sulcus or capsular bag
 Material:
 PMMA with UV absorber


storz Model P010UV
Length 13.75 mm Power +04.0
SN: x00401

**MODEL
P010UV
DIOPTER
+04.0**

LENTH
13.75mm
OPTIC
6.50mm

EXP. DATE
2003-04
SN
x00401

REF # **1C040** LOT **0TEST**

SN
x00401

"A" Factor : P.C. 118.0
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA with UV absorber
 Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.

REF # 1C040

This certifies that
 received the following intraocular lens implant in the
 Left eye Right eye Left eye
P010UV Posterior Chamber IOL
Diopter: +04.0
 Length: 13.75 mm
 SN: x00401

Haptic angulation: 10 degrees

Optic size: 6.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Ciliary sulcus or capsular bag
 Material: PMMA with UV absorber

Physician:

[01] 10757770083883 [21] x00401

50050 Rev A

P359UV storz +04.0
Posterior Chamber IOL
Control#x00101 Length:12.25

Optic Size: 5.50mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Capsular bag
 Material:
 PMMA with UV absorber


storz Model P359UV
 Length 12.25 mm Power +04.0
 Control # x00101

MODEL
P359UV
+04.0
 DIOPTER

LENGTH
12.25mm
 OPTIC
5.50mm

EXP. DATE
2000-01
 CONTROL #
x00101

Part # **2R040** Batch # **OTEST**

Control #
x00101

"A" Factor : **P.C. 118.0**
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA with UV absorber

This certifies that _____ received the following intraocular lens implant in the _____
 [] Left eye [] Right eye
P359UV Posterior Chamber IOL
Length: 12.25 mm
Control # x00101

Haptic angulation: Step-vaulted
 Optic size: 5.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Capsular bag
 Material: PMMA with UV absorber

Physician: _____

PART # **2R040**

P329UV storz +04.0
Posterior Chamber IOL
Control#x00201 Length:12.00

Optic Size: 5.00 x 6.00mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Capsular bag
 Material:
 PMMA with UV absorber

storz Model P329UV
 Length 12.00 mm Power +04.0
 Control # x00201

MODEL
P329UV
 DIOPTER
+04.0

l-x mm-l LENGTH
12.00mm
 OPTIC
5.00 x 6.00mm

EXP. DATE
2003-04
 SN
x00601

REF # **2C040** LOT **OTEST**

SN
x00601

"A" Factor : **P.C. 118.4**
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA with UV absorber

**Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.**

REF # **2C040**

This certifies that _____ received the following intraocular lens implant in the _____
 [] Left eye [] Right eye
P329UV Posterior Chamber IOL
Length: 12.00 mm
SN: x00601

Haptic angulation: Step-vaulted
 Optic size: 5.00 x 6.00mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Capsular bag
 Material: PMMA with UV absorber

Physician: _____

[01] 10757770090300 [21] x00601

50063 Rev A

storz +04.0
S122UV Anterior Chamber IOL
SN: x00101 Length:12.50
 Optic Size: 6.00mm
 Optic Design: Equiconvex
 Haptic Design:
 Modified S
 Placement: Anterior Chamber
 Material:
 PMMA with UV absorber

storz Model S122UV
 Length 12.50 mm Power +04.0
 SN: x00101

This certifies that _____ received the following intraocular lens implant in the _____
 Left eye
 Right eye

S122UV Anterior Chamber IOL
Diopter: +04.0

Optic size: 6.00mm
 Optic design: Equiconvex
 Haptic design: Modified S
 Placement: Anterior Chamber
 Material: PMMA with UV absorber

Physician: _____

storz +04.0
BV359 Posterior Chamber IOL
SN: x00401 Length:12.25
 Optic Size: 5.50mm
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Capsular bag
 Material:
 One piece PMMA with UV absorbing optic and blue haptic

storz Model BV359
 Length 12.25 mm Power +04.0
 SN: x00401

This certifies that _____ received the following intraocular lens implant in the _____
 Left eye
 Right eye

BV359 Posterior Chamber IOL
Diopter: +04.0

Optic size: 5.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Capsular bag
 Material: One piece PMMA with UV absorbing optic and blue haptic


Physician: _____

BVR-150S storz +04.0
Posterior Chamber IOL
SN: x00701 Length:12.25
 Optic Size: 5.0mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Capsular bag
 Material:
 One piece PMMA with UV absorbing optic and blue haptics

storz Model BVR-150S
Length 12.25 mm Power +04.0
SN: x00501

MODEL BVR-150S	LENGTH 12.25mm	EXP. DATE 2003-02
DIOPTER +04.0	OPTIC 5.00mm	SN x00701
REF # KY040 LOT OTEST		
SN x00701		
"A" Factor : P.C. 118.1 Optic Design : Equiconvex Haptic Design : Modified C Material : One piece PMMA with UV absorbing optic and blue haptics Caution: Federal [USA] Law restricts this device to sale by or on the order of a physician.		
REF # KY040		
[01] 10757770068699 [21] x00701		

This certifies that _____ received the following intraocular lens implant in the _____ eye.
 Left eye
 Right eye
BVR-150S Posterior Chamber IOL
Diopter: +04.0
SN: x00701
 Optic size: 5.0mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Capsular bag
 Material: One piece PMMA with UV absorbing optic and blue haptics

Physician: _____

BVR-155S storz +04.0
Posterior Chamber IOL
SN: x01001 Length:12.25
 Optic Size: 5.50mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Capsular bag
 Material:
 One piece PMMA with UV absorbing optic and blue haptics

storz Model BVR-155S
Length 12.25 mm Power +04.0
SN: x01001

MODEL BVR-155S	LENGTH 12.25mm	EXP. DATE 2003-02
DIOPTER +04.0	OPTIC 5.50mm	SN x01001
REF # K1040 LOT OTEST		
"A" Factor : P.C. 118.1 Optic Design : Equiconvex Haptic Design : Modified C Material : One piece PMMA with UV absorbing optic and blue haptics Caution: Federal [USA] Law restricts this device to sale by or on the order of a physician.		
REF # K1040		
[01] 10757770069306 [21] x01001		

This certifies that _____ received the following intraocular lens implant in the _____ eye.
 Left eye
 Right eye
BVR-155S Posterior Chamber IOL
Diopter: +04.0
SN: x01001
 Optic size: 5.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Capsular bag
 Material: One piece PMMA with UV absorbing optic and blue haptics

Physician: _____

BVR-160M Storz +04.0
Posterior Chamber IOL
SN: x00501 Length:12.75

Optic Size: 6.00mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Capsular bag
 Material:
 One piece PMMA with UV absorbing optic and blue haptics

storz Model BVR-160M
 Length 12.75 mm Power +04.0
 SN: x00501

MODEL BVR-160M	LENGTH 12.75mm	OPTIC 6.00mm	EXP. DATE 2003-02
DIOPTER +04.0	OPHTHALMICS		SN x00501
SN x00501		REF # K5040 LOT OTEST	
"A" Factor : P.C. 118.1 Optic Design : Equiconvex Haptic Design : Modified C Material : One piece PMMA with UV absorbing optic and blue haptics Caution: Federal [USA] Law restricts this device to sale by or on the order of a physician.			
REF # K5040			
[01] 10757770071743 [21] x00501			

50005 Rev A

BVR-165L Storz +04.0
Posterior Chamber IOL
SN: x00201 Length:13.25

Optic Size: 6.50mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Ciliary sulcus or capsular bag
 Material:
 One piece PMMA with UV absorbing optic and blue haptics

storz Model BVR-165L
 Length 13.25 mm Power +04.0
 SN: x00201

MODEL BVR-165L	LENGTH 13.25mm	OPTIC 6.50mm	EXP. DATE 2003-02
DIOPTER +04.0	OPHTHALMICS		SN x00201
SN x00201		REF # YA040 LOT OTEST	
"A" Factor : P.C. 118.1 Optic Design : Equiconvex Haptic Design : Modified C Material : One piece PMMA with UV absorbing optic and blue haptics Caution: Federal [USA] Law restricts this device to sale by or on the order of a physician.			
REF # YA040			
[01] 10757770080868 [21] x00201			

50007 Rev A

This certifies that _____ received the following intraocular lens implant in the _____
 [] Left eye [] Right eye
BVR-160M Posterior Chamber IOL
Length: 12.75 mm
SN: x00501
 Optic size: 6.00mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Capsular bag
 Material: One piece PMMA with UV absorbing optic and blue haptics

Physician: _____

This certifies that _____ received the following intraocular lens implant in the _____
 [] Left eye [] Right eye
BVR-165L Posterior Chamber IOL
Length: 13.25 mm
SN: x00201
 Optic size: 6.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Ciliary sulcus or capsular bag
 Material: One piece PMMA with UV absorbing optic and blue haptics

Physician: _____

EZE-50 Storz +04.0
Posterior Chamber IOL
SN: x00801 Length:12.75

Optic Size: 5.00mm
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Capsular bag
 Material:
 PMMA with UV absorber

storz Model EZE-50
 Length 12.75 mm Power +04.0
 SN: x00901

MODEL
EZE-50
 DIOPTER
+04.0

I-x mm-L LENGTH
12.75mm
 OPTIC
5.00mm

EXP. DATE
2003-02
 SN
x00901

SN
x00901

REF # **[LOT]**
8M040 0TEST

"A" Factor : P.C. 118.1
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA with UV absorber

**Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.**

REF # **8M040**

[01] 10757770148346 [21] x00901

50025 Rev A


EZE-55 Storz +04.0
Posterior Chamber IOL
SN: x01101 Length:12.75

Optic Size: 5.50mm
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Capsular bag
 Material:
 PMMA with UV absorber

storz Model EZE-55
 Length 12.75 mm Power +04.0
 SN: x01101

MODEL
EZE-55
 DIOPTER
+04.0

I-x mm-L LENGTH
12.75mm
 OPTIC
5.50mm

EXP. DATE
2003-02
 SN
x01101

SN
x01101

REF # **[LOT]**
7W040 0TEST

"A" Factor : P.C. 118.1
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA with UV absorber

**Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.**

REF # **7W040**

[01] 10757770140227 [21] x01101

50027 Rev A

This certifies that _____ received the following intraocular lens implant in the _____
 Right eye Left eye
EZE-55 Posterior Chamber IOL
Diopter: +04.0
SN: **x01101**
 Optic size: 5.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Capsular bag
 Material: PMMA with UV absorber
 Physician: _____

EZE-65 Storz +04.0
Posterior Chamber IOL
SN: x00301 Length:13.25

Optic Size: 6.50mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Ciliary sulcus or capsular bag
 Material:
 PMMA with UV absorber

storz Model EZE-65
Length 13.25 mm Power +04.0
SN: x00301

MODEL EZE-65	LENGTH 13.25 mm	EXP. DATE 2003-02
DIOPTER +04.0	OPTIC 6.50mm	SN x00301
OPHTHALMICS		
REF # 8N040		LOT 8N040 0TEST
SN: x00301 "A" Factor : P.C. 118.1 Optic Design : Equiconvex Haptic Design : Modified C Material : PMMA with UV absorber Caution: Federal [USA] Law restricts this device to sale by or on the order of a physician. REF # 8N040		
 [01] 10757770148957 [21] x00301		

This certifies that _____ received the following intraocular lens implant in the _____
 Right eye Left eye
EZE-65 Posterior Chamber IOL
Length: 13.25 mm
SN: x00301

Optic size: 6.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Ciliary sulcus or capsular bag
 Material: PMMA with UV absorber

Physician: _____

P484UV Storz +04.0
Posterior Chamber IOL
SN: x00601 Length:13.00

Optic Size: 6.50mm Haptic angulation: 10 degrees
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Ciliary sulcus or capsular bag
 Material:
 PMMA with UV absorber

storz Model P484UV
Length 13.00 mm Power +04.0
SN: x01101

MODEL P484UV	LENGTH 13.00 mm	EXP. DATE 2003-02
DIOPTER +04.0	OPTIC 6.50mm	SN x01101
OPHTHALMICS		
REF # 5Q040		LOT 5Q040 0TEST
SN: x01101 "A" Factor : P.C. 117.9 Optic Design : Equiconvex Haptic Design : Modified C Material : PMMA with UV absorber Caution: Federal [USA] Law restricts this device to sale by or on the order of a physician. REF # 5Q040		
 [01] 10757770117939 [21] x01101		

This certifies that _____ received the following intraocular lens implant in the _____
 Right eye Left eye
P484UV Posterior Chamber IOL
Length: 13.00 mm
SN: x01101

Optic size: 6.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Ciliary sulcus or capsular bag
 Material: PMMA with UV absorber

Physician: _____


P517UV Storz +04.0
Posterior Chamber IOL
SN: x00901 Length:13.25

Optic Size: 6.50mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Ciliary sulcus or capsular bag
 Material:
 PMMA with UV absorber

Storz Model P517UV
 Length 13.25 mm Power +04.0
 SN: x01201

MODEL
P517UV
 DIOPTER
+04.0

EXP. DATE
2003-02
 SN
x00901

REF # **6L040** LOT **OTEST**

SN
x00901

"A" Factor : P.C. 118.0
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA with UV absorber

**Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.**

REF # **6L040**

[01] 10757770127709 [21] x00901

50109 Rev A


P574UV Storz +04.0
Posterior Chamber IOL
SN: x01201 Length:13.25

Optic Size: 6.00mm Haptic angulation: Step-vaulted
 Optic Design: Equiconvex
 Haptic Design:
 Modified C
 Placement: Ciliary sulcus or capsular bag
 Material:
 PMMA with UV absorber

Storz Model P574UV
 Length 13.25 mm Power +04.0
 SN: x01201

MODEL
P574UV
 DIOPTER
+04.0

EXP. DATE
2003-02
 SN
x01201

REF # **8J040** LOT **OTEST**

SN
x01201

"A" Factor : P.C. 118.1
 Optic Design : Equiconvex
 Haptic Design : Modified C
 Material : PMMA with UV absorber

**Caution: Federal [USA] Law restricts this device
 to sale by or on the order of a physician.**

REF # **8J040**

[01] 10757770147554 [21] x01201

50125 Rev A

This certifies that _____
 received the following intraocular lens implant in the
 Left eye
 Right eye

**P517UV Posterior Chamber IOL
 Length: 13.25 mm
 SN: x00901**

Haptic angulation: Step-vaulted
 Optic size: 6.50mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Ciliary sulcus or capsular bag
 Material: PMMA with UV absorber

Distribution:

This certifies that _____
 received the following intraocular lens implant in the
 Left eye
 Right eye

**P574UV Posterior Chamber IOL
 Length: 13.25 mm
 SN: x01201**

Haptic angulation: Step-vaulted
 Optic size: 6.00mm
 Optic design: Equiconvex
 Haptic design: Modified C
 Placement: Ciliary sulcus or capsular bag
 Material: PMMA with UV absorber

Physician initials: _____