

Qualification & Maintenance for LC System

Scott Liu Application Specialist Solution Center, Taiwan

© 2017 Waters Corporation COMPANY CONFIDENTIAL

USP 1058 Terminology



According to USP 1058.....

- Instruments are QUALIFIED, Was the instrument built right?
- Processes are VALIDATED, Was the right instrument built?
- AIQ (Analytical Instrument Qualification)
 - "AIQ Is documented evidence that an instrument performs suitably for its intended purpose and that it is properly maintained and calibrated"
 - Does not include people (training), processes performed on instruments (analytical methods)
 - Helps to justify the continued use of instrument, but it alone does not ensure quality of the data
- It is easiest to qualify instruments using standardized tests
- Procedures and methods are validated for each specific analysis

Performing AIQ



- Design Qualification (DQ)
- Installation Qualification (IQ)
- Operational Qualification (OQ)
- Performance Qualification (PQ)
- Change Control

Chosen because they are the most widely understood terms

THE SCIENCE OF

©2017 Waters Corporation COMPANY CONFIDENTIAL

Operational Qualification - Current

- OQ Parameters: Test operation of instrument as per specifications in the user's environment
 - Test critical parameters to assure required performance
 - Parameters based on manufacturer's recommendation and on user's intended use
 - Secure data storage, backup and archive process
- Performance of OQ
 - Use non-method specific testing
 - Repeat relevant OQ tests when instrument undergoes major repairs or modifications
 - Repeat tests on a regular basis

Operational Qualification – Proposed



THE SCIE

- OQ Parameters additions
 - Parameters to qualify described in the general chapters for the analytical technique used
 - Test critical functions under actual operating conditions
 - Demonstrate that the entire system, including software, works as intended
- Performance of OQ addition
 - OQ performed using the same software configuration as that used for routine analysis
- The best and easiest way to meet this requirement is using SystemsQT

©2017 Waters Corporation COMPANY CONFIDENTIAL

Performance Qualification - Current

- PQ tests are performed on a periodic basis to ensure that the instrument remains in qualified state after IQ/OQ
 - Run tests to check and verify satisfactory performance of the instrument
 - These tests can verify suitability of the instrument for your specific intended use and configuration
- Performance of PQ
 - Perform at specified intervals
 - Tests based on good science
 - Tests reflect general use of the instrument
 - Same tests can be used repeatedly to generate history of instrument performance
 - Specifications for PQ tests can be different than for OQ if required
 - Perform after repairs and maintenance
 - Maintain an SOP for operation/calibration/maintenance



- Instrumentation that fails must be investigated and an explanation provided with the reason for failure
- Addition of periodic review for critical instruments and computer systems

© 2017 Waters Corporation COMPANY CONFIDENTIAL

Change Control

- Changes are inevitable as new features are added and corrections are made
- Follow DQ/IQ/OQ/PQ classification process
 - DQ: Review the change; adopt only useful or necessary changes
 - $_{\odot}$ Not making changes for too long can also cause problems
 - Ex: Running unsupported versions of software or firmware
 - IQ: Install the changes
 - $_{\circ}\,$ Log any changes in the system log
 - OQ/PQ:
 - Revise tests and specifications if necessitated by the change
 - Change SOPs as necessary
 - Perform changed OQ or PQ tests



8







Annex 1 Specification Examples

RID detector

Examples of requirements for HPLC instruments and detectors								
Instrument module	Parameter to be checked	Typical tolerance limits						
Solvent delivery system	Flow rate Proportioning accuracy and precision (gradient test) Proportioning ripple	 ±5% ±2 ≤02% 						
Injector	Volume precision Carry-over	RSD ≤ 1.0 % see Annex I						
Autosampler	 Thermostatting accuracy and precision 	• ±3 °C						
Oven or cooling device	 Thermostatting accuracy 	• ±2°C						
UV/DAD detector	LinearityWavelength accuracy	 r² ≥ 0.999 ± 2 nm 						
Fluorescence detector	 Wavelength accuracy excitation Wavelength accuracy emission Sensitivity 	 ± 3 nm ± 3 nm see Annex I 						
Electrochemical detector	 Accuracy of the signal 	 see Annex I 						

see Annex I

see Annex I

± 0.1 mV/min

٠

Stability of the signal

Signal/Noise ratio

Drift over time

٠

•

Typical Qualification Tests for LC system

- Solvent Delivery System
 - Flow rate accuracy
 - Flow rate precision
 - Flow rate linearity
- Injector System
 - Injection accuracy
 - Injection linearity
 - Injector precision
 - Injector carryover
 - Sample position accuracy
- Detector System
 - Detector wavelength accuracy
 - Detector absorbance linearity
 - Detector noise and drift

COMPANY CONFIDENTIAL ©2017 Waters Corporation

Example 1 - HPLC Flow Linearity and Accuracy



	Linearity R^2 Specification	Linearity R*2 Linearity R*2 Accuracy Pass Fail Specification Accuracy		Accuracy Pass Fail	00		
1	R'2 > or = 0.90	1.00	PASS	< or = 100 ul	1	PASS	System Nam Projaci Nam
Г							Prosessing Report Meth User Name: Sample Sett
ſ	Accuracy Cale	ulation :	Flow rate Accur	acv = IY interd	ent/Slope	I * 1000	Sample

(to convert from mL to uL). The flow rate accuracy passes when the result is < or = 100 uL.

Y intercept = -0.001423	Slope = 1.240509		
Flow rate Accuracy =	-0.001423/ 1.240509 * 1	000 = 1.1	



THE SCIENCE OF WHAT'S POSSIBLE.

COMPANY CONFIDENTIAL ©2017 Waters Corporation

13

Waters

THE SCIEN







What do really know after these tests?



- Solvent Delivery System
 - motors work correctly
 - The gradient are working properly and the mixing valve (low pressure mixing) all actuate without sticking
- Injector System
 - syringe drive motors, valves and gears all work as intended.
- Detector System
 - gratings, mirrors, cell and lamp are functioning as intended in terms of directing light through the cell and that the correct wavelengths are used

©2017 Waters Corporation COMPANY CONFIDENTIAL



21

Recommended Qualification Procedure -Pump

Repair Type and Affected Device	Procedure
Major repair	
Pump motor driver PCBs	
• Gears	Requalify each parts
Piston drives	
Pump castings	
Pump control modules	
Minor repair	
Check valves	Static leak test
Check valve cartridges	Static leak test
• Plungers	Static leak test
• Seals	Static leak test
Pressure transducer	Verify proper zero and pressure increase
Gradient proportioning valves	Verify gradient valve proportioning or gradient operation
Column heater	Verify temperature accuracy
Column heater/cooler	Verify temperature accuracy

©2017 Waters Corporation COMPANY CONFIDENTIAL



23

'aters

Recommended Qualification Procedure -Injector

Repair Type and Affected Device Procedure Major repair Injector/valve driver PCBs Requalify each parts Injector assemblies Valve assemblies Minor repair · Sample positioning drives Verify sample positioing accuracy Pressure transducers Verify zero and deflection at pressure · Seal pack and needle Verify injector accuracy Valve seals Verify injector accuracy Sample heater/cooler Verify temperature accuracy Column heater Verify temperature accuracy Column heater/cooler Verify temperature accuracy

© 2017 Waters Corporation COMPANY CONFIDENTIAL



25

aters

Recommended Qualification Procedure – Optical detectors

Repair Type and Affected Device	Procedure
Major repair • Cell assembliers • Analog PCBs • Preamp PDBs • Photodiodes • Optics benches • Gratings • Grating drive devices • Mirrors • Beam splitters • Photodiode arrays • Photodiode arrays • Photomultiplier tubes • Optical slits • Filter wheel replacement • Drift tube • Nebulizer • Nebulizer heater/cooler	Requalify each parts
Minor repair • Lamps • Cell windows	Verify lamp calibration via start-up diagnostics Verify lamp calibration via start-up diagnostics

©2017 Waters Corporation COMPANY CONFIDENTIAL

Advantages of using CDS Managed Qualification



27

- Accurate Qualification Testing and Analysis
 - Less opportunity for human error
 - Measures peak areas, peak heights and retention times accurately and consistently
 - Custom field calculations and regression analysis
 - Testing consistent from system to system
 - Reduces time that system is off-line by about half
 - Multiple systems can be qualified at once
 - Qualifies software and systems in their analytical configuration
 - Qualifies using same peak processing and quantitation algorithms as during use on CDS system of record
 - Demonstrates system level fitness for chromatographic use

Empower SQT (System Qualification Tool)

Waters THE SCIENCE OF WHAT'S POSSIBLE!

29

Built into Empower

Secure and Auditable Qualification Data

- 21 CFR Part 11 Compliant Ready Qualification Data
- All of your data is maintained on your CDS
- On-Line Qualification Documentation for Easy Inspection
- Easy tracking and trending of qualification results
- Audit trails and method change control part of the data system
- No need for external spreadsheets or third party software
- Secure data environment
- On-line review and approval available

©2017 Waters Corporation COMPANY CONFIDENTIAL



Empower 3 Audit Trial Overview

Michelle Ho Informatics BD Manager – APAC June 2017

© 2017 Waters Corporation COMPANY CONFIDENTIAL

Disclaimer



- This presentation is for informational purposes only and should not be taken as advice regarding any particular course of action to be followed.
- Waters does not make any representations or warranties, express or implied, to any party, regarding use of the information contained in this presentation to make decisions regarding the implementation and maintenance of effective quality control systems and quality assurance testing programs, including but not limited to the applicable Good Manufacturing Regulations that apply to the manufacture of regulated products.

Empower Audit Trails

- Understand the Information in the Audit Trails and what causes an entry to be made in these Audit Trails
- Differentiate between the various Audit Trails
 - System
 - Project
 - Method
 - ∘ Sample
 - Sample Set
 - Result

System Audit Trail=specific to installation to Empower All data related Audit Trails are in the project

COMPANY CONFIDENTIAL ©2017 Waters Corporation

Vaters **Empower Audit Trails ID numbers** THE SCIEN Empower is built into an Oracle Database This database gives each object or result a Unique Identifier for tracking the values and records This identifier is unique within each project. Modification of any data base object results in a NEW record with NEW identifiers (Nothing is ever over written but versioned) Many users of Empower use these ID number to prove and identify results to auditors

Also to track for their own purposes



THE SCIENCE OF

Empower Audit Trail

- The built in Empower Database
 - Enables every object to be uniquely referenced
 - Can never overwrite data
 - Can never mistake which data went with which method
 - Ensures easy and accurate data review
- Automatic versioning for results / methods
 - With full computer generated audit trail
 - WHO changed WHAT (before and after values) WHEN.... And WHY?)

©2017 Waters Corporation COMPANY CONFIDENTIAL

dit Troil Sv	votom Audi	it Troil Water
uit Trail – Sy	Sterri Audi	IL II dili THE SCIENCE OF WHAT'S
島 EMPOWER3 WORLD as System/Admini	strator - Configuration Manager	
File Edit View Records Tools He	lo	
	Filter By:	Default Edit View Update Max Rows 1000 J4 4
Empower 3 Configuration	Action	D
Projects	1 Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
Aldol_Condensation	2 Unsuccessful Logon Attempt	User: system Node: Democlientfr1 - Authentication Failed
Bio Clinical	3 Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
EIN Project	4 Successfully Logged Off	User: System Node: Democlientfr1
Gimeniride PATROL	5 Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
GLYCAN March 2013	6 Modified User	User: System Reason: Set password (system generated)
Impurity Assay	7 Successfully Logged Off	User: System Node: Democlientfr1
LIMS_LabSystems	8 Altered System Policy	Default Full Audit Trail Confirm Identity on Deletion Changes - False> True Reason: added policy
🕂 🔄 Mass Spectrometry	9 Altered System Policy	Disable Copy from non Full Audit Trail project to Full Audit Trail project - False> True Reason: added policy
🕀 🕞 Pharmaceutical	10 Altered System Policy	Passwords Expire - False> True Reason: added policy
SQT BEAE8EBFDX41	11 Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
SQT_ACQUITY_PDA	12 Unsuccessful Logon Attempt	User: system Node: Democlientfr1 - Authentication Failed
SQT_ACQ_DIAIpha1_FLK	13 Successfully Logged Off	User: System Node: Democlientfr1
All SWI Demo	14 Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
TTAcetaminophen	15 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
- (국) Valid_OQ1	16 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
Valid_OQApex	17 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
Valid_OQMVM	18 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
	19 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
	20 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
E & Libraries	21 Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
eCord	22 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
V Users	23 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
©? User Types	24 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
All Disto Tumor	25 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
System Audit Trail	26 Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
🖏 Offline System Audit Trail	27 Successfully Logged Off	User: System Node: Democlientfr1
	28 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
	29 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
	30 Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
	24 0	Hann Conten Made Demotionality Johnstone Ten 8/2



What is System Audit Trail

- System Audit Trail keeps track of actions taken at the system level.
- The audit trail is divided into the following fields:
 - Action
 - Details
 - Date and time stamp as the when the action took place
 - User Who took the action
- The system audit trail shows changes to system objects and system policies
 - details archive activity
 - notes all changes to security (users, user types etc)
 - documents all successful and unsuccessful logins
 - $_{\odot}\,$ you have a history of who was logged into the application at any time
 - 。 you have information about system break in attempts
 - o includes the client the login/login attempt occurred at

©2017 Waters Corporation COMPANY CONFIDENTIAL

System Audit Trail examples...



8

THE SCIENCE OF

131	Action	
1	Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
2	Unsuccessful Logon Attempt	User: system Node: Democlientfr1 - Authentication Failed
3	Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
4	Successfully Logged Off	User: System Node: Democlientfr1
5	Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
6	Modified User	User: System Reason: Set password (system generated)
7	Successfully Logged Off	User: System Node: Democlientfr1
8	Altered System Policy	Default Full Audit Trail Confirm Identity on Deletion Changes - False> True Reason: added policy
9	Altered System Policy	Disable Copy from non Full Audit Trail project to Full Audit Trail project - False> True Reason: added policy
10	Altered System Policy	Passwords Expire - False> True Reason: added policy
11	Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
12	Unsuccessful Logon Attempt	User: system Node: Democlientfr1 - Authentication Failed
13	Successfully Logged Off	User: System Node: Democlientfr1
14	Successfully Logged On	User: System Node: Democlientfr1 Interface: Pro
15	Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
16	Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit
17	Successfully Logged On	User: System Node: Democlientfr1 Interface: ToolKit

Empower Project Audit Trails

The Project audit trail keeps track of action taken within the project. Each project has its own unique audit trail.

Filte	Fiter By: Default Edit View Update Max Rows: 1000								
•	L Sample Sets Injections Channels Methods Result Sets Results Peaks Fractions Sign Offs Curves View Filters Custom Field								
E	Action	Details	Change Date						
1	Created Manual Result	Sample Name: S160920_00126 Vial: 5 Injection No.: 1 Channet 2487Channel 1 Method: Acetaminophen PM Result ID: 2611 Channel ID: 1248 Result Set ID: 2587	Wednesday, November 23, 2016 4						
2	Updated Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2595 Calibration Source: Auto	Wednesday, November 23, 2016 4						
3	Created Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2593 Calibration Source: Auto	Wednesday, November 23, 2016 4						
4	Created Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2595 Calibration Source: Auto	Wednesday, November 23, 2016 4						
5	Created Result Set	Result Set: ELN_Assay Sample Set Method: 8536_02 071402 Method: Acetaminophen PM Processed How: Processing Method Result Set ID: 2587 Sample Set ID: 1804	Wednesday, November 23, 2016 4						
6	Created Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2588 Calibration Source: Auto	Wednesday, November 23, 2016 4						
7	Modified Method	Method: Acetaminophen PM Type: Processing Version: 5	Wednesday, November 23, 2016 4						
8	Created Manual Result	Sample Name: S160920_00126 Vial: 5 Injection No.: 1 Channel: 2487Channel 1 Method: Acetaminophen PM Result ID: 2561 Channel ID: 1248 Result Set ID: 2539	Tuesday, November 22, 2016 6:32						
9	Updated Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2545 Calibration Source: Auto	Tuesday, November 22, 2016 6:24						
10	Created Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2543 Calibration Source: Auto	Tuesday, November 22, 2016 6:24						
11	Created Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2545 Calibration Source: Auto	Tuesday, November 22, 2016 6:24						
12	Created Result Set	Result Set: ELN_Assay Sample Set Method: 8536_02 071402 Method: Acetaminophen PM Processed How: Processing Method Result Set ID: 2539 Sample Set ID: 1804	Tuesday, November 22, 2016 6:24						
13	Created Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2537 Calibration Source: Auto	Tuesday, November 22, 2016 6:23						
14	Created Manual Result	Sample Name: STD-3 Vial: 9 Injection No.: 1 Channel: 2487Channel 1 Method: Acetaminophen PM Result ID: 2535 Channel ID: 1260 Result Set ID: 0	Tuesday, November 22, 2016 6:23						
15	Modified Method	Method: Acetaminophen PM Type: Processing Version: 4	Tuesday, November 22, 2016 6:23						
16	Created Manual Result	Sample Name: S160920_00126 Vial: 5 Injection No.: 1 Channel: 2487Channel 1 Method: Acetaminophen PM Result ID: 2507 Channel ID: 1248 Result Set ID: 2467	Thursday, November 17, 2016 5:3						
17	Created Manual Result	Sample Name: STD-2 Vial: 3 Injection No.: 1 Channel: 2487Channel 1 Method: Acetaminophen PM Result ID: 2500 Channel ID: 1238 Result Set ID: 0	Thursday, November 17, 2016 5:3						
18	Created Manual Result	Sample Name: Blind Vial: 1 Injection No.: 1 Channet: 2487Channel 1 Method; Acetaminophen PM Result ID: 2503 Channel ID: 1226 Result Set ID: 2467	Thursday, November 17, 2016 5:3						
19	Updated Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2501 Calibration Source: Auto	Thursday, November 17, 2016 5:3						
20	Updated Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2473 Calibration Source: Auto	Tuesday, September 20, 2016 9:3						
21	Created Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration D: 2471 Calibration Source: Auto	Tuesday, September 20, 2016 9:3						
22	Created Calibration	System: B101R319_Waters1 Method: Acetaminophen PM Channel: 2487Channel 1 Calibration ID: 2473 Calibration Source: Auto	Tuesday, September 20, 2016 9:3						
23	Created Result Set	Result Set: ELN_Assay Sample Set Method: 8536_02 071402 Method: Acetaminophen PM Processed How: Processing Method Result Set D: 2467 Sample Set D: 1804	Tuesday, September 20, 2016 9:3						
24	Created Method	Method: Michele0920 Type: Sample Set Version: 1	Tuesday, September 20, 2016 9:30						

©2017 Waters Corporation COMPANY CONFIDENTIAL





THE SCIENCE OF

- Sample Audit Trail
 - Tracks changes to entered data about each sample
- Result Audit Trail
 - Links results to instruments, samplesets, methods, calibration curves and standards used in calibration.
 - Also traces any manual manipulation of data
- Method Audit Trail
 - Keeps all versions of method for recreation of results
 - Audit Trail monitors each change, before and after values, who when and why
 - Different versions can be compared to identify the differences

Built in Audit Trails in Empower

- All user actions are logged in various audit trails and associated with the logged in USERNAME
 - Assumes all users have unique User Account
- It is not possible to create, manipulate, modify or delete data inside Empower without creating an audit trail entry
- Multiple "modes" of audit trail
 - Silent
 - Full Includes the requirement to enter a reason "Why?"
 - $_{\circ}~$ With free form reasons
 - $_{\circ}~$ With predefined reasons only
 - Reauthentication (re entry of password to confirm identity)
- Empower Audit trails are not editable or modifiable by ANY USER

© 2017 Waters Corporation COMPANY CONFIDENTIAL

	Empower Sample Aud	it Trail		Waters THE SCIENCE OF WHAT'S POSSIBLE!					
	C Alliance UV System in Pharmaceutical_Project as CWeis/Manager - QuickStart - [Find Data]								
	File Edit View Tools Database Manage Help Run Samples Sample Queue Control Panel Sample History	Detailed Unicetion Info ceutical_Project as (20172) 2010 100 100 100 100 100 100 100 100 10	CWeis/Manager	- Alter Sample					
	Browse Proje File Help View Data I User: mharnois Date: 5/3/2002 12:49:16 // 2 Modified Vial(Lot_No): <no value=""> - 3 User: TBrown Date: 5/3/2002 12:52:31 A 4 Modified Vial(SampleVveight): 1.0000</no>	Sample History AM Reason: All sample > SS 4Comp M Reason: Wrong wie 0 -> 23.40000	e information not ki ight values entere	A Fhan A Phan A Phan					
/	Mathad Set			A Phan					
3	User: TBrown Date: 5/3/2002 12:52:31	AM Reason: \	Wrong wei	ght values entered at time of analysis.					
4	Modified Vial(SampleVVeight): 1.000	000 -> 23.4000	00						
	49 SystemSuitability SS 4 50 SystemSuitability SS 4 View Acculisition Image: Signal State	how All Columns 5 rint Table 6 able Properties	Standard Standard	Friday, November 30, 2001 3:17:21 PM Phar Friday, November 30, 2001 3:14:26 PM Phar					



Audit Trial summary



- There are a variety of Audit Trails through out Empower keeping track of critical activities.
- As a Data reviewer it is important to understand the information tracked in these audit trails and what caused entries to be made.
- Comments entered should reflect why a change was made. The audit trails already track what was changed.

Result Audit Viewer

- Understand how to use the Result Audit Viewer. This tool is designed to help the user review all the history that contributed to the generation of one Empower result.
- Access audit trails information for results in one window.
- As more companies are reviewing data electronically the workflow would be to use the Result audit reviewer before sign offs.
- Electronic Sign-off would be used to confirm that the results had been reviewed.





15

Result Audit	Viewer Too	l					THE S	Waters	
Di Inform2012_Tutorial2_Process as Employee 00789/Administrator - OuickStart - (Review D									
File Edit View Tools Plot Proces	Analgesics in Inform2012_Tutorial2 File View Help	Process as Employee_00789/	Administr	ator - Resu	ilt Audit V	iewer			
Alter Sample	Alter Sample 關								
Sample Queue Result Audit View Results									
Control Panel	Result Sample Name	Manual Result Comments		ummary	Result	Result		A .	
	L3 Id Semple Name	Manual Result Comments	duits	Faults	#	Superseded			
Browse Project	8 1160 AG Standard 3				10		njection Volume = 2.00	Acetaminophen Value = 31.250000	
View bain	9 1161 AG Standard 4 10 1162 AG Standard 5				10 10		njection Volume = 2.00 njection Volume = 2.00	Acetaminophen Value = 34.400000 Acetaminophen Value = 37.500000	
View Method B- 0 / Method Set	Result History Result Differences Pr	ocessing Method Sample Set N	lethod In	strument M	ethod Me	Afte	er This Date: 1/ 9	>/2000 ▼ Update	
Processing R A		la durante Malla a Uda 4000	Reason	User	7050044	Date	Action Type	Source	
	1 uto Additions : Injection Id : 1087	Instrument Method Id : 1063	N/A	System	7/25/2011	2:21:45 PM CEs 2:28:35 PM	ST N/A	Sample Set Method Properties	
One	Stop Solution	:		System	7/22/2011	1:50:32 PM	N/A	Instrument Method Properties	
One	Stop Solution	-		Rune	6/17/2011	7:57:13 AM	N/A	Processing Method Properties	
Project Aug	dit Trails			Rune	6/16/2011	11:07:47 AM	N/A	Processing Method Properties	
		-		Rune	6/16/2011	10:08:34 AM	N/A	Processing Method Properties	
 Method His 	story and Differ	ences		Rune	6/16/2011	10:01:33 AM	N/A	Instrument Method Properties	
 Sample His 	story	-		System	0/13/2011	3.13.30 PM	108	Processing Method Properties	
Injection L ©2017 Waters Corporation COMPA	OG ANY CONFIDENTIAL			_	_	٨	lew in Fea	ature Release 2	
Result Audit	Viewer						THE S	Waters	
Image: Collog_2695_996 in ICH Impurity Processing as Syzeries File Edit View Tools Plot Process Mavigate Image: Control Panel Browse Project Oze Oze Oze Deak List: Acid Degradants Glimepiride Impurity 1 Impurity 2 Peak in blank Total Impurities Show Me Oze	ten/Administrator - QuickStart - (Review Data Ma Qptions Window Spectrum Review Library P O Result 14000+002 [140000 Result Name 3984 Impurity 2 3983 Impurity 2	in Window) Manage Help D D D C Q @ Protocol QPM0879 in ICH Impu File View Help Result Differen Area (µV*se	ity Processing ces c) 4740 7168	g as System/A Results	dministrator	- Result Audit View 052 Rb 079 Vv * Supersi	er D 20 Int Type	Differences in the results are in blue. Results outside limits are in red Manual integration	
View Acquisition 0.04 + 0.04 + 0.01 + 0	10 3.60 3.80 4.00 4.20 Minules — Time Pt: 28 Amount Time minimation Area % Area % Adjust 0.003 1.764 937 0.09 2.072 2517 0.24 2.199 2865 0.27 2.289 4828 0.46 2.312 9068 0.86 2.2754 2311 0.22 2.814 1303 0.12 2.860 4740 0.45 2052 2.860 4740 0.052 2.860 4740 20 Channl •	Image: second	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	Res lame 2 y y 2	··· I ↓ □ ↓ □ ↓ □ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓	30 Image: Constraint of the second seco	Injecton volum Injecton Volum Injecton Volum Injecton Volum Int Type 2 Rb 2 Rb 9 Vv	 noted in Integration Type field 	

Result Audit Viewer Summary

- Result Audit Viewer brings together all the audit history related to the result set in view. Audit history is pulled from various locations in the software which were explored in Chapter 4
- Result Audit Viewer additionally aids the user in comparing result changes and identify manual manipulations/faults in a result set.

©2017 Waters Corporation COMPANY CONFIDENTIAL



19

