

# Method of Test for Cu-pyropheophytin A in Edible Oils

## 1. Scope

This method is applicable to the determination of Cu-pyropheophytin a, the major copper chlorophyll, in edible oils.

## 2. Method

After extraction and purification, Cu-pyropheophytin a is determined by high-performance liquid chromatography (HPLC) and confirmed by liquid chromatography/high resolution tandem mass spectrometry (LC/HRMS<sup>2</sup>) or liquid chromatography/ tandem mass spectrometry (LC/MS/MS).

### 2.1. Equipments

#### 2.1.1. High-performance liquid chromatograph

2.1.1.1. Detector: photodiode array detector.

2.1.1.2. Column: GL Sciences InertSustain C18, 2- $\mu$ m, 2.1-mm i.d.  $\times$  10-cm, or an equivalent product.

#### 2.1.2. Liquid chromatograph/high resolution tandem mass spectrometer

2.1.1.1. Ion source: negative ion atmospheric pressure chemical ionization, APCI<sup>-</sup>.

2.1.1.2. Column: HALO C18, 2.7- $\mu$ m, 4.6-mm i.d.  $\times$  10-cm, or an equivalent product.

#### 2.1.3. Liquid chromatograph/tandem mass spectrometer

2.1.1.1. Ion source: negative ion atmospheric pressure chemical ionization, APCI<sup>-</sup>.

2.1.1.2. Column: HALO C18, 2.7- $\mu$ m, 4.6-mm i.d.  $\times$  7.5-cm, or an equivalent product.

#### 2.1.4. Solid phase extraction vacuum manifolds.

#### 2.1.5. Nitrogen evaporator.

### 2.2. Chemicals

Petroleum ether, HPLC grade;

Ethyl ether, HPLC grade;

Acetone, HPLC grade;

Methanol, HPLC grade;

Anhydrous ethanol, HPLC grade;

Ammonium acetate, GR grade;

Milli-Q water, resistivity  $\geq$  18 M $\Omega$   $\cdot$  cm (25  $^{\circ}$ C);

Cu-pyropheophytin a, reference standard.

### 2.3. Apparatus

2.3.1. Centrifuge tube: 15-mL, PP.

2.3.2. Membrane filter: 0.22- $\mu$ m, Nylon.

2.3.3. Solid phase extraction cartridge: Sep-Pak<sup>®</sup> silica, 1-g, 6-mL, or an equivalent product.

## 2.4. Reagents

### 2.4.1. 1 M ammonium acetate solution

Dissolve 77 g of ammonium acetate in Milli-Q water and dilute with Milli-Q water to 1000 mL.

### 2.4.2. Mobile phase for HPLC

Solvent A: methanol: 1 M ammonium acetate (8:2, v/v).

Solvent B: acetone.

Solvent C: methanol.

Solvent D: Milli-Q water.

### 2.4.3. Mobile phase for LC/HRMS<sup>2</sup>

Solvent A: acetone.

Solvent B: methanol.

### 2.4.4. Mobile phase for LC/MS/MS

Solvent A: methanol.

Solvent B: anhydrous ethanol.

## 2.5. Standard solution preparation

Transfer about 10 mg of Cu-pyropheophytin a reference standard accurately weighed into 10-mL volumetric flask, dissolve and dilute to volume with acetone as a stock solution. When to use, dilute the stock solution with acetone to 0.05 - 5.0 µg/mL as standard solutions.

## 2.6. Sample solution preparation

Transfer about 1 g of the sample accurately weighed into a centrifuge tube and transfer the sample from the tube into a Sep-Pak<sup>®</sup> silica cartridge by a plastic dropper. Wash the centrifuge tube with 3 mL of petroleum ether and transfer the washing solution into the cartridge. Repeat the-washing step twice. Discard the eluents. Wash the cartridge with 9 mL of petroleum ether: ethyl ether (9:1, v/v) and discard the eluent. Add 6 mL of acetone to the cartridge, collect the eluent and evaporate to dryness by gently flushing with a stream of nitrogen at 25°C. Dissolve the residue with 1 mL of acetone, then filter the solution with a membrane filter, and use it as the sample solution.

## 2.7. Identification and quantification

Separately inject about 20 µL of the sample solution and the standard solutions into the HPLC and perform HPLC analysis. Identify Cu-pyropheophytin a by retention times and absorption spectra. Calculate the amount of Cu-pyropheophytin a in the sample by the following formula:

$$\text{The amount of Cu-pyropheophytin a in the sample (ppm)} = \frac{C \times V}{M}$$

where

C: the concentration of Cu-pyropheophytin a in the sample solution calculated by the standard curve ( $\mu\text{g/mL}$ )

V: the make up volume of sample (mL)

M: the weight of sample (g)

HPLC operating conditions:

Photodiode array detector: 430 nm.

Column: GL Sciences InertSustain C18, 2  $\mu\text{m}$ , 2.1-mm i.d.  $\times$  10-cm.

Column temperature: 30  $^{\circ}\text{C}$ .

Mobile phase: gradient.

Time (min)	Solvent A (%)	Solvent B (%)	Solvent C (%)	Solvent D (%)
0.0 $\rightarrow$ 3.0	30 $\rightarrow$ 30	0 $\rightarrow$ 0	0 $\rightarrow$ 0	70 $\rightarrow$ 70
3.0 $\rightarrow$ 7.0	30 $\rightarrow$ 100	0 $\rightarrow$ 0	0 $\rightarrow$ 0	70 $\rightarrow$ 0
7.0 $\rightarrow$ 9.0	100 $\rightarrow$ 100	0 $\rightarrow$ 0	0 $\rightarrow$ 0	0 $\rightarrow$ 0
9.0 $\rightarrow$ 11.0	100 $\rightarrow$ 50	0 $\rightarrow$ 25	0 $\rightarrow$ 25	0 $\rightarrow$ 0
11.0 $\rightarrow$ 13.0	50 $\rightarrow$ 50	25 $\rightarrow$ 25	25 $\rightarrow$ 25	0 $\rightarrow$ 0
13.0 $\rightarrow$ 15.0	50 $\rightarrow$ 0	25 $\rightarrow$ 50	25 $\rightarrow$ 50	0 $\rightarrow$ 0
15.0 $\rightarrow$ 25.0	0 $\rightarrow$ 0	50 $\rightarrow$ 50	50 $\rightarrow$ 50	0 $\rightarrow$ 0
25.0 $\rightarrow$ 26.0	0 $\rightarrow$ 0	50 $\rightarrow$ 70	50 $\rightarrow$ 30	0 $\rightarrow$ 0
26.0 $\rightarrow$ 35.0	0 $\rightarrow$ 0	70 $\rightarrow$ 90	30 $\rightarrow$ 10	0 $\rightarrow$ 0
35.0 $\rightarrow$ 35.1	0 $\rightarrow$ 30	90 $\rightarrow$ 0	10 $\rightarrow$ 0	0 $\rightarrow$ 70
35.1 $\rightarrow$ 40.0	30 $\rightarrow$ 30	0 $\rightarrow$ 0	0 $\rightarrow$ 0	70 $\rightarrow$ 70

Injection volume: 20  $\mu\text{L}$ .

Flow rate: 0.25 mL/min.

## 2.8. Confirmation

### 2.8.1. LC/HRMS<sup>2</sup>

Separately inject about 40  $\mu\text{L}$  of the HPLC positive sample solution and the standard solution into the LC/HRMS<sup>2</sup> and operate according to the following LC/HRMS<sup>2</sup> conditions. Identify Cu-pyropheophytin a based on retention times, mass accuracy ( $< 5$  ppm) and relative ion intensities.

LC/HRMS<sup>2</sup> operating conditions:

Column: HALO C18, 2.7- $\mu\text{m}$ , 4.6-mm i.d.  $\times$  10-cm.

Column temperature: 30  $^{\circ}\text{C}$ .

Mobile phase: gradient.

Time (min)	Solvent A (%)	Solvent B (%)
0.0 $\rightarrow$ 10.0	5 $\rightarrow$ 80	95 $\rightarrow$ 20
10.0 $\rightarrow$ 11.0	80 $\rightarrow$ 100	20 $\rightarrow$ 0

11.0 → 14.0	100 → 100	0 → 0
14.0 → 14.1	100 → 5	0 → 95
14.1 → 21.0	5 → 95	95 → 95

Injection volume: 40  $\mu$ L.

Flow rate: 1 mL/min.

Ion source: negative ion atmospheric pressure chemical ionization, APCI<sup>-</sup>.

Collision energy: 20 eV.

Detection mode: product ion scan.

Resolution: 70000.

Analyte	Precursor ion ( <i>m/z</i> )	Product ion ( <i>m/z</i> )
		522.1468
Cu-pyropheophytin a	873.4749	550.1799
		594.1697

### 2.8.2. LC/MS/MS

Separately inject about 40  $\mu$ L of the HPLC positive sample solution and the standard solution into the LC/MS/MS and operate according to the following LC/MS/MS conditions. Identify Cu-pyropheophytin a based on retention times and relative ion intensities.

LC/MS/MS operating conditions:

Column: HALO C18, 2.7- $\mu$ m, 4.6-mm i.d.  $\times$  7.5-cm.

Column temperature: 30  $^{\circ}$ C.

Mobile phase: gradient.

Time (min)	Solvent A (%)	Solvent B (%)
0.0 → 1.0	100 → 70	0 → 30
1.0 → 6.0	70 → 40	30 → 60
6.0 → 14.0	40 → 0	60 → 100
14.0 → 15.5	0 → 0	100 → 100
15.5 → 16.0	0 → 100	100 → 0
16.0 → 20.0	100 → 100	0 → 0

Injection volume: 40  $\mu$ L.

Flow rate: 1 mL/min.

Ion source: negative ion atmospheric pressure chemical ionization, APCI<sup>-</sup>.

Curtain gas: 20 psi.

Collision gas: high.

Gas 1: 55 psi.

Gas 2: 0 psi.

Temperature: 400  $^{\circ}$ C.

Detection mode: multiple reaction monitoring (MRM).

Analyte	Precursor ion ( <i>m/z</i> ) > product ion ( <i>m/z</i> )	Declustering potential (V)	Collision energy (eV)
Cu-pyropheophytin a	873.5 > 522	-68	-50
	873.5 > 535	-68	-57
	873.5 > 550	-68	-50
	873.5 > 594	-68	-38

Notes:

1. All the parameters can be adjusted depending on the instruments used if the above conditions are not applicable.
2. Maximum permitted tolerances for relative ion intensities by LC/HRMS<sup>2</sup> or LC/MS/MS are as follows:

Relative ion intensity (% of base peak)	Tolerance (%)
> 50	± 20
> 20~50	± 25
> 10~20	± 30
≤ 10	± 50

3. Further validation is necessary when interference compounds appear in samples.