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E-mail: sales@macroitd.com.tw

Compiled by Sean Kuo SEP. 2011, version A

MEDICAL GAS ANALYSIS & QUALITY CONTROL SYSTEM (MGAQ)

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Compiled by Sean Kuo Oct. 2010, version A

MEDICAL GAS ANALYSIS & QUALITY CONTROL (MGAQ)

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MEDICAL OXYGEN (O₂) MGAQ MAINTENANCE & CALIBRATION

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
Figure 1.1 Key features of Xentra

Key:	A	FRONT VIEW	8	Sample inlet(s)
	B	REAR VIEW	9	Mains power connector
	1	Sample filter (optional)	10	Fan and filter
	2	Flowmeter(s) (optional)	11	Sample outlet(s)
	3	Display	12	Functional earth
	4	Keypad	13	Serial output port
	5	Display adjustment	14	Signal terminals
	6	Needle valve(s) (optional)	15	Screen
	7	Rack mounting brackets		

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
GfX Bench - exploded view (Showing location of Fuses)

On older units F3 was F1.6A.

If working on one of these earlier units fuse F3 should be changed to F630mA

- Fuse F1 (T3.15A) Sample cell heater.
- Fuse F2 (T1.6A) Chopper wheel.
- Fuse F3 (F630mA) IR source power*.
- Fuse F4 (F2.0A) Chopper box heater.

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
LED D13 (Motor OK)

LED D16 (Bench OK)

(This LED will turn off if any other fault LED extinguishes).

Picture of LED locations on housekeeping P.C.B. (01210905).

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- If open circuit, replace fuse.


If OK,

- Switch on bench, connect voltmeter to inner pins (2 & 3) of PL4 (source voltage sense). Check voltage is 2.0volts dc +/- 0.02v

If not OK, replace housekeeping pcb.

If OK,

- Switch off power, disconnect IR source from PL4. Check source resistance between inner pins (2 & 3). Resistance should be <2 ohms. If outside of these limits, replace source.

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
Check fuse F2 (T1.6A).

If OK,

Check the chopper wheel is free to rotate, not hitting or grazing any part of the chopper box.

Tip.

If ambient noise is low, listen carefully near chopper box to see if you can hear wheel rotation without any knocks or scrapping sounds.

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




01210902 P.C.B. details

Note...

Be cautious if removing spacer, as damage can occur to capacitors C29 & C30. Use a 5.5mm socket to loosen/tighten spacer.

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
Gas signal adjustments

*(Under zero conditions
i.e. Nitrogen in sample Cell).*

- Gas Sig (V_{gas}) TP4(+) wrt TP5 (0VA)
- Gas Sig (V_{N_2}) TP3 (+) wrt TP5 (0VA)
- Diff Sig (V_{diff}) TP2 (+) wrt TP5 (0VA)

Signal Processing P.C.B. test points

Location of RV1 – Detector pre-amp Gain adjustment

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
Detector signal from TP1 of signal processing P.C.B

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1. SPARES KIT, 2 YEARS OF OPERATION
2. FILTER ELEMENT, FAN,
3. SENSOR INTERFACE PCB,
4. INVERTER 24 VDC TO 600 V
5. CCFL BACKLIGHT ASSEMBLY, FOR BLUE/WHITE DISPLAY
6. SWITCHED MODE POWER SUPPLY
7. TUBING / FITTINGS REFURB. KIT

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1. SCRUBBER ASSEMBLY
2. KIT SCRUBBER SACHET
3. KIT, I.R. SOURCE
4. KIT, 1210 FUSES

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
Principle of Operation

Disadvantages 缺點

- ◆ P₂O₅ coating can be destroyed when on UHP gases, need to wet up periodically, 五氧化二磷感測器在長時間使用於高純度氣體中需要定期對感測電極塗料
- ◆ Problems when H₂ and O₂ are present in sample gas due to recombination, 若氫氣和氧氣重組時會造成錯誤的水份讀值
- ◆ Needs precise flow control, 流量控制的誤差會造成準確度大幅降低
- ◆ High moisture levels will overload sensor, 高水含量會造成感測器受損

MEECO AQUAVOLT+ 水份分析儀

MEECO AQUAVOLT+ 規格

Range 量測範圍	0-20 ppmV
Accuracy 精確度	± 2% of reading or 20 ppb, which ever is greater
Alarms 警報	(1) Cell Failure (2) User Adjustable Moisture Level Alarms
Flow Control 流量	Recommended Inlet Flow: 1.1 L/min combined sample & bypass
Lower Detection Limit 最低偵測極限	35 ppbV
Environmental Condition 周界操作環境	0° to +60° Celsius +20% to 80% RH, non-condensing, non-corrosive atmosphere
Power 電源	100 - 240 VAC, 50/60 Hz
Sample Gas 樣品氣組成	AIR, ARGON, CARBON DIOXIDE, HELIUM, HYDROGEN, METHANE, NEON, NITROGEN, OXYGEN, XENON, OTHER
Sample Inlet Pressure 樣品氣壓力	10 – 3,000 PSIG (0.68 - 206.8 BarG)
Signal Outputs 訊號輸出	One isolated 0 – 5 V DC analog output or One isolated 4-20 mA output; user defined RS-232 serial communication port
Units of Measure 可換算單位	ppbV, ppmV, ppmW, dewpoint in °C or °F

MEECO AQUAVOLT+ 操作介面介紹


Figure 1

- | | |
|--------------------------|---|
| (A) Alphanumeric Display | Vacuum-fluorescent display, 2 characters high x 20 characters wide; displays moisture reading, units of measure, set-up selections and alarm status. |
| (B) Power Switch | Main Power on/off function |
| (C) "Menu" Key | Use to access the various set-up menus or to enter a selected option on the alphanumeric display |
| (D) Arrow Keys | Left and right arrow keys are used to display menu parameters for review or editing. The up and down arrow keys enable the operator to adjust the value of the parameter displayed. |

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Rear Panel


Figure 2a (also see figure 2B and 2C page 13)

- | | |
|------------------------|--|
| (A) Sample Inlet | 1/4" VCR male bulkhead fitting; 10-3000 psig (0.7-204 Bar) inlet pressure |
| (B) Alarms | Connector for alarm outputs |
| (C) Bypass Outlet | 1/8" male compression bulkhead fitting |
| (D) Sample Outlet | 1/8" male compression bulkhead fitting |
| (E) Cell Access Panel | For removal and replacement of electrolytic cell; also for access to internal pressure regulator for pressure adjustment |
| (F) Power Entry Module | Connector for main power |
| (G) Fuses (2) | Holds two 2.5A main power fuses |
| (H) RS-232 Connector | Connector for serial interface |
| (J) I/O Connector | Connector for 0-5 VDC / 4-20 mA / 0-20mA / 0-24 mA output |

MEECO AQUAVOLT+ 內部元件介紹

Internal View


Figure 3

Internal View (components)

(A) Electrolytic Cell	Senses moisture
(B) Sample Inlet Pressure Regulator	Do not allow inlet pressures to exceed 3000 psig (204 Bar) or go lower than 10 psig (0.7 Bar)
(C) Sample Mass Flow Controller	0–200 sccm control range (100 sccm standard for N ₂)
(D) Vacuum Fluorescent Display	Alphanumeric display; 2 characters high x 20 characters wide
(E) Power Supply	±15 VDC and +5 VDC; power supply for the mass flow controller, display and main board.
(F) Main Board	Single controller board for all logic, A/D, relay and communication functions
(G) 67.5V Battery	Provides power to cell during shipment or power loss

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
MEECO PRESSURE MEASUREMENT KIT

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- (A) Pressure gauge, 0-10 psig readout
- (B) One female fitting (reduction fitting)
- (C) One rotameter, 0-3 scfh readout
- (D) One male run tee
- (E) 1/8" OD stainless steel or PTFE tube (approximately 24")
- (F) Cap nut and ferrule

TOLERANCES SHOWN IN NOTES			
MEECO			
DESIGN	SCALE	DRAWN BY	EMP
#	1:2	APPROVED BY	
FUNCTIONAL	TYPE	EXTERNAL PRESSURE MEASUREMENT KIT	
#			
ANALOG	DATE	DRAWING NUMBER	
	3-19-90	103029	

MEECO AQUAVOLT+ STATUS DISPLAY

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Status Display


Figure 4

(A) Moisture Reading

This reading represents the measured moisture content of the sample gas.

(B) Mode Indicator

S - Service
M - Measure

(C) Flow Rate Indication

This area of the display will show the sample flow rate in sccm.

(D) Alarm State

1 - User Set
2 - User Set
C - Cell Alarm
F - Flow Alarm

MEECO AQUAVOLT+ ADJUSTING PARAMETERS

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Display Text 參數	Range of Value 選項	Definition or Usage 功能
Operating Mode: 操作模式	MEASURE, SERVICE SHUTDOWN	Measure Mode measures moisture under default parameters, whereas Service Mode parameters may be altered.
Moisture Units 濕度單位	PPMv, PPBv, DewPoint ° C, # / MMSCF DewPoint ° F,	
Sample Gas Type: 樣品氣種類	AIR, ARGON, CARBON DIOXIDE, HELIUM, HYDROGEN, METHANE, NEON, NITROGEN, OXYGEN, XENON, OTHER	The type of gas being analyzed. This is required so that the mass flow controller can adapt for differing molecular weights.
K-fac for User Gas: 樣品氣種類代數	0.50 to 3.00	K factor, for mass flow controllers, when using a gas not listed in the Gas Type parameter. This menu is shown only when gas type = OTHER. (see Appendix 3)
Communications ID#: 通訊代號	0 to 15	Aquavolts adress for serial communications.

MEECO AQUAVOLT+ ADJUSTING PARAMETERS

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Recorder Output MODE: 類比輸出訊號模式	0-5 Vdc, 4-20 mA, 0-20 mA, 0-24 mA	The selection of recorder output signal.
Recorder Scale: 類比輸出範圍	1 - any value, up to the upper limit of the unit. (.02ppm increments)	The range of values to provide full scale output of the recorder output..
Level Alarm Type: 濃度警報驅動RELAY輸出功能	1 or 2 (ENABLED or DISABLED)	The effect of the derived moisture level on the alarm output relay. Alarms will always be displayed in the Status Display regardless of this setting.
Level Alarm Value: 濃度警報設定值	Upper limit of unit (user set) to lower detection limit of unit.	The value at which point the derived moisture level is considered to be an alarm.
Cell Alarm: 水分感測器警報	Always Enabled	Alarm will trigger when the reading is below the lower detection limit.
Flow Alarm Type: 流量警報驅動RELAY輸出功能	ENABLED, DISABLED	The effect of flow deviation on the alarm output relay. Alarms will always be displayed in the Status Display regardless of this setting.
Flow Alarm Deviation: 流量誤差容許	3 to 25 (% of scale deviation)	The amount of deviation from set point at which to be considered a flow alarm.
Software Revision: 軟體版本	Read only, not adjustable	Software version and release date of the installed EPROM.
Uptime: 距離上次開機已運轉時間	Read only, not adjustable	Duration that the unit has been powered up.

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Alarm Relay Connector

1. 0-5Vdc Signal Out
2. 0-5 Vdc Common
3. 4-20 / 0-20 / 0-24 mA Common
4. 4-20 / 0-20 / 0-24 mA Signal Output

nc - Normally Closed
com - Common
no - Normally Open

5. com - Low Level Alarm
4. nc - Low Level Alarm
9. no - Low Level Alarm
1. com - Level Alarm1
2. nc - Level Alarm 1
6. no - Level Alarm1
8. com - Level Alarm 2
3. nc - Level Alarm 2
7. no - Level Alarm 2

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



MEECO AQUAVOLT+ , SPARE PARTS 備品清單

FUSE, 2.5 A	P/N: L3380	2 UNITS
AILR (Inert gases only) cell with battery (shelf life-3 months)	P/N: T7000	1 UNIT
AOLR (Oxygen Cell with battery (shelf life-3 months)	P/N: T7001	1 UNIT
AHLR (Hydrogen Cell with battery (shelf life-3 months)	P/N: T7002	1 UNIT
Battery	P/N: L0027	1 UNIT
+5/ +/- 15 Volt Power Supply	P/N: L3240	1 UNIT
PCB Assembly	P/N: L3795	1 UNIT

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







G.E. PANAMETRICS MMS35 露點控制器

G.E. PANAMETRICS MMS35 控制器

Functions: 功能	Dew Point 顯示露點 PPMv in gases at constant pressure 顯示 PPM水份 (pressure by programmable constant)
Inputs: 使用 M 系列 sensor	Single input via M-Series probe connected to a terminal strip. Probe may be remotely located up to 600 m (2,000 ft) from electronic console. Auxiliary: 4 to 20 mA input.
Recorder Output: 類比輸出	0 to 20-mA, 4 to 20-mA or 0 to 2-V analog, linear in parameter chosen.
Alarm Relays: 警報繼電器	Standard: 8A @ 250 VAC, 8A @ 30 VDC Hermetically Sealed: 0.3A @ 115 VAC, 2A @ 28 VDC Optional Form: C
Display Units: 顯示單位	Dew/frost point temperature: °C, °F PPMv Pressure: psi(g), bar, kPa(g), kg/cm2(g) MH: raw signal.
Power Requirements: 電源	AC power supply: Specify as 100, 120, 220, or 240 VAC, 50/60 Hz Optional DC power supply: 24 VDC Input power: 12 watts, maximum.
Temperature: 周界操作溫度	Operating: 0° to 60°C (32° to 140°F) Storage: -20° to 70°C (-22° to 158°F)

G.E. PANAMETRICS MMS35 電器接端

G.E. PANAMETRICS MMS35 電器接端

Table 2-5: High & Low Alarm Connections

Connect Alarm A:	To ALARM A Terminal Block:
NC Contact	pin #4
NO Contact	pin #5
A Contact	pin #6
Connect Alarm B:	To ALARM B Terminal Block:
NC Contact	pin #7
NO Contact	pin #8
A Contact	pin #9

Table 2-8: Recorder Output Connections

Connect:	To REC OUT Terminal Block:
Out (+)	pin #10
Return (-)	pin #11

Table 2-9: Auxiliary Input Connections

Connect:	To IN Terminal Block:
Out (+)	pin #7
Return (-)	pin #8
Power (+V)	pin #9

Table 2-6: Fault Alarm Connections

Connect Fault Alarm:	To FAULT ALARM Terminal Block:
NC Contact	pin #1
NO Contact	pin #2
A Contact	pin #3




Figure 2-7: S1 Location on 703-1175 Output Board

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CONNECTING SENSOR PROBE TO MMS35 感測器接線

Connect:	To PROBE Terminal Block:
Red (H2) wire	pin #1
Shield	pin #2
Green (H1) wire	pin #3

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

GE PANAMETRICS MMS TROUBLESHOOTING 故障排除

症狀	原因	系統反應	解決方案
SENSOR 準確度 有疑慮	等待分析儀穩定時間不夠	SENSOR 在乾燥製成露點讀值過高， 在濕製成露點讀值過低	確認系統purge時間足夠，改變流量後若 露點改變代表管線有漏。
	取樣點露點與 製成露點不同	SENSOR 露點讀值過高或過低	Readings may be correct if the sampling point and main stream do not run under the same process conditions. The different process conditions cause readings to vary. Refer to Appendix A for more information. If sampling point and main stream conditions are the same, check sample system pipes, and any pipe between the sample system and main stream for leaks. Also, check sample system for adsorbing water surfaces, such as rubber or plastic tubing, paper-type filters, or condensed water traps. Remove or replace contaminating parts with stainless steel parts.
	SENSOR 受汙染	Probe reads too wet or too dry.	Clean the sensor and the sensor shield as described in Appendix A. Then reinstall the sensor.
	SENSOR 受汙染	Probe reads high dew point.	Clean the sensor and the sensor shield as described in Appendix A, then reinstall the sensor. Also, install a proper filter (i.e. sintered or coalescing element).

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


NOTE 附註:

- Used standards:
A. Analyzer: Chilled Mirror Hygrometer, model #: RH 373LX, S/N: 12-0748
B. Analyzer calibrated by: Thunder Scientific Corporation, report #: 10552, NIST certified
- 95% confidence interval with uncertainties considering intrinsic errors with this method:
DP range -80°C to -70°C: $\pm 0.6^\circ\text{C}$ DP range -70°C to 10°C: $\pm 0.2^\circ\text{C}$
- NIST traceable operating range: -80°C to +20°C, analyzer displays from: -110°C to +20°C
- Ambient temperature at calibration: 24.0°C, ambient humidity at calibration: 60.0%

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


COSA XENTAUR LPDT 操作模式

Appendix A: Operating State User Interface flowchart


LEGEND: 1. Arrows leading out from below a button depict the flow as a result of pressing the button.
 2. Arrows leading out of other locations depict flow as a result of completing a function or a time-out.
 3. A slash '/' in the display area, is used to depict two alternately shown (flashing back and forth) messages.

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
MICHELL EASIDEW 露點計規範

Transmitter

Type	Michell Easidew
Transmitter Torque Loading	Minimum 30.5 Nm (22.5 lbf-ft)
Calibration Range	-100 to +20°C (-148 to +68°F) dew point
Output	4-20 mA current source configurable over the entire dew-point range
Interchangeability	Fully interchangeable
Dew-Point Accuracy	± 2°C (± 3.6°F)
Operating Temperature	-40 to +60°C (-40 to +140°F)
Storage Temperature	-40 to +75°C (-40 to +167°F)
Temperature Coefficients	Temperature compensated across operating temperature range
Operating Pressure	10-6 bara vacuum to 45 MPa (450 barg / 6500 PSI)
Flow Rate	1 to 5 l/min (2.1 to 10.6 SCFH) mounted in standard sampling block 0 to 10 m/sec (0 to 353 SCFS) direct insertion (80 µm sintered guard)
Traceable Certification	-90 to +20°C (-130 to +68°F) dew point traceable to the National Physical Laboratory; -75 to +20°C (-103 to +68°F) dew point traceable to NIST (USA) [For dew points < -90°C (< -130°F): Direct reference to a fundamental cooled mirror dew-point meter]
Environmental Protection	IP66 and NEMA 4
Weight	150g (5.29oz)

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MICHELL CERMAX 露點計

MICHELL CERMAX 露點計規範

Sensor	Kahn Ceramic Moisture Sensor
Gas wetted components	316 stainless steel Components
Gas connections	1/8" Swagelok inlet and outlet fittings.
Display	Standard version: 240 x 64 dot matrix LCD
Range	Calibrated from -100 to +20 oC dewpoint; readings to -120 and +30 oC dewpoint
Accuracy	+/- 1 oC from -60 to +20 oC dewpoint; +/- 2 oC from -100 to -60 oC dewpoint
Resolution	0.1 oC dewpoint, 3 significant figures for other units
Units	oC, oF, K dewpoint, PPM(V), PPM(W) for air, N2, H2, CO2, natural gas, %RH, gm-3 (natural gas), #/MMSCF (natural gas)
Secondary input	0-20 mA or 4-20 mA for temperature or pressure
Data storage	Up to 10,000 samples of primary and secondary variable, time and date stamp and identification tag
Communications	RS232C for stored data
Power	Internal re-chargeable battery pack (removable) charged by external a.c. powered 6V charger (supplied). 24 hours normal operation between charges. Charge time is up to 12 hours but can be up to 14 hours if the instrument is in use at the same time.
Case	Custom polyurethane case with integral padded carry handle
Dimensions	250W x 300D x 150H mm approx.
Weight	8.8 lbs (3 kg)
Ingress Protection	IP65
Operating temperature	-20 to +50 OC

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