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* CE, collision energy; CXP, collision cell exit potential.

** transitions used for quantification.

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- chromatography-tandem mass spectrometry using $^{13}\text{C}15$ -deoxynivalenol as internal standard. *Food Addit. Contam. Part A* 27: 1701-1713.
8. Desmarchelier, A. and Seefelder, W. 2011. Survey of deoxynivalenol and deoxynivalenol-3-glucoside in cereal-based products by liquid chromatography electrospray ionization tandem mass spectrometry. *World Mycotoxin J.* 4: 29-35.
 9. Vendl, O., Crews, C., MacDonald, S., Krska, R. and Berthiller, F. 2010. Occurrence of free and conjugated Fusarium mycotoxins in cereal based food. *Food Addit. Contam. Part A* 27: 1148-1152.
 10. Beyer, M., Danicke, S., Rohweder, D. and Humpf, H. U. 2010. Determination of deoxynivalenol-sulfonate (DONS) in cereals by hydrophilic interaction chromatography coupled to tandem mass spectrometry. *Mycotoxin Res.* 26: 109-117.
 11. Iqbal, S., Paterson, R., Bhatti, I. and Asi, M. 2010. Survey of aflatoxins in chilies from Pakistan produced in rural, semi-rural and urban environments. *Food Addit. Contam. Part B* 3: 268-274.
 12. Ahn, J., Kim, D., Kim, H. and Jahng, K. Y. 2010. Quantitative determination of mycotoxins in urine by LC-MS/MS. *Food Addit. Contam. Part A* 27: 1674-1682.
 13. Han, Z., Zheng, Y., Luan, L., Ren, Y. and Wu, Y. 2010. Analysis of ochratoxin A and ochratoxin B in traditional Chinese medicines by ultra-high-performance liquid chromatography-tandem mass spectrometry using [$^{13}\text{C}20$]-ochratoxin A as an internal standard. *J. Chromatogr. A* 1217: 4365-4374.
 14. Yue, Y., Zhang, X., Yang, M., Ou-Yang, Z. and Liu, H. 2010. Simultaneous determination of deoxynivalenol and nivalenol in traditional Chinese medicine by SPE and LC. *Chromatographia*. 72: 551-555.
 15. Monbaliu, S., Poucke, C. V., Detavernier, C., Dumoulin, F., Mario, V. D. V., Schoeters, E., Stefaan, V. D., Averkieva, O., van Peteghem, C. and de Saeger, S. 2010. Occurrence of mycotoxins in feed as analyzed by a multi-mycotoxin LC-MS/MS method. *J. Agric. Food Chem.* 58(1): 66-71.
 16. Beltrán, E., Ibáñez, M., Sancho, J. V., Cortés, M. Á., Yusà, V. and Hernández, F. 2011. UHPLC-MS/MS highly sensitive determination of aflatoxins, the aflatoxin metabolite M1 and ochratoxin A in baby food and milk. *Food Chem.* 126(2): 734-744.
 17. Ediage, E. N., Di Mavungu, J. D., Monbaliu, S., van Peteghem, C. and de Saeger, S. 2011. A validated multianalyte LC-MS/MS method for quantification of 25 mycotoxins in cassava flour, peanut cake and maize samples. *J. Agric. Food Chem.* 59: 5173-5180.
 18. Tanaka, K., Sago, Y., Zheng, Y., Nakagawa, H. and Kushiro, M. 2007. Mycotoxin in rice. *Int. J. Food Microbiol.* 119: 59-66.
 19. Tutelyan, V. A. 2004. Deoxynivalenol in cereal in Russia. *Toxicol. Lett.* 53: 173-179.

A Survey of Mycotoxins in Cereal Products

YING-CHUN CHEN, WEI-GUANG WU, CHIA-DING LIAO, HSU-YANG LIN,
LIH-CHING CHIUEH AND DANIEL YANG-CHIH SHIH

Division of Research and Analysis, FDA

ABSTRACT

Mycotoxins are secondary metabolites produced by *Aspergillus*, *Fusarium*, *Penicillium* and *Alternaria* that are known as carcinogenic and harmful to the health of humans and animals, and even caused huge economic loss. Some molds are capable of producing more than one mycotoxin and some mycotoxins are produced by several fungal species. Often many mycotoxins are found on a contaminated substrate. Factors influencing the presence of mycotoxins in foods include temperature and environmental conditions. A multi-mycotoxin analytical method was employed for simultaneous determination of aflatoxin B₁ (AFB₁), aflatoxin B₂ (AFB₂), aflatoxin G₁ (AFG₁), aflatoxin G₂ (AFG₂), fumonisin B₁ (FB₁), fumonisin B₂ (FB₂), deoxynivalenol (DON), zearalenone (ZON), ochratoxin A (OTA), T-2 and HT-2 toxin in cereals. The regulations have not been set for the FB₁, FB₂, DON, ZON, T-2 and HT-2 in Taiwan. In order to understand the contamination situation in the food, a survey of multi-mycotoxin contents in market products was conducted by TFDA. One hundred products including rice, corn, nut, wheat and peanut products were collected from convenient stores, supermarkets and grocery stores from February to March, 2012. The results showed that mycotoxins were detected in 15 samples ranging from 0.57 ppb (AFB₂) to 281 ppb (DON), in 1 rice, 3 corn, 7 wheat, and 4 peanut products. The result also showed that none of sample was contaminated with mycotoxin levels above the Taiwan and international regulation limits. Except for total aflatoxin and ochratoxin A, there is no official regulation in Taiwan. This survey result may offer government a reference to set up the regulation limits for the sake of national health and safety.

Key words: mycotoxin, cereals