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Arsenic compound	Intercept	Slope	R ²
As ³⁺	19809	49450	0.9995
As ⁵⁺	28808	55454	0.9997
MMA	3701.9	68406	1.0000
DMA	6106.9	69820	1.0000
AsB	12028	73322	0.9999
AsC	28287	76974	0.9997

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
Inorganic arsenic	Spiked level (mg/kg)	S/N ratio ^a
As ³⁺	0.01	4.97 ± 0.95
	0.02	10.1 ± 0.85
	0.05	21.2 ± 3.68
As ⁵⁺	0.01	7.67 ± 1.01
	0.02	14.2 ± 1.23
	0.05	36.8 ± 5.86

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Inorganic arsenic	Intra-day precision (n=5)		
	Spiked level (mg/kg)	Recovery ^a (%)	CV (%)
As ³⁺	0.02	98.0 \pm 2.87	2.93
	0.1	89.5 \pm 3.48	3.89
	0.2	92.8 \pm 5.60	6.04
As ⁵⁺	0.02	103.3 \pm 1.65	1.59
	0.1	91.5 \pm 3.67	4.01
	0.2	93.0 \pm 5.04	5.42

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1. International Programme on Chemical Safety. 2001. Arsenic and arsenic compounds. Environmental Health Criteria 224. [<http://www.inchem.org/documents/ehc/ehc/ehc224.htm>].
2. Sahoo, P. K. and Kim, K. 2013. A review of the arsenic concentration in paddy rice from the perspective of geoscience. Geosci. J. 17(1): 107-122.
3. Cullen, W. R. and Reimer, K. J. 1989. Arsenic speciation in the environment. Chem. Rev. 89: 713-764.
4. Francesconi, K. A. and Edmonds, K. J. 2013. Arsenic and marine organisms. Adv. Inorg. Chem. 44: 147-189.
5. Siscovick, D. S., Barringer, T. A., Fretts, A. M. and *et al.* 2017. Omega-3 polyunsaturated fatty acid (fish oil) supplementation and the prevention of clinical cardiovascular disease: a science advisory from the American heart association. Circulation 135(15): e867-e884.
6. 謝俊明。1999。砷生物偵測物種分析研究。行政院勞工委員會勞工安全衛生研究所，台北。
7. International Agency for Research on Cancer. 2020. Agents Classified by the IARC Monographs, Volumes 1-127. [<http://monographs.iarc.fr/agents-classified-by-the-iarc/>].
8. Hata, A., Yamanaka, K., Endo, G. and *et al.* 2012. Arsenic metabolites in humans after ingestion of wakame seaweed. 16th International Conference on Heavy Metals in the Environment (ICHMET 2012), Angelicum Conference Centre, Rome, Italy.
9. Crebelli, R. and Leopardi, P. 2012. Long-term risks of metal contaminants in drinking water: a critical appraisal of guideline values for arsenic and vanadium. Ann. Ist. Super. Sanita. 48(4): 354-361.
10. Ma, L., Yang, Z., Kong, Q. and *et al.* 2017. Extraction and determination of arsenic compound in leafy vegetables: method development and application. Food Chem. 217: 524-530.
11. Yuan, C. G., Jiang, G. B. and He, B. 2005. Evaluation of the extraction methods for arsenic speciation in rice straw, *Oryza sativa* L., and analysis by HPLC-HG-AFS. J. Anal. At. Spectrom. 20: 103-110.
12. Baranowska, I. 2015. Fundamentals and applications cavitation. Handbook of Trace Analysis. pp. 136-137. Springer International Publishing, Switzerland.
13. Bruno, E. S., Luciana, M., Cleide, S. T. and *et al.* 2016. Analytical strategies for the determination of arsenic in rice. J. Chem. 2016: 1-11.
14. Romina, P., Juan, A., Silvia, S. and *et al.* 2011. Organic solvent-free reversed-phase ion-pairing liquid chromatography coupled to atomic fluorescence spectrometry for organoarsenic compound determination in

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Development of the Analytical Method for Inorganic Arsenic in Aquatic Animal Oil

YU-CHENG LAI, YU-NING SHIH, YING-RU SHEN, NU-CHING LIN,
YA-MIN KAO, SU-HSIANG TSENG AND DER-YUAN WANG

Division of Research and Analysis, TFDA

ABSTRACT

Arsenic can be classified into organic and inorganic compounds. Inorganic arsenic is much more toxic than the organic ones. Aquatic animals may accumulate arsenic due to the contaminated environment. To safeguard the public health, TFDA actively developed the testing method on the speciation of arsenic in fish oil by high performance liquid chromatography coupled with inductively coupled plasma mass spectrometry (HPLC-ICP/MS). Sample was homogenized and extracted with 0.28 M HNO₃ in an ultrasonic bath at 80°C for 30 min twice. Analytes such as inorganic arsenic (arsenite, arsenate) and organic arsenic (DMA, MMA, AsB, AsC) in the sample solution were separated effectively on an anion exchange column, and detected by ICP-MS in less than 15 min. The calibration curves observed coefficients of determination higher than 0.9995. The limits of quantification for arsenite and arsenate in fish oil were 0.02 mg/kg. The average recoveries of arsenite and arsenate spiked into fish oil were in the range of 89.5-98.0% and 91.5-103.3%, respectively. Surveillance studies contained 16 fish oil samples purchased from markets were conducted. The results showed that inorganic arsenic was not found in all samples. The results also complied with the maximum limits of inorganic arsenic in fish oil listed in "Sanitation Standards for Contaminants and Toxins in Food" preannounced by the Ministry of Health and Welfare, Taiwan. The testing method offered fast analysis and simple sample pretreatment, and can be applied into routine testing for the determination of inorganic arsenic in fish oil products for the purpose of public health protection.

Key words: fish oil, inorganic arsenic, HPLC-ICP/MS.