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Interference substances	Concentration
Endogenous substances	
Mucin	100 µg/mL
Biotin	100 µg/mL
Exogenous substances - Relevant medicines	
Lumefantrine (Malaria)	50 µM
Artemether (Malaria)	50 µM
Zanamivir (Influenza)	5 mg/mL
Oseltamivir (Influenza)	10 mg/mL
Chlorpheniramine (Influenza)	0.08 mg/dL
Ivermectin (Scabies)	4.4 µg/mL
Exogenous substances - Antiviral and antibacterial	
Amoxicillin	5.4 mg/dL
Ritonavir (HCV)	16.4 µg/mL
Exogenous substances - Anti-inflammatory	
Ibuprofen	21.9 mg/dL
Acetaminophen	200 µM
Acetylsalicylic acid	3 mg/dL
Triamcinolone (Stomatitis)	1.18 ng/mL
Exogenous substances - Homeopathic allergy relief medicine	
Olopatadine Hydrochloride	10 mg/mL
Cromolyn Sodium salt	20 mg/mL
Exogenous substances - Nasal sprays	
Oxymetazoline hydrochloride	10% v/v
Exogenous substances - Anti-hypertensive	
Chlorothiazide	2.7 mg/dL
Indapamide	140 ng/mL
Exogenous substances - Anti-diabetic	
Glimepiride	0.164 mg/dL
Acarbose	0.03 mg/dL
Exogenous substances - Other	
Benzocaine (Oral anesthetic)	0.13 mg/mL
Budesonide (Bronchial Asthma)	2.76 ng/mL
Guaiacol glyceryl ether (Common cold)	1 µg/mL
Diphenhydramine (Allergic rhinitis)	0.08mg/dL
Sulfur	9.23 µg/mL
Methanol	0.15 mg/mL

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TTREND COVID-19 Antigen Test	Real-time PCR method: LabTurbo AIO COVID-19 RNA Testing Kit		
	+	-	Total
+	7	0	7
-	0	52	52
Total	7	52	59

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- Lu R, Zhao X, Li J, *et al.* Genomic characterization and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. Lancet 2020; 395:565.
- Zhu N, Zhang D, Wang W, *et al.* A Novel Coronavirus from Patients with Pneumonia in China, 2019. N Engl J Med 2020; 382(8): 727–733.
- Li Q, Guan X, Wu P, *et al.* Early transmission dynamics in Wuhan, China, of novel coronavirus-infected pneumonia. N Engl J Med 2020; 382:1199.

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