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

The diagram illustrates the process of DNA replication in three main stages:

- Initiation:** The original DNA double helix (green strands) is unwound. A DNA primer (red) is synthesized at the 3' end of one template strand. Nucleotides (blue) are added to the 5' end of the primer.
- Elongation:** The DNA polymerase (orange) moves along the template strand, synthesizing a new strand (blue) in the 5' to 3' direction. The template strand is labeled with 5' and 3' ends. The new strand is labeled with 3' and 5' ends. The process is labeled with a blue circle containing the number 1.
- Termination:** The DNA polymerase reaches the end of the template strand, and the new strand is released. The process is labeled with a blue circle containing the number 2.


The final result is two identical DNA double helices, each consisting of one original template strand (green) and one newly synthesized strand (blue). The process is labeled with a blue circle containing the number 3.

Labels in the diagram include: original DNA to be replicated, DNA primer, nucleotide, 5', 3', R, F, double stranded DNA, Hydrogen bond, and a chemical structure of a nucleotide base (G) with its sugar and phosphate groups.

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



BioEdit

UTR-HAV68 (F)

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
BioEdit Sequence Alignment Editor

File Edit view Zoom Horizontal Scale Accessory Application RNA Window Help


ABI Chromatogram: C:\Users\shingyuan\AppData\Local\Temp\15\Temp529781e-b201-498b-aac7-89c082aee5b0_SD230506004_230508K08D.zip\UTR-560_HAV68.ab1

Selected: none Sample: UTR-560_HAV68 File: C:\Users\shingyuan\AppData\Local\Temp\15\Temp529781e-b201-498b-aac7-89c082aee5b0_SD230506004_230508K08D.zip\UTR-560_HAV68.ab1

ABI Chromatogram: C:\Users\shingyuan\AppData\Local\Temp\15\Tempf9873d56-0ca6-4d28-af06-2e572253ef27_SD230506004_230508K08D.zip\UTR-560_HAV240.ab1

Selected: none Sample: UTR-560_HAV240 File: C:\Users\shingyuan\AppData\Local\Temp\15\Tempf9873d56-0ca6-4d28-af06-2e572253ef27_SD230506004_230508K08D.zip\UTR-560_HAV240.ab1

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
 An official website of the United States government [Here's how you know](#) ✓

NIH National Library of Medicine
National Center for Biotechnology Information

Log in

All Databases Search

NCBI Home
Resource List (A-Z)
All Resources
Chemicals & Bioassays
Data & Software
DNA & RNA
Domains & Structures
Genes & Expression
Genetics & Medicine
Genomes & Maps
Homology
Literature
Proteins
Sequence Analysis
Taxonomy
Training & Tutorials
Variation

Welcome to NCBI
The National Center for Biotechnology Information advances science and health by providing access to biomedical and genomic information.
[About the NCBI](#) | [Mission](#) | [Organization](#) | [NCBI News & Blog](#)

Submit
Deposit data or manuscripts into NCBI databases

Download
Transfer NCBI data to your computer

Learn
Find help documents, attend a class or watch a tutorial

Develop
Use NCBI APIs and code libraries to build applications

Analyze
Identify an NCBI tool for your data analysis task

Research
Explore NCBI research and collaborative projects

Popular Resources
PubMed
Bookshelf
PubMed Central
BLAST
Nucleotide
Genome
SNP
Gene
Protein
PubChem

NCBI News & Blog
Comparing Yeast Species Used in Beer Brewing and Bread Making 29 Sep 2023
Using the NIH Comparative Genomics Resource (CGR) to gain knowledge
Join NCBI at ASHG 2023 27 Sep 2023
November 1-5 in Washington, D.C. We look forward to seeing you in person at the American Society for Human

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NCBI BLAST

National Library of Medicine
National Center for Biotechnology Information

Log in

BLAST®

[Home](#) [Recent Results](#) [Saved Strategies](#) [Help](#)

Check out the ClusteredNR database on BLAST+

[Learn more](#)

[Give us feedback](#)

Basic Local Alignment Search Tool

BLAST finds regions of similarity between biological sequences. The program compares nucleotide or protein sequences to sequence databases and calculates the statistical significance. [Learn more](#)

NEWS

ClusteredNR database on BLAST+

The ClusteredNR database is now available for BLAST+

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
Standard Nucleotide BLAST

blastn blastp blastx tblastn tblastx

BLASTN programs search nucleotide databases using a nucleotide query

Enter Query Sequence

Enter accession number(s), gi(s), or FASTA sequence(s) ? Clear

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Job Title

Enter a descriptive title for your BLAST search ?

☐ Align two or more sequences ?

Choose Search Set

Database ☒ Standard databases (nr etc.): ☐ rRNA/ITS databases ☐ Genomic + transcript databases ☐ Betacoronavirus

New ☐ Experimental databases [Try experimental taxonomic nt databases](#)
For more info see [What are taxonomic nt databases?](#)

Nucleotide collection (nr/nt) ?

Organism Optional ☐ exclude

Enter organism common name, binomial, or tax id. Only 20 top taxa will be shown ?

Exclude Optional ☐ Models (XM/XP) ☐ Uncultured/environmental sample sequences

Limit to Optional ☐ Sequences from type material

Entrez Query Optional [YouTube](#) [Create custom database](#)

Program Selection

Optimize for ☒ Highly similar sequences (megablast)
☐ More dissimilar sequences (discontiguous megablast)
☐ Somewhat similar sequences (blastn)
Choose a BLAST algorithm ?

BLAST Search **database Nucleotide collection (nr/nt)** using **Megablast (Optimize for highly similar sequences)**
☐ Show results in a new window

UTR-HAV68 (序列長度144 bp)

Job Title	Nucleotide Sequence
RID	HR4A59AK013 Search expires on 10-04 11:57 am Download All ▾
Program	BLASTN ? Citation ▾
Database	nt See details ▾
Query ID	lcl Query_30177
Description	None
Molecule type	dna
Query Length	144
Other reports	Distance tree of results MSA viewer ?

Filter Results

Organism only top 20 will appear ☐ exclude

Type common name, binomial, taxid or group name

[+ Add organism](#)

Percent Identity

E value

Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len
228	228	87%	3e-55	99.21%	675
222	222	87%	1e-53	98.43%	7477

Descriptions

Graphic Summary

Alignments

Taxonomy

Sequences producing significant alignments

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UTR-HAV68 (序列長度125bp)

Job Title	Nucleotide Sequence
RID	HRFEXZWZ01R Search expires on 10-04 15:08 pm Download All ▾
Program	BLASTN ? Citation ▾
Database	nt See details ▾
Query ID	lcl Query_185205
Description	None
Molecule type	dna
Query Length	125
Other reports	Distance tree of results MSA viewer ?

Filter Results

Organism *only top 20 will appear* ☐ exclude

Type common name, binomial, taxid or group name

[+ Add organism](#)

Percent Identity to **E value** to **Query Coverage** to

[Filter](#) [Reset](#)

Descriptions

[Graphic Summary](#)[Alignments](#)[Taxonomy](#)

Sequences producing significant alignments

[Download](#) ▾[Select columns](#) ▾Show [?](#)☒ select all 100 sequences selected[GenBank](#)[Graphics](#)[Distance tree of results](#)[MSA Viewer](#)

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
☒	Hepatitis A virus isolate G1B3-IRES 5' UTR	Hepatovirus A	226	226	97%	8e-55	100.00%	675	FJ829475.1
☒	Hepatitis A virus isolate G1B1-IRES 5' UTR	Hepatovirus A	226	226	97%	8e-55	100.00%	675	FJ829473.1
☒	Hepatovirus A strain HepA/USA_NY/51657/2019_partial genome	Hepatovirus A	220	220	97%	4e-53	99.18%	7361	ON524431.1
☒	Hepatitis A virus isolate A3_complete genome	Hepatovirus A	220	220	97%	4e-53	99.18%	7339	KC182589.1
☒	Hepatitis A virus isolate A2_complete genome	Hepatovirus A	220	220	97%	4e-53	99.18%	7339	KC182587.1
☒	Hepatitis A virus isolate HAV-Arg/06_complete genome	Hepatovirus A	220	220	97%	4e-53	99.18%	7477	HM769724.1
☒	Hepatitis A virus isolate G1B4-IRES 5' UTR	Hepatovirus A	220	220	97%	4e-53	99.18%	675	FJ829476.1
☒	Hepatitis A virus isolate G1B2-IRES 5' UTR	Hepatovirus A	220	220	97%	4e-53	99.18%	675	FJ829474.1
☒	Hepatitis A virus gene_5' UTR_partial sequence_isolate: HA-JNG02-91	Hepatovirus A	220	220	97%	4e-53	99.18%	395	AB245895.1

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序列比對Alignments

Hepatitis A virus isolate G1B3-IRES 5' UTR

Sequence ID: [FJ829475.1](#) Length: 675 Number of Matches: 1

Range 1: 58 to 179 [GenBank](#) [Graphics](#)

▼ Next Match

Score	Expect	Identities	Gaps	Strand
226 bits(122)	8e-55	122/122(100%)	0/122(0%)	Plus/Plus
Query 3	CCTATCTCCCTTTGTTTTGTTTGTAATATTAATTCCTGCAGGTTCAAGGTTCTTTAATC	62		
Sbjct 58	CCTATCTCCCTTTGTTTTGTTTGTAATATTAATTCCTGCAGGTTCAAGGTTCTTTAATC	117		
Query 63	TGTTTCTCTATAAGAACTCAATTTTCACGCTTCTGTCTTCTTTCTTCCAGGGCTCTC	122		
Sbjct 118	TGTTTCTCTATAAGAACTCAATTTTCACGCTTCTGTCTTCTTTCTTCCAGGGCTCTC	177		
Query 123	CC 124			
Sbjct 178	CC 179			

Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
226	226	97%	8e-55	100.00%	675	FJ829475.1
220	220	97%	4e-53	99.18%	7477	HM769724.1

Query coverage= 122/125

Hepatitis A virus isolate HAV-Arg/06, complete genome

Sequence ID: [HM769724.1](#) Length: 7477 Number of Matches: 1

Range 1: 120 to 241 [GenBank](#) [Graphics](#)


▼ Next Match

Score	Expect	Identities	Gaps	Strand
220 bits(119)	4e-53	121/122(99%)	0/122(0%)	Plus/Plus
Query 3	CCTATCTCCCTTTGTTTTGTTTGTAATATTAATTCCTGCAGGTTCAAGGTTCTTTAATC	62		
Sbjct 120	CCTATTTCCCTTTGTTTTGTTTGTAATATTAATTCCTGCAGGTTCAAGGTTCTTTAATC	179		
Query 63	TGTTTCTCTATAAGAACTCAATTTTCACGCTTCTGTCTTCTTTCTTCCAGGGCTCTC	122		
Sbjct 180	TGTTTCTCTATAAGAACTCAATTTTCACGCTTCTGTCTTCTTTCTTCCAGGGCTCTC	239		
Query 123	CC 124			
Sbjct 240	CC 241			

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Graphic Summary and Taxonomy

UTR-HAV240 (序列長度151bp)

Job Title	Nucleotide Sequence		
RID	HRF0KNY3016	Search expires on 10-04 15:00 pm	Download All ▾
Program	BLASTN ?	Citation ▾	
Database	nt	See details ▾	
Query ID	cl Query_409437		
Description	None		
Molecule type	dna		
Query Length	151		
Other reports	Distance tree of results	MSA viewer	?

Filter Results

Organism only top 20 will appear ☐ exclude

Type common name, binomial, taxid or group name

[+ Add organism](#)

Percent Identity to **E value** to **Query Coverage** to

[Filter](#) [Reset](#)

Descriptions

[Graphic Summary](#)[Alignments](#)[Taxonomy](#)

Sequences producing significant alignments

[Download](#) ▾[Select columns](#) ▾Show [?](#)☒ select all 100 sequences selected[GenBank](#)[Graphics](#)[Distance tree of results](#)[MSA Viewer](#)

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
☒	Hepatitis A virus isolate HAV-Arg/06, complete genome	Hepatovirus A	230	230	86%	8e-56	98.47%	7477	HM769724.1
☒	Hepatitis A virus isolate G1B3-IRES 5' UTR	Hepatovirus A	230	230	86%	8e-56	98.47%	675	FJ829475.1
☒	Hepatitis A virus isolate G1B1-IRES 5' UTR	Hepatovirus A	230	230	86%	8e-56	98.47%	675	FJ829473.1
☒	Hepatovirus A strain HepA/USA_PA/42382/2019, partial genome	Hepatovirus A	224	224	86%	4e-54	97.71%	7361	ON524425.1
☒	Hepatitis A virus isolate G1B2-IRES 5' UTR	Hepatovirus A	224	224	86%	4e-54	97.71%	675	FJ829474.1
☒	Hepatitis A virus GBM/HFS RNA	Hepatovirus A	224	224	86%	4e-54	97.71%	7446	X75216.1
☒	Hepatitis A virus genome, 5'UTR	Hepatovirus A	224	224	86%	4e-54	97.71%	709	X70846.1
☒	Hepatovirus A genomic RNA, complete genome isolate: MNA06-2148	Hepatovirus A	219	219	86%	2e-52	96.95%	7492	LC049342.1

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
Query Coverage

 to [Filter](#)[Reset](#)

No significant similarity found. For reasons why, [click here](#)

ABI Chromatogram: C:\Users\shingyuan\AppData\Local\Temp\15\Temp7582c2ba-742b-40aa-beed-9e9b61b2ea66_SD230506004_230508K08D.zip\D-716-UTR_HAV68.ab1

Selected: none Sample: D-716-UTR_HAV68 File: C:\Users\shingyuan\AppData\Local\Temp\15\Temp7582c2ba-742b-40aa-beed-9e9b61b2ea66_SD230506004_230508K08D.zip\D-716-UTR_HAV68.ab1

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Job Title	Nucleotide Sequence
RID	HRNYMKA401R Search expires on 10-04 16:59 pm Download All ▼
Program	BLASTN ? Citation ▼
Database	nt See details ▼
Query ID	lc Query_372299
Description	None
Molecule type	dna
Query Length	137
Other reports	Distance tree of results MSA viewer ?

Filter Results

Organism only top 20 will appear ☐ exclude

Type common name, binomial, taxid or group name

+ [Add organism](#)

Percent Identity

E value

Query Coverage

to

to

to

Filter

Reset

Descriptions

Graphic Summary

Alignments

Taxonomy

Sequences producing significant alignments

Download ▼

Select columns ▼

Show 100 ▼

☒ select all 1 sequences selected

[GenBank](#)

[Graphics](#)

[Distance tree of results](#)

[MSA Viewer](#)

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
☒	PREDICTED: <i>Fragaria vesca</i> subsp. <i>vesca</i> acetylornithine aminotransferase, mitochondrial-like (LOC101295...	<i>Fragaria vesca</i> s...	126	126	54%	1e-24	97.30%	1695	XM_004294056.2

PREDICTED: *Fragaria vesca* subsp. *vesca* acetylornithine aminotransferase, mitochondrial-like (LOC101295312), mRNA

Sequence ID: [XM_004294056.2](#) Length: 1695 Number of Matches: 1

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The image shows the NCBI BLAST search interface. It has tabs for different BLAST programs: **blastn**, blastp, blastx, tblastn, and tblastx. The **blastn** tab is selected.

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Job Title

Enter a descriptive title for your BLAST search [?](#)

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Program Selection

Optimize for

☒ Highly similar sequences (megablast)

☐ More dissimilar sequences (discontiguous megablast)

☐ Somewhat similar sequences (blastn)

Choose a BLAST algorithm [?](#)

BLAST

Search **nucleotide sequence** using **Megablast (Optimize for highly similar sequences)**

☐ Show results in a new window

Annotations: A pink arrow points to the query sequence input field. A green arrow points to the 'Align two or more sequences' checkbox. A pink arrow points to the subject sequence input field. A pink arrow points to the 'BLAST' button.

序列比對Alignments

Job Title	Nucleotide Sequence
RID	HRWRBUJU114 Search expires on 10-04 18:54 pm Download All ▼
Program	Blast 2 sequences Citation ▼
Query ID	lcl Query_45117 (dna)
Query Descr	None
Query Length	342
Subject ID	lcl Query_45119 (dna)
Subject Descr	None
Subject Length	342
Other reports	MSA viewer ?

Filter Results

Percent Identity	E value	Query Coverage
to	to	to
Filter		Reset

Descriptions		Graphic Summary	Alignments	Dot Plot
Sequences producing significant alignments				
Download ▼ Select columns ▼ Show 100 ?				
☒ select all 1 sequences selected				
Description		Scientific Name		
☒ None provided				
Range 1: 21 to 342 Graphics		▼ Next Match ▲		
Score	Expect	Identities	Gaps	Strand
584 bits(316)	1e-171	321/323(99%)	1/323(0%)	Plus/Plus
Query 20	GGATTCTCCAGGTTTCAATTCAGGAAATGCTCAGGTACTTTCTTTGCTAAACTGGATC	79		
Sbjct 21	GGATTCTCCAGGTTTC-ATTCAGGAAATGCTCAGGTACTTTCTTTGCTAAACTGGATC	79		
Query 80	CTCAATTGTTGTGATAGCTCCACAGGTGCTTGACTCCTGAAACATCCATTTCCCTCT	139		
Sbjct 80	CTCAATTGTTGTGATAGCTCCACAGGTGCTTGACTCCTGAAACATCCATTTCCCTCT	139		
Query 140	GTTAGCTTTTCCTTTCAAATCTTTTCATGGTGTGTATACCAACTTGGGGATCTGGAACATT	199		
Sbjct 140	GTTAGCTTTTCCTTTCAAATCTTTTCATGGTGTGTATACCAACTTGGGGATCTGGAACATT	199		
Query 200	CTGTTCTGTAGAAACTGTTGTTGAAAAACCTCCAGAATCATCTCCAACCTGTGTAGTAAC	259		
Sbjct 200	CTGTTCTGTAGAAACTGTTGTTGAAAAACCTCCAGAATCATCTCCAACCTGTGTAGTAAC	259		
Query 260	ATCCATAGCATGATAAAGAGGAGCAAAACATTCCTCAAGTTAATTGCTGAAAGATAAACATT	319		
Sbjct 260	ATCCATAGCATGATAAAGAGGAGCAAAACATTCCTCAAGTTAATTGCTGAAAGATAAACATT	319		
Query 320	CACTCTGACATGGCAAGCAACGA	342		
Sbjct 320	CACTCTGACATGGGAAGCAACGA	342		

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