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
1. Groves, M.L. 1960. The isolation of a red protein from Milk. *J. Am. Chem. Soc.* 82(13): 3345-3350.
2. Johanson, B. 1960. Isolation of an iron-containing red protein from human milk. *Acta Chem. Scand.* 14(2): 510-512.
3. Baker, E.N. 1994. Structure and reactivity of transferrins. *Adv. Inorg. Chem.* 41: 389-464.
4. Lambert, L.A., Perri, H. and Meehan, T.J.

2005. Evolution of duplications in the transferrin family of proteins. *Comp. Biochem. Physiol. B* 140(1): 11-25.
5. Baveye, S., Elass, E., Mazurier, J., Spik, G. and Legrand, D. 1999. Lactoferrin: a multifunctional glycoprotein involved in the modulation of the inflammatory process. *Clin. Chem. Lab. Med.* 37(3): 281-286.
  6. 衛生福利部。2022。食品添加物使用範圍及限量暨規格標準。111年3月10日衛授食字第1101360128號公告修正。
  7. Ellingson, D.J., Shippar, J.J., Vennard, T.R., Moloney, C. and *et al.* 2019. Analytical method for lactoferrin in milk-based infant formulas by signature peptide quantification with ultra-high performance LC-tandem mass spectrometry. *J. AOAC Int.* 102(3): 915-925.
  8. Zhang, J., Lai, S., Cai, Z., Chen, Q. and *et al.* 2014. Determination of bovine lactoferrin in dairy products by ultra-high performance liquid chromatography–tandem mass spectrometry based on tryptic signature peptides employing an isotope-labeled winged peptide as internal standard. *Anal. Chim. Acta* 829: 33-39.
  9. Yuan, M., Feng, C., Wang, S., Zhang, W. and *et al.* 2017. Selection of possible signature peptides for the detection of bovine lactoferrin in infant formulas by LC-MS/MS. *PLoS One* 12(9): e0184152.
  10. Ke, X., Chen, Q., Pan, X., Zhang, J. and *et al.* 2016. Quantification of lactoferrin in breast milk by ultra-high performance liquid chromatography-tandem mass spectrometry with isotopic dilution. *RSC Adv.* 6(15): 12280-12285.
  11. van den Broek, I., Niessen, W.M. and van Dongen, W.D. 2013. Bioanalytical LC-MS/MS of protein-based biopharmaceuticals. *J. Chromatogr. B* 929: 161-179.
  12. Lange, V., Picotti, P., Domon, B. and Aebersold, R. 2008. Selected reaction monitoring for quantitative proteomics: a tutorial. *Mol. Syst. Biol.* 4(1): 222.
  13. 衛生福利部食品藥物管理署。2013。食品化學檢驗方法之確效規範。102年09月09日第二次修正。
  14. MacLean, B., Tomazela, D.M., Shulman, N., Chambers, M. and *et al.* 2010. Skyline: an open source document editor for creating and analyzing targeted proteomics experiments. *Bioinformatics* 26(7): 966-968.
  15. EU Reference Laboratories for Residues of Pesticides. SANTE/11312/2021. Analytical quality control and method validation procedures for pesticide residues analysis in food and feed. [[https://www.eurl-pesticides.eu/userfiles/file/EurlALL/SANTE\\_11312\\_2021.pdf](https://www.eurl-pesticides.eu/userfiles/file/EurlALL/SANTE_11312_2021.pdf)]
  16. Jenkins, R., Duggan, J. X., Aubry, A. F., Zeng, J. and *et al.* 2015. Recommendations for validation of LC-MS/MS bioanalytical methods for protein biotherapeutics. *AAPS J.* 17: 1-16.
  17. AOAC International. 2020. Standard Method Performance Requirements (SMPR®) for determination of bovine lactoferrin in infant and adult/pediatric nutritional formula. [<https://www.aoac.org/resources/smpr-2020-005>]
  18. 衛生福利部。2021。包裝食品營養標示應遵行事項。110年4月27日衛授食字第1101300478號公告修正。
  19. General European OMCL Network. 2016. Interpretation of screening results for unknown peptides and proteins by mass spectrometry based methods. PA/PH/OMCL (15) 04 2R.

# Method Development of Bovine Lactoferrin Analysis in Infant Formula

CHENG-HONG QIN, PAI-WEN WU, YUNG-LUN HSIEH, SHU-HAN CHANG,  
YA-MIN KAO, SU-HSIANG TSENG AND DER-YUAN WANG

Division of Research and Analysis, TFDA, MOHW

## ABSTRACT

Lactoferrin is an iron-binding glycoprotein commonly found in mammalian milk. To comply with the amended “Standards for Specification, Scope, Application and Limitation of Food Additives” announced by the Ministry of Health and Welfare of March 10, 2022, which stipulates that “Lactoferrin can be used in foods labeled with daily intake limits, and the total daily intake should not exceed 300 mg” and “This product can be used in special nutritional foods as needed”. A test method for bovine lactoferrin in infant formula products was developed. The protocol in brief was that bovine lactoferrin was reduced by dithiothreitol, alkylated by iodoacetamide, hydrolyzed by trypsin, and then analyzed by an ultra-high performance liquid chromatography/tandem mass spectrometer (UPLC-MS/MS). The average recoveries of bovine lactoferrin spiked at 4, 10 and 40 mg/100 g into the blank matrix made of camel milk nutrition powder were 91.0, 95.1 and 105.7%, and the coefficients of variance were 6.7, 3.2 and 6.6%, respectively in intra-day test (n=5); the average recoveries were 93.0, 99.4 and 112.2%, and the coefficients of variance were 7.8, 7.8 and 7.4%, respectively in inter-day test (n=10). We tested commercial nutrition supplement products for the application of the testing method, the results showed the detected values of bovine lactoferrin in 3 commercial products with bovine lactoferrin labeling were in the range of 120-253 mg/100 g. One hydrolyzed formula powder and 2 commercial products constituted of goat milk were not detected of bovine lactoferrin.

Key words: infant formula, bovine lactoferrin, UPLC-MS/MS