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Parameter	Condition		
LC system	ACQUITY UPLC I-Class with Sample manager (FL)		
Column	Atlantis Premier BEH C18 AX (100 × 2.1 mm, 1.7 μm)		
Delay column	Atlantis Premier BEH C18 AX (50 × 2.1 mm, 5 μm)		
Mobile phase	A: 2 mM Ammonium acetate in H ₂ O B: 0.1% (v/v) Ammonium hydroxide in methanol		
Gradient program	Time (min)	A (%)	B (%)
	0 → 2	99 → 99	1 → 1
	2 → 3	99 → 75	1 → 25
	3 → 8	75 → 50	25 → 50
	8 → 15	50 → 15	50 → 85
	15 → 16	15 → 0	85 → 100
	16 → 20	0 → 0	100 → 100
	20 → 20.1	0 → 100	100 → 0
	20.1 → 23.5	100 → 100	0 → 0
	23.5 → 24	100 → 99	0 → 1
Flow rate	0.3 mL/min		
Injection volume	20 μL		
Column temperature	40°C		
MS instrument	Xevo TQ-XS		
Ionization mode	ESI negative		
Capillary voltage	0.5 kV		
Desolvation temperature	400° C		
Desolvation gas flow	1,200 L/hr		
Cone gas flow	250 L/hr		
Ion source temp.	100° C		

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1. Brunn, H., Arnold, G., Körner, W., Rippen, G. *et al.* 2023. PFAS : forever chemicals - persistent, bioaccumulative and mobile. Reviewing the status and the need for their phase out and remediation of contaminated sites. *Environ. Sci. Eur.* 35: 1-50.
2. Androulakis, A., Alygizakis, N., Bizani, E. and Thomaidis, N.S. 2022. Current progress in the environmental analysis of poly-and perfluoroalkyl substances (PFAS). *Environ. Sci. Adv.* 1: 705-724.
3. Consumer Product Safety Commission. 2023. Characterizing PFAS chemistries, sources, uses, and regulatory trends in U.S. and international markets. [<https://www.cpsc.gov/s3fs-public/CPSC-PFAS-WhitePaper.pdf>]
4. Vorst, K.L., Saab, N., Silva, P., Curtzwiler, G. *et al.* 2021. Risk assessment of per-and polyfluoroalkyl substances (PFAS) in food: Symposium proceedings. *Trends Food Sci. Technol.* 116: 1203-1211.
5. European Food Safety Authority. 2020. Risk to human health related to the presence of perfluoroalkyl substances in food. *EFSA J.* 18: 1-391.
6. Seyyedsalehi, M.S. and Boffetta, P. 2023. Per-and poly-fluoroalkyl substances (PFAS) exposure and risk of kidney, liver, and testicular cancers: a systematic review and meta-analysis. *Med. Lav.* 114: 1-19.
7. United States Environmental Protection Agency. 2024. Final PFAS National Primary Drinking Water Regulation. [<https://www.epa.gov/sdwa/and-polyfluoroalkyl-substances-pfas>]
8. European Commission. 2019. Regulation (EC) No 2019/1021 of the European Parliament and of the Council of 20 June 2019 on persistent organic pollutants (recast). *Off. J. Eur. Union* L169: 45-77.
9. European Commission. 2020. Directive (EC) No 2020/2184 of the European Parliament and of the Council of 16 December 2020 on the quality of water intended for human consumption (recast). *Off. J. Eur. Union* L435: 1-62.
10. European Commission. 2023. Regulation (EC) No 2023/915 of 25 April 2023 on maximum levels for certain contaminants in food and repealing Regulation (EC) No 1881/2006. *Off. J. Eur. Union* L119: 103-157.
11. 環境部。2024。飲用水水質標準。113年11月25日環部水字第1131072826號修正公告。 [<https://law.moj.gov.tw/LawClass/LawAll>].

- aspx?pcode=O0040019]
12. United States Food and Drug Administration. 2024. Determination of 30 per and polyfluoroalkyl substances (PFAS) in food and feed using liquid chromatography-tandem mass spectrometry (LC-MS/MS). C-010.03. [https://www.fda.gov/media/131510/download]
13. 衛生福利部食品藥物管理署。2021。食品化學檢驗方法之確效規範。110年11月01日修正公布。

Development of an Analytical Method for Per- and Polyfluoroalkyl Substances in Foods

CHIA-HSIN LIU, MIN-CHIH YUAN, YU WU, GUAN-JHIH PENG,
SHU-HAN CHANG, YA-MIN KAO, MEI-CHIH LIN
AND SU-HSIANG TSENG

Division of Research and Analysis, TFDA, MOHW

ABSTRACT

Per- and polyfluoroalkyl substances (PFAS) are widely used across industries due to their exceptional stability as well as hydrophobic and oleophobic properties. However, these compounds are highly persistent to decompose which carrying risk of bioaccumulation. The European Food Safety Authority (EFSA) proposed in 2020 that PFAS may cause developmental effects and adverse impacts on serum cholesterol, liver function and immune system. EFSA proposed a tolerable weekly intake (TWI) of 4.4 ng/kg body weight/week for the combined exposure to PFOA, PFNA, PFHxS and PFOS. The European Union set the maximum levels for PFAS in food stuffs in 2023. 0.2 to 50 µg/kg for individual PFAS and 1.3 to 50 µg/kg for the sum of the four regulated compounds. This study established a method for the determination of 19 PFAS in fish, meat, eggs, milk, and edible offal. The approach combined QuEChERS pretreatment, an EMR lipid solid-phase extraction cartridge, and a delay column to eliminate background interference from the analytical system in liquid chromatograph-tandem mass spectrometry (LC-MS/MS) analysis. The results of the spiking tests showed that the average recoveries of 19 PFAS in fish were 77.2-97.9% with the coefficients of variation of 2.8-16.3%; those in eggs were 73.7-91.1% with the coefficients of variation of 4.6-18.3%; those in milk were 63.7-97.7% with the coefficients of variation of 3.2-14.9%; those in meat were 71.7-96.3% with the coefficients of variation of 3.3-15.1%; those in edible offal were 70.0-102.2% with the coefficients of variation of 2.7-20.6%. All results met the “Validation guideline of the food chemical testing method” published by the Taiwan Food and Drug Administration. The limits of quantification were 0.02-0.25 µg/kg, and the linear range of the standard calibration curves calibrated by the internal standards were 0.005-0.5 ng/mL with the coefficients of correlation higher than 0.99. The developed method is sensitive, stable, suitable for various applications across diverse food matrices, offering good accuracy and precision suitable for routine monitoring and public health protection.

Key words: per- and polyfluoroalkyl substances, QuEChERS pretreatment, EMR lipid SPE cartridge, LC-MS/MS