

Method of Test for Synthetic Phenethylamines in Urine (1)

1. Scope

This method is applicable to the determination of 74 phenethylamines (amphetamine etc. listed as the attached table) in urine.

2. Method

After dilution, analytes are determined by liquid chromatography/tandem mass spectrometry (LC-MS/MS).

2.1. Equipment

2.1.1. Liquid chromatograph/tandem mass spectrometer

2.1.1.1. Ion source: electrospray ionization (ESI).

2.1.1.2. Column: Kinetex® Phenyl-hexyl, 1.7 μm , 2.1 mm i.d. $\times$ 10 cm, or an equivalent product.

2.1.2. Vortex mixer.

2.2. Chemicals

Methanol, HPLC grade;

Formic acid and ammonium acetate, reagent grade;

Artificial urine (UTAK 88121-CDF(L) or an equivalent product), reagent grade;

Deionized water, resistivity $\geq 18 \text{ M}\Omega\cdot\text{cm}$ (at 25°C);

Amphetamine etc. listed in the attached table, reference standards;

Amphetamine- d_8 and other isotope-labeled internal standards (listed in the attached table).

2.3. Apparatus

2.3.1. Volumetric flask: 1 mL and 10 mL.

2.3.2. Membrane filter: 0.22 μm , PVDF.

2.4. 50% Methanol:

Mix methanol and deionized water at the ratio of 1:1 (v/v).

2.5. Mobile phase

2.5.1. Solvent A

Dissolve and dilute 0.4 g of ammonium acetate with deionized water to 1000 mL, add 1 mL of formic acid, and filter with a membrane filter.

2.5.2. Solvent B:

Dilute 1 mL of formic acid with methanol to 1000 mL, and filter with a membrane filter.

2.6. Internal standard solution preparation

Transfer 1 mg of internal standards accurately weighed into each 10-mL

volumetric flask, dissolve and dilute with methanol to volume as the internal standard stock solutions. Store at -20°C in the dark. Prior to use, mix appropriate volume of the internal standard stock solutions, and dilute with 50% methanol to 500 ng/mL as the internal standard solution.

2.7. Standard solution preparation

Transfer 1 mg of reference standards accurately weighed into each 10-mL volumetric flask, dissolve and dilute with methanol to volume as the standard stock solution. Store at -20°C in the dark. Prior to use, mix appropriate volume of the standard stock solutions, and dilute with 50% methanol to 500 ng/mL as the standard solution.

2.8. Sample solution preparation

Transfer 20 µL of the homogenized sample and 20 µL of the internal standard solution into a 1-mL volumetric flask, and dilute to volume with 50% methanol. Filter with a membrane filter, and take the filtrate as the sample solution.

2.9. Calibration curve

Use the artificial urine as the blank sample. Separately take 20 µL of the artificial urine, add 2-100 µL of the standard solution and 20 µL of the internal standard solution, dilute with 50% methanol to 1 mL, and filter with a membrane filter. Take the filtrates as the calibration standard solutions. Operate LC-MS/MS according to the following conditions. Establish the calibration curve of each phenethylamine by the ratios of the peak area of each phenethylamine to that of the respective internal standard vs. the added concentrations (1-50 ng/mL).

LC-MS/MS operating conditions^(note):

Column: Kinetex® Phenyl-hexyl, 1.7 µm, 2.1 mm i.d. × 10 cm.

Column temperature: 40°C.

Injection volume: 3 µL.

Mobile phase: a gradient program of solvent A and solvent B is as follows.

Time (min)	A (%)	B (%)
0.0 → 0.5	95 → 95	5 → 5
0.5 → 1.0	95 → 70	5 → 30
1.0 → 1.5	70 → 70	30 → 30
1.5 → 2.0	70 → 63	30 → 37
2.0 → 2.5	63 → 63	37 → 37
2.5 → 2.6	63 → 60	37 → 40

2.6 → 3.0	60 → 60	40 → 40
3.0 → 5.0	60 → 54	40 → 46
5.0 → 5.5	54 → 54	46 → 46
5.5 → 8.0	54 → 50	46 → 50
8.0 → 9.5	50 → 35	50 → 65
9.5 → 10.5	35 → 0	65 → 100
10.5 → 11.0	0 → 0	100 → 100
11.0 → 11.1	0 → 95	100 → 5
11.1 → 14.0	95 → 95	5 → 5

Flow rate: 0.3 mL/min.

Ion spray voltage: 5.5 kV.

Ionization mode: ESI⁺.

Turbo heater temperature: 550°C.

Nebulizer gas, GS1: 50 psi.

Heated gas, GS2, 60 psi.

Curtain gas: 30 psi.

Collision gas: Medium.

Detection mode: multiple reaction monitoring (MRM). Detection of ion pair, declustering potential and collision energy are shown in the attached table.

Note: All the parameters can be adjusted depending on the instruments used if the above conditions are not applicable.

2.10. Identification and quantification

Accurately inject 3 µL of the sample solution and the calibration standard solutions into LC-MS/MS separately. Operate according to the conditions in section 2.9. Identify each phenethylamine based on the retention time and the relative ion intensities^(note). Calculate the amount of each phenethylamine in the sample by the following formula:

$$\text{The amount of each phenethylamine in the sample (ng/mL)} = \frac{C \times V}{M}$$

Where,

C: the concentration of each phenethylamine in the sample solution calculated by the calibration curve (ng/mL)

V: the final make-up volume of the sample (mL)

M: the volume of the sample (mL)

Note: Relative ion intensities are calculated by peak areas of quantitative ions divided by peak areas of qualitative ions ($\leq 100\%$). Maximum permitted tolerances of relative ion intensities are as the following:

Relative ion intensity (%)	Tolerance (%)
> 50	± 20
> 20~50	± 25
> 10~20	± 30
≤ 10	± 50

Remark

1. Limit of quantification (LOQ) for each phenethylamine is 50 ng/mL.
2. Further validation should be performed when interference compounds appear in the samples.

Reference

Dong, Y., Yan, K., Ma, Y., Wang, S., He, G., Deng, J. and Yang, Z. 2015. A sensitive dilute-and-shoot approach for the simultaneous screening of 71 stimulants and 7 metabolites in human urine by LC–MS-MS with dynamic MRM. J. Chromatogr. Sci. 53: 1528-1535.

Reference chromatogram




Figure. The MRM chromatograms of 74 phenethylamines and 10 isotope-labeled internal standards analyzed by LC-MS/MS.




Figure. The MRM chromatograms of 74 phenethylamines and 10 isotope-labeled internal standards analyzed by LC-MS/MS (continued).




Figure. The MRM chromatograms of 74 phenethylamines and 10 isotope-labeled internal standards analyzed by LC-MS/MS (continued).




Figure. The MRM chromatograms of 74 phenethylamines and 10 isotope-labeled internal standards analyzed by LC-MS/MS (continued).




Figure. The MRM chromatograms of 74 phenethylamines and 10 isotope-labeled internal standards analyzed by LC-MS/MS (continued).




Figure. The MRM chromatograms of 74 phenethylamines and 10 isotope-labeled internal standards analyzed by LC-MS/MS (continued).

Table. MRM parameters of 74 phenethylamines and 10 isotope-labeled internal standards

Analyte	Ion pair	Declustering potential (V)	Collision energy (eV)	Internal standard
	Precursor ion(<i>m/z</i>) >product ion (<i>m/z</i>)			
2-(2,5-Dimethoxy-4-nitrophenyl)- <i>N</i> -(2-methoxybenzyl)ethanamine (25N-NBOMe)	347 > 121*	30	23	DOB-d ₅
	347 > 91	30	59	
2-(2,5-Dimethoxyphenyl)- <i>N</i> -(2-methoxybenzyl)ethanamine (25H-NBOMe)	302 > 121*	33	23	DOB-d ₅
	302 > 91	33	51	
2,5-Dimethoxy-4- <i>n</i> -propylthiophenethylamine (2C-T-7)	256 > 167*	67	35	DOB-d ₅
	256 > 224	67	24	
<i>N,N</i> -Dimethylamphetamine (<i>N,N</i> -DMA)	164 > 91*	57	30	MDA-d ₅
	164 > 119	57	18	
<i>N,N</i> -Dimethyl-3,4-methylenedioxyamphetamine (MDDMA)	208 > 163*	21	19	MDMA-d ₅
	208 > 135	21	29	
Escaline	226 > 181*	60	22	Mescaline-d ₉
	226 > 91	60	44	
5-(2-Ethylaminopropyl)benzofuran (5-EAPB)	204 > 131*	34	29	5-APB-d ₅
	204 > 159	34	9	
<i>N</i> -Ethylamphetamine	164 > 119*	68	16	MDA-d ₅
	164 > 65	68	56	
4-Ethyl- <i>N</i> -(<i>o</i> -methoxybenzyl)amphetamine (4-EA-NBOMe)	284 > 121*	44	24	25I-NBOMe-d ₃
	284 > 91	44	59	
4-(Ethylthio)-2,5-dimethoxy- <i>N</i> -[(2-methoxyphenyl)methyl]benzeneethanamine (25T2-NBOMe)	362 > 121*	23	27	25I-NBOMe-d ₃
	362 > 91	23	65	
Fenproporex	189 > 119*	40	15	Amphetamine-d ₈
	189 > 65	40	59	
4-Fluoroamphetamine (4-FA)	154 > 109*	40	28	Amphetamine-d ₈
	154 > 137	40	15	
<i>N</i> -(2-Fluorobenzyl)-4-iodo-2,5-dimethoxyphenethylamine (25I-NBF)	416 > 291*	47	25	25I-NBOMe-d ₃
	416 > 276	47	39	
4-Fluoroethamphetamine (4-FEA)	182 > 109*	65	31	MDEA-d ₅
	182 > 137	65	17	
4-Fluoromethamphetamine (4-FMA)	168 > 109*	50	30	Amphetamine-d ₈
	168 > 83	50	52	

* Quantitative ion pair.

Table. MRM parameters of 74 phenethylamines and 10 isotope-labeled internal standards (continued)

Analyte	Ion pair Precursor ion(<i>m/z</i>) >product ion (<i>m/z</i>)	Declustering potential (V)	Collision energy (eV)	Internal standard
3-Fluoro-4-methoxyamphetamine (3-F-4-MOA)	184 > 139* 184 > 167	36 36	24 35	Amphetamine-d ₈
5-Fluoro-2-methoxyamphetamine (5-F-2-MOA)	184 > 139* 184 > 109	40 40	25 14	Amphetamine-d ₈
1-(4-Fluorophenyl)butan-2-amine (FPBA)	168 > 109* 168 > 83	47 47	28 51	5-APB-d ₅
<i>N</i> -Hydroxy-3,4- methylenedioxyamphetamine (<i>N</i> -Hydroxy-MDA)	196 > 163* 196 > 135	34 34	15 28	Amphetamine-d ₈
Lefetamine	226 > 103* 226 > 165	61 61	41 38	DOB-d ₅
Mescaline	212 > 165* 212 > 180	40 40	30 24	Mescaline-d ₉
Methamphetamine	150 > 119* 150 > 65	40 40	14.5 55.5	Methamphetamine-d ₈
4-Methoxyamphetamine (PMA)	166 > 121* 166 > 91	39 39	25 42	Amphetamine-d ₈
<i>N</i> -(2-Methoxybenzyl)-4-bromo-2,5- dimethoxyphenethylamine (25B-NBOMe)	380 > 65* 380 > 93	72 72	105 44	25I-NBOMe-d ₃
<i>N</i> -(<i>o</i> -Methoxybenzyl)-3,4- dimethoxyamphetamine (3,4-DMA-NBOMe)	316 > 121* 316 > 91	17 17	24 55	25I-NBOMe-d ₃
<i>N</i> -(2-Methoxybenzyl)-4-iodo-2,5- dimethoxyphenethylamine (25I-NBOMe)	428 > 65* 428 > 272	117 117	113 23	25I-NBOMe-d ₃
Methoxyethylamphetamine (PMEA)	194 > 149* 194 > 121	42 42	18 27	Methamphetamine-d ₈
Methoxymethamphetamine (PMMA)	180 > 78* 180 > 65	33 33	55 62	Amphetamine-d ₈
5-Methoxy-3,4-methylenedioxyamphetamine (MMDA)	210 > 135* 210 > 165	51 51	27 26	Amphetamine-d ₈
5-(2-Methylaminopropyl)benzofuran (5-MAPB)	190 > 131* 190 > 159	27 27	26 16	5-APB-d ₅

* Quantitative ion pair.

Table. MRM parameters of 74 phenethylamines and 10 isotope-labeled internal standards (continued)

Analyte	Ion pair	Declustering potential (V)	Collision energy (eV)	Internal standard
	Precursor ion(<i>m/z</i>) >product ion (<i>m/z</i>)			
5-(2-Methylaminopropyl)-2,3- dihydrobenzofuran (5-MAPDB)	192 > 161*	15	19	5-APB-d ₅
	192 > 133	15	33	
4-Methylamphetamine (4-MA)	150 > 105*	42	26	Amphetamine-d ₈
	150 > 133	42	12	
4-Methyl-2,5-dimethoxyamphetamine (DOM)	210 > 165*	63	24	DOB-d ₅
	210 > 178	63	25	
3,4-Methylenedioxyamphetamine (MDA)	180 > 105*	51	27	MDA-d ₅
	180 > 135	51	25	
3,4-Methylenedioxy- <i>N</i> - ethylamphetamine (MDEA)	208 > 163*	57	17	MDEA-d ₅
	208 > 135	57	30	
3,4-Methylenedioxymethamphetamine (MDMA)	194 > 163*	54	16	MDMA-d ₅
	194 > 105	54	36	
<i>N</i> -Methyl- α -ethyl-3,4- methylenedioxyphenethylamine (MBDB)	208 > 135*	60	25	5-APB-d ₅
	208 > 51	60	83	
4-Methylthioamphetamine (4-MTA)	182 > 137*	44	24	DOB-d ₅
	182 > 115	44	50	
Phentermine	150 > 65*	38	52	Phentermine-d ₅
	150 > 133	38	14	
4-Propoxy-3,5-dimethoxyamphetamin (3C-P)	254 > 237*	54	12	DOB-d ₅
	254 > 107	54	38	
Proscaline	240 > 181*	23	21	DOB-d ₅
	240 > 121	23	33	
2,4,5-Trimethoxyamphetamine (TMA-2)	226 > 194*	47	28	25I-NBOMe-d ₃
	226 > 179	47	34	
2,4,6-Trimethoxyamphetamine (TMA-6)	226 > 209*	72	16	25I-NBOMe-d ₃
	226 > 181	72	28	
3,4,5-Trimethoxyamphetamine (TMA)	226 > 209*	64	14	25I-NBOMe-d ₃
	226 > 181	64	25	

* Quantitative ion pair.

Table. MRM parameters of 74 phenethylamines and 10 isotope-labeled internal standards (continued)

Analyte	Ion pair	Declustering potential (V)	Collision energy (eV)	Internal standard
	Precursor ion(<i>m/z</i>) >product ion (<i>m/z</i>)			
5-(2-Aminopropyl)benzofuran-d ₅ (5-APB-d ₅) (I.S.)	181 > 164	27	12	-
Amphetamine-d ₈ (I.S.)	144 > 97	50	24	-
4-Bromo-2,5-dimethoxyamphetamine-d ₅ (DOB-d ₅) (I.S.)	279 > 262	67	16	-
4-Iodo-2,5-dimethoxy- <i>N</i> -[(2-methoxyphenyl)methyl] benzeneethanamine-d ₃ (25I-NBOMe-d ₃) (I.S.)	431 > 124	26	26	-
Mescaline-d ₉ (I.S.)	221 > 204	59	16	-
Methamphetamine-d ₈ (I.S.)	158 > 93	36	33	-
3,4-Methylenedioxyamphetamine-d ₅ (MDA-d ₅) (I.S.)	185 > 168	58	14	-
3,4-Methylenedioxy- <i>N</i> -ethylamphetamine-d ₅ (MDEA-d ₅) (I.S.)	213 > 163	61	18	-
3,4-Methylenedioxymethamphetamine-d ₅ (MDMA-d ₅) (I.S.)	199 > 165	22	19	-
Phentermine-d ₅ (I.S.)	155 > 96	42	32	-