

Method of Test for Synthetic Cannabinoids in Urine (1)

1. Scope

This method is applicable to the determination of 39 synthetic cannabinoids (AB-CHMINACA etc. listed as the attached table) in urine.

2. Method

After hydrolysis and purification, analytes are determined by liquid chromatography/tandem mass spectrometry (LC-MS/MS).

2.1. Equipment

2.1.1. Liquid chromatograph/tandem mass spectrometer

2.1.1.1. Ion source: electrospray ionization (ESI)

2.1.1.2. Column: Sunshell® RP-AQUA, 2.6 µm, 2.1 mm i.d. × 10 cm, or an equivalent product

2.1.2. Vortex mixer

2.1.3. Ultrasonicator

2.1.4. Reciprocal shaking water bath

2.1.5. Nitrogen evaporator

2.1.6. Solid phase vacuum extraction manifold

2.2. Chemicals

Methanol and acetonitrile, HPLC grade;

β-glucuronidase (100,000 units/mL);

Acetic acid, analytical grade;

Ethyl acetate, sodium acetate (CH₃COONa), reagent grade;

Formic acid, mass spectrometry grade;

Artificial urine (UTAK 88121-CDF(L) or an equivalent product), reagent grade;

Deionized water, resistivity ≥ 18 MΩ·cm (at 25°C);

AB-CHMINACA etc. listed in the attached table, reference standards;

AB-CHMINACA-d₄ and other isotope-labeled internal standards (listed in the attached table).

2.3. Apparatus

2.3.1. Volumetric flask: 1 mL and 10 mL

2.3.2. Centrifuge tube: 15 mL, PP

2.3.3. Membrane filter: 0.22 µm, PVDF

2.3.4. Solid-supported liquid-liquid extraction cartridge: Novum™ SLE cartridge, 3 mL, or an equivalent product

2.4. Reagent solution preparation

2.4.1. 0.1 M sodium acetate solution

Dissolve and dilute 0.82 g of sodium acetate with deionized water to 100 mL.

2.4.2. 20% acetic acid solution

Dilute 20 mL of acetic acetate with 80 mL of deionized water to 100 mL.

2.5. Mobile phase

2.5.1. Solvent A

Dilute 1 mL of formic acid with deionized water to 1000 mL, mix well and filter with a membrane filter.

2.5.2. Solvent B

Dilute 1 mL of formic acid with methanol/acetonitrile solution (1:1, v/v) to 1000 mL, mix well and filter with a membrane filter.

2.6. Internal standard solution preparation

Transfer 1 mg of the 18 internal standards accurately weighed into each 10-mL volumetric flask, dissolve and dilute with methanol to volume as the internal standard stock solutions. Store at -20°C in the dark. Upon use, mix adequate volume of the internal standard stock solutions, and dilute with methanol to 1 µg/mL as the internal standard solution.

2.7. Standard solution preparation

Transfer 1 mg of the 39 reference standards accurately weighed into each 10-mL volumetric flask, dissolve and dilute with methanol to volume as the standard stock solutions. Store at -20°C in the dark. Upon use, mix adequate volume of the standard stock solutions, and dilute with methanol to 1 µg/mL as the standard solution.

2.8. Sample solution preparation

2.8.1. Hydrolysis

Mix 1 mL of the homogenized sample accurately with 20 µL of the internal standard solution, and 1 mL of 0.1 M sodium acetate solution in a 15 mL centrifuge tube. Add 500 µL of β-glucuronidase and mix well. Transfer the above solution into the reciprocal shaking water bath and incubate at 56 °C with 80 rpm for 15 min. Allow the solution to cool to room temperature and add 500 µL of 20% acetic acid solution for subsequent purification.

2.8.2. Purification

Transfer accurate 400 µL of the hydrolyzed sample solution mentioned in

section 2.8.1. into the solid-supported liquid-liquid extraction cartridge and hold for 5 min. Elute the cartridge with 0.9 mL of ethyl acetate for twice. Collect the eluent and evaporate to dryness by gently flushing with a stream of nitrogen at 40°C. Dissolve the residue with 1 mL of methanol and filter with a membrane filter as the sample solution.

2.9. Calibration curve

Use the artificial urine as the blank sample. Separately take 5-100 µL of the standard solution and mix with artificial urine to 1 mL. Prepare the calibration solutions following the procedure in section 2.8. Operate LC-MS/MS according to the following conditions. Establish the calibration curve of each synthetic cannabinoid by the ratios of the peak area of each synthetic cannabinoid to that of the internal standard vs. the added concentrations (5-100 ng/mL).

LC-MS/MS operating conditions⁽¹⁾:

Column: SunShell® RP-AQUA, 2.6 µm, 2.1 mm i.d. × 10 cm

Column temperature: 40°C

Injection volume: 3 µL

Mobile phase: a gradient program of solvent A and solvent B is as follows.

Time (min)	A (%)	B (%)
0.0 → 0.5	60 → 60	40 → 40
0.5 → 12.0	60 → 10	40 → 90
12.0 → 12.5	10 → 10	90 → 90
12.5 → 12.6	10 → 60	90 → 40
12.6 → 14.5	60 → 60	40 → 40

Flow rate: 0.4 mL/min

Ion spray voltage:

Electrospray ionization positive mode (ESI⁺): 5.5 KV

Electrospray ionization negative mode (ESI⁻): -4.5 KV

Turbo heater temperature: 550°C

Nebulizer gas (GS1): 50 psi

Heated gas (GS2): 60 psi

Curtain gas: 30 psi

Collision gas: Medium

Detection mode: multiple reaction monitoring (MRM). Selected ion pairs, declustering potential and collision energy are shown in

the attached table.

Note 1: All the parameters can be adjusted depending on the instruments used if the above conditions are not applicable.

2.10. Identification and quantification

Accurately inject 3 µL of the sample solution and the calibration standard solutions into LC-MS/MS separately. Operate according to the conditions in section 2.9. Identify each synthetic cannabinoid based on the retention time and the relative ion intensities⁽²⁾. Calculate the amount (ng/mL) of each synthetic cannabinoid in the sample by the following formula:

$$\text{The amount of each synthetic cannabinoid in the sample (ng/mL)} = \frac{C \times V}{M}$$

Where:

C: the concentration of each synthetic cannabinoid in the sample solution calculated by the calibration curve (ng/mL)

V: the final make-up volume of the sample (mL)

M: the volume of the sample (mL)

Note 2: Relative ion intensities are calculated by peak areas of quantitative ions divided by peak areas of qualitative ions ($\leq 100\%$). Maximum permitted tolerances of relative ion intensities are as the following:

Relative ion intensity (%)	Tolerance (%)
> 50	± 20
> 20-50	± 25
> 10-20	± 30
≤ 10	± 50

Remark

1. Limit of quantification (LOQ) for each synthetic cannabinoid is 50 ng/mL.
2. Further validation should be performed when interference compounds appear in the samples.

Reference

Scheidweiler, K. B. and Huestis, M. A. 2014. Simultaneous quantification of 20 synthetic cannabinoids and 21 metabolites, and semi-quantification of 12 alkyl hydroxy metabolites in human urine by liquid chromatography–tandem mass spectrometry. J. Chromatogr. A 1327: 105-117.

Reference chromatograms




Figure. The MRM chromatograms of 39 synthetic cannabinoids and 18 internal standards analyzed by LC/MS/MS.




Figure. The MRM chromatograms of 39 synthetic cannabinoids and 18 internal standards analyzed by LC/MS/MS (continued).




Figure. The MRM chromatograms of 39 synthetic cannabinoids and 18 internal standards analyzed by LC/MS/MS (continued).




Figure. The MRM chromatograms of 39 synthetic cannabinoids and 18 internal standards analyzed by LC/MS/MS (continued).




Figure. The MRM chromatograms of 39 synthetic cannabinoids and 18 internal standards analyzed by LC/MS/MS (continued).

Table. MRM parameters of 39 synthetic cannabinoids and 18 internal standards.

Analyte	Ion pairs	Declustering potential (V)	Collision energy (eV)	Internal Standard
	Precursor ion (<i>m/z</i>)			
	> product ion (<i>m/z</i>)			
AB-CHMINACA	357 > 241*	59	35	AB-CHMINACA-d ₄
	357 > 312	59	23	
AB-CHMINACA metabolite M4	259 > 241*	53	21	JWH-018 NPA-d ₅
	259 > 145	53	34	
AB-FUBINACA	369 > 253*	83	32	AB-FUBINACA-d ₄
	369 > 324	83	21	
AB-FUBINACA metabolite 3	370 > 109*	76	70	UR-144 NPA-d ₅
	370 > 253	76	31	
AB-PINACA	331 > 215*	40	36	AB-PINACA-d ₉
	331 > 286	40	20	
AB-PINACA <i>N</i> -(4-hydroxypentyl) metabolite (AB-PINACA N4HP)	347 > 330*	65	14	UR-144 NPA-d ₅
	347 > 213	65	38	
AB-PINACA 5-pentanoic acid metabolite (AB-PINACA NPA)	361 > 316*	64	22	UR-144 NPA-d ₅
	361 > 217	64	42	
ADB-FUBINACA	383 > 338*	72	21	AB-FUBINACA-d ₄
	383 > 253	72	36	
ADBICA <i>N</i> -pentanoic acid metabolite (ADBICA NPA)	377 > 244*	58	30	UR-144 NPA-d ₅
	377 > 144	58	52	
AKB48 <i>N</i> -(4-hydroxypentyl) metabolite (AKB48 N4HP)	382 > 135*	56	30	JWH-018 NPA-d ₅
	382 > 93	56	75	
AKB48 <i>N</i> -pentanoic acid metabolite (AKB48 NPA)	396 > 135*	75	29	UR-144 NPA-d ₅
	396 > 93	75	73	
AM-2201	360 > 155*	117	37	AM-2201-d ₅
	360 > 127	117	70	
AM-2201 <i>N</i> -(4-hydroxypentyl) metabolite (AM-2201 N4HP)	376 > 155*	76	34	AM-2201 N4HP-d ₅
	376 > 127	76	61	
CP 47,497-C8-homolog C-8-hydroxy metabolite (8-HydroxyCP-47497)	349 > 175*	59	19	UR-144 NPA-d ₅
	349 > 331	59	12	
4F-MDMB-BINACA	364 > 219*	60	33	AB-FUBINACA-d ₄
	364 > 304	60	23	

*Quantitative ion

Table. MRM parameters of 39 synthetic cannabinoids and 18 i internal standards (continued).

Analyte	Ion pairs	Declustering potential (V)	Collision energy (eV)	Internal Standard
	Precursor ion (m/z) > product ion (m/z)			
5F-MDMB-PICA	377 > 232*	89	32	AM-2201-d ₅
	377 > 144	89	48	
5-fluoro MDMB-PINACA metabolite 7 (5F-ADB metabolite 7)	364 > 233*	89	32	UR-144 NPA-d ₅
	364 > 318	89	22	
5-fluoro AKB48 <i>N</i> -(4-hydroxypentyl) metabolite (5F-AKB48 N4HP)	400 > 135*	136	28	UR-144 NPA-d ₅
	400 > 93	136	74	
(S)-5F-AMB acid metabolite (5F-AMB acid metabolite)	350 > 233*	120	28	UR-144 N5HP-d ₅
	350 > 145	120	53	
FUB-PB-22 3-carboxyindole metabolite (FUB-PB-22 3CI)	270 > 109*	68	25	JWH-018 NPA-d ₅
	270 > 83	68	81	
JWH-022	340 > 155*	142	32	JWH-018-d ₉
	340 > 127	142	60	
JWH-073	328 > 155*	151	32	JWH-073-d ₇
	328 > 127	151	65	
JWH-250	336 > 121*	129	27	JWH-250-d ₅
	336 > 91	129	65	
JWH-073 <i>N</i> -butanoic acid metabolite (JWH-073 NBA)	358 > 155*	55	30	JWH-203 NPA-d ₅
	358 > 127	55	64	
JWH-018 <i>N</i> -(5-hydroxypentyl) metabolite (JWH-018 N5HP)	358 > 155*	75	29	JWH-018 NPA-d ₅
	358 > 127	75	67	
JWH-019 <i>N</i> -(6-hydroxyhexyl) metabolite (JWH-019 N6HH)	372 > 155*	91	32	UR-144 NPA-d ₅
	372 > 127	91	72	
JWH-073 <i>N</i> -(3-hydroxybutyl) metabolite (JWH-073 N3HB)	344 > 155*	40	34	JWH-203 NPA-d ₅
	344 > 127	40	60	
JWH-081 <i>N</i> -(5-hydroxypentyl) metabolite (JWH-081 N5HP)	388 > 185*	152	33	JWH-081 N5HP-d ₅
	388 > 157	152	58	
JWH-122 <i>N</i> -(5-hydroxypentyl) metabolite (JWH-122 N5HP)	372 > 169*	85	32	JWH-122 N5HP-d ₅
	372 > 115	85	97	
JWH-210 <i>N</i> -(5-hydroxypentyl) metabolite (JWH-210 N5HP)	386 > 183*	145	34	JWH-210 N4HP-d ₅
	386 > 153	145	65	

*Quantitative ion.

Table. MRM parameters of 39 synthetic cannabinoids and 18 internal standards (continued).

Analyte	Ion pairs	Declustering potential (V)	Collision energy (eV)	Internal Standard
	Precursor ion (<i>m/z</i>) > product ion (<i>m/z</i>)			
JWH-250 <i>N</i> -(4-hydroxypentyl) metabolite (JWH-250 N4HP)	352 > 121* 352 > 91	86 86	28 67	JWH-250 N4HP-d ₅
JWH-018 <i>N</i> -pentanoic acid metabolite (JWH-018 NPA)	372 > 155* 372 > 127	72 72	33 74	JWH-018 NPA-d ₅
JWH-203 <i>N</i> -pentanoic acid metabolite (JWH-203 NPA)	370 > 125* 370 > 200	135 135	34 27	JWH-203 NPA-d ₅
JWH-210 <i>N</i> -pentanoic acid metabolite (JWH-210 NPA)	400 > 183* 400 > 153	70 70	34 70	JWH-210 N4HP-d ₅
MDMB-CHMICA metabolite M2	371 > 240* 371 > 144	79 79	28 48	UR-144 NPA-d ₅
(±)-11-nor-9-carboxy-Δ ⁹ -THC (THCA)	343 > 245* 343 > 191	-140 -140	-38 -42	THCA-d ₉
UR-144 <i>N</i> -(4-hydroxypentyl) metabolite (UR-144 N4HP)	328 > 125* 328 > 97	94 94	27 34	UR-144 N5HP-d ₅
UR-144 <i>N</i> -pentanoic acid metabolite (UR-144 NPA)	342 > 125* 342 > 244	47 47	30 32	UR-144 NPA-d ₅
XLR-11 <i>N</i> -(4-hydroxypentyl) metabolite (XLR-11 N4HP)	346 > 125* 346 > 248	94 94	28 33	XLR-11 N4HP-d ₅
AB-CHMINACA-d ₄ (I.S.)	361 > 245	114	35	-
AB-FUBINACA-d ₄ (I.S.)	373 > 328	68	21	-
AB-PINACA-d ₉ (I.S.)	340 > 224	40	36	-
AM-2201-d ₅ (I.S.)	365 > 127	31	73	-
AM-2201 <i>N</i> -(4-hydroxypentyl) metabolite-d ₅ (AM-2201 N4HP-d ₅) (I.S.)	381 > 155	88	35	-
JWH-018-d ₉ (I.S.)	351 > 155	36	33	-
JWH-073-d ₇ (I.S.)	335 > 155	76	33	-
JWH-250-d ₅ (I.S.)	341 > 121	46	27	-
JWH-081 <i>N</i> -(5-hydroxypentyl) metabolite-d ₅ (JWH-081 N5HP-d ₅) (I.S.)	393 > 185	118	32	-
JWH-122 <i>N</i> -(5-hydroxypentyl) metabolite-d ₅ (JWH-122 N5HP-d ₅) (I.S.)	377 > 169	91	33	-

*Quantitative ion

Table. MRM parameters of 39 synthetic cannabinoids and 18 internal standards (continued).

Analyte	Ion pairs	Declustering potential (V)	Collision energy (eV)	Internal Standard
	Precursor ion (<i>m/z</i>) > product ion (<i>m/z</i>)			
JWH-210 <i>N</i> -(4-hydroxypentyl) metabolite- <i>d</i> ₅ (JWH-210 N4HP- <i>d</i> ₅) (I.S.)	391 > 183	111	32	—
JWH-250 <i>N</i> -(4-hydroxypentyl) metabolite- <i>d</i> ₅ (JWH-250 N4HP- <i>d</i> ₅) (I.S.)	357 > 121	76	29	—
JWH-018 <i>N</i> -pentanoic acid metabolite- <i>d</i> ₅ (JWH-018 NPA- <i>d</i> ₅) (I.S.)	377 > 155	64	34	—
JWH-203 <i>N</i> -pentanoic acid metabolite- <i>d</i> ₅ (JWH-203 NPA- <i>d</i> ₅) (I.S.)	375 > 125	73	35	—
(±)-11-nor-9-carboxy-Δ ⁹ -THC- <i>d</i> ₉ (THCA- <i>d</i> ₉) (I.S.)	352 > 254	-160	-38	—
UR-144 <i>N</i> -(5-hydroxypentyl) metabolite- <i>d</i> ₅ (UR-144 N5HP- <i>d</i> ₅) (I.S.)	333 > 125	74	25	—
UR-144 5-pentanoic acid metabolite- <i>d</i> ₅ (UR-144 NPA- <i>d</i> ₅) (I.S.)	347 > 125	64	31	—
XLR-11 <i>N</i> -(4-hydroxypentyl) metabolite- <i>d</i> ₅ (XLR-11 N4HP- <i>d</i> ₅) (I.S.)	351 > 125	83	29	—