

## Method of Test for 2'-Fucosyllactose in Foods

### 1. Scope

This method is applicable to the determination of 2'-fucosyllactose in infant formula.

### 2. Method

After extraction, 2'-fucosyllactose is determined by high performance ion chromatograph (HPIC).

#### 2.1. Equipment

2.1.1. High performance ion chromatograph.

2.1.1.1. Detector: pulsed electrochemical detector.

2.1.1.1.1. Gold working electrode.

2.1.1.1.2. Ag/AgCl reference electrode.

2.1.1.2. Column: CarboPac PA20, 3 mm × 15 cm, or an equivalent product.

2.1.1.3. Guard column: CarboPac PA20 Guard, 3 mm × 3 cm, or an equivalent product.

2.1.2. Centrifuge: centrifugal force > 9000 ×g.

2.1.3. Ultrasonicator.

2.1.4. Vortex mixer.

2.1.5. Shaker.

#### 2.2. Chemicals

50% Sodium hydroxide, HPIC grade;

Ethanol, reagent grade;

Deionized water, resistivity ≥ 18 MΩ·cm (at 25°C);

2'-Fucosyllactose, reference standard.

#### 2.3. Apparatus

2.3.1. Volumetric flask: 1 mL and 20 mL.

2.3.2. Centrifuge tube: 50 mL, PP.

2.3.3. Membrane filter: 0.22 μm, PVDF.

#### 2.4. 50% ethanol

Dilute 500 mL of ethanol with deionized water to 1000 mL.

#### 2.5. Mobile phase

2.5.1. Solvent A, deionized water.

2.5.2. Solvent B

Dilute 10.5 mL of 50% sodium hydroxide with deionized water to 1000 mL. Filter with a membrane filter, and take the filtrate as the mobile phase B.

## 2.6. Standard solution preparation

Transfer about 5 mg of 2'-fucosyllactose reference standard accurately weighed to a 1-mL volumetric flask. Dissolve and dilute to volume with 50% ethanol as the standard stock solution. When to use, dilute appropriate volume of the standard stock solution with deionized water to 1~10 µg/mL as the standard solutions.

## 2.7. Sample solution preparation

Transfer about 1 g of the homogenized sample accurately weighed into a 10-mL volumetric flask. Add 10 mL of 50% ethanol, vortex-mix, ultrasonicate for 10 min, shake for 10 min, and dilute with 50% ethanol to 20 mL. Centrifuge at 9000 ×g for 30 min. Take 1 mL of the supernatant, and dilute 20 times with deionized water. Filter with a membrane filter, and take the filtrate as the sample solution.

## 2.8. Identification and quantification

Accurately inject 10 µL of the sample solution and the standard solutions into the HPIC separately, and operate according to the following conditions. Identify 2'-fucosyllactose based on the retention time. Calculate the amount of 2'-fucosyllactose in the sample by the following formula:

The amount of 2'-fucosyllactose in the sample (g/100 g)

$$= \frac{C \times V \times F}{M \times 10000}$$

Where,

C: the concentration of 2'-fucosyllactose in the sample solution calculated by the standard curve (µg/mL)

V: the final make-up volume of the sample (20 mL)

M: the weight of the sample (g)

F: dilution factor (20)

HPIC operating conditions<sup>(Note)</sup>:

Detector: pulsed electrochemical detector.

Gold working electrode.

Ag/AgCl reference electrode.

Column: CarboPac PA20, 3 mm × 15 cm.

Guard column: CarboPac PA20 Guard, 3 mm × 3 cm.

Column temperature: 25°C.

Mobile phase: a gradient program of solvent A and solvent B is as follows:

| Time (min) | A (%)   | B (%)     |
|------------|---------|-----------|
| 0 → 15     | 90 → 35 | 10 → 65   |
| 15 → 15.1  | 35 → 0  | 65 → 100  |
| 15.1 → 25  | 0 → 0   | 100 → 100 |
| 25 → 25.1  | 0 → 90  | 100 → 10  |
| 25.1 → 35  | 90 → 90 | 10 → 10   |

Flow rate: 0.5 mL/min.

Injection volume: 10  $\mu$ L.

Note: All the parameters can be adjusted depending on the instruments used if the above conditions are not applicable.

### Remark

1. Limit of quantitation (LOQ) for 2'-fucosyllactose is 0.04 g/100 g.
2. Further validation should be performed when interference compounds appear in samples.

### References

1. Auer, F., Jarvas, G. and Guttman, A. 2021. Recent advances in the analysis of human milk oligosaccharides by liquid phase separation methods. J. Chromatogr. B 1162: 112497.
2. Taiwan Food and Drug Administration. 2015. Method of test for sugars in foods (TFDAO0022.01). Published on December 12, 2015.

### Reference chromatogram




Figure. The HPIC chromatogram of 2'-fucosyllactose standard.