

Method of Test for Synthetic Cannabinoids in Urine (2)

1. Scope

This method is applicable to the determination of 34 synthetic cannabinoids (AB-001 etc. listed as the attached table) in urine.

2. Method

After dilution, analytes are determined by liquid chromatography/tandem mass spectrometry (LC-MS/MS).

2.1. Equipment

2.1.1. Liquid chromatograph/tandem mass spectrometer

2.1.1.1. Ion source: electrospray ionization, ESI.

2.1.1.2. Column: Sunshell® RP-AQUA, 2.6 µm, 2.1 mm i.d. × 10 cm, or an equivalent product.

2.1.2. Vortex mixer

2.2. Chemicals

Methanol, HPLC grade;

Acetonitrile, HPLC grade;

Formic acid, reagent grade;

Artificial urine (UTAK 88121-CDF(L)) or an equivalent product;

Deionized water, resistivity $\geq 18 \text{ M}\Omega \cdot \text{cm}$ (at 25°C);

AB-001 and other 33 synthetic cannabinoids listed in the attached table, reference standards.

AB-CHMINACA-d₄ and 12 other internal standards listed in the attached table.

2.3. Apparatus

2.3.1. Volumetric flask: 1 mL and 10 mL.

2.3.2. Membrane filter: 0.22 µm, PVDF.

2.4. Reagent solution preparation

2.4.1. 50% Methanol

Mix methanol and deionized water at the ratio of 1:1 (v/v).

2.4.2. 50% Acetonitrile in Methanol

Mix acetonitrile and methanol at the ratio of 1:1 (v/v).

2.5. Mobile phase

2.5.1. Solvent A

Dilute 1 mL of formic acid with deionized water to 1000 mL, mix well and filter with a membrane filter.

2.5.2. Solvent B

Dilute 1 mL of formic acid with 50% acetonitrile in methanol to 1000 mL, mix well and filter with a membrane filter.

2.6. Internal standard solution preparation

Transfer 1 mg of the 13 internal standards accurately weighed into each 10-mL volumetric flask, dissolve and dilute with methanol to volume as the internal standard stock solutions. Store at -20°C in the dark. Upon use, mix appropriate volume of the internal standard stock solutions, and dilute with 50% methanol to 1 µg/mL as the internal standard solution.

2.7. Standard solution preparation

Transfer 1 mg of the 34 reference standards accurately weighed into each 10-mL volumetric flask, dissolve and dilute with methanol to volume as the standard stock solution. Store at -20 °C in the dark. Upon use, mix appropriate volume of the standard stock solutions, and dilute with 50% methanol to 500 ng/mL as the standard solution.

2.8. Sample solution preparation

Transfer 100 µL of the homogenized sample and 20 µL of the internal standard solution into a 1-mL volumetric flask, and dilute to volume with 50% methanol. Filter with a membrane filter, and take the filtrate as the sample solution.

2.9. Calibration curve

Use the artificial urine as the blank sample. Transfer 100 µL of the artificial urine, 10-200 µL of the standard solution, and 20 µL of the internal standard solution to a 1-mL volumetric flask and dilute with 50% methanol to volume. Filter with a membrane filter, and take the filtrates as the calibration standard solutions. Operate LC-MS/MS according to the following conditions. Establish the calibration curve of each synthetic cannabinoid by the ratios of the peak area of each synthetic cannabinoid to that of the respective internal standard vs. the added concentrations (5-100 ng/mL).

LC-MS/MS operating conditions^(note):

Column: SunShell® RP-AQUA, 2.6 µm, 2.1 mm i.d. × 10 cm.

Column temperature: 40°C.

Injection volume: 3 µL.

Mobile phase: a gradient program of solvent A and solvent B is as follows.

Time (min)	A (%)	B (%)
0.0 → 0.5	60 → 60	40 → 40
0.5 → 12.0	60 → 10	40 → 90
12.0 → 12.5	10 → 10	90 → 90
12.5 → 12.6	10 → 60	90 → 40
12.6 → 14.5	60 → 60	40 → 40

Flow rate: 0.4 mL/min

Ion spray voltage:

5.5 kV for ESI⁺ Ionization mode;

4.5 kV for ESI⁻ Ionization mode.

Turbo heater temperature: 550°C.

Nebulizer gas, GS1: 50 psi.

Heated gas, GS2: 60 psi.

Curtain gas: 30 psi.

Collision gas: Medium.

Detection mode: multiple reaction monitoring (MRM). Selected ion pair, declustering potential and collision energy are shown in the attached table.

Note: All the parameters can be adjusted depending on the instruments used if the above conditions are not applicable.

2.10. Identification and quantification

Accurately inject 3 µL of the sample solution and the calibration standard solutions into LC-MS/MS separately. Operate according to the conditions in section 2.9. Identify each synthetic cannabinoid based on the retention time and the relative ion intensities^(note). Calculate the amount of each synthetic cannabinoid in the sample by the following formula:

The amount of each synthetic cannabinoid in the sample (ng/mL)

$$= \frac{C \times V}{M}$$

Where,

C: the concentration of each synthetic cannabinoid in the sample solution calculated by the calibration curve (ng/mL)

V: the final make-up volume of the sample (mL)

M: the volume of the sample (mL)

Note: Relative ion intensities are calculated by peak areas of quantitative ions divided by peak areas of qualitative ions ($\leq 100\%$). Maximum permitted tolerances of relative ion intensities are as the following:

Relative ion intensity (%)	Tolerance (%)
> 50	± 20
> 20-50	± 25
> 10-20	± 30
≤ 10	± 50

Remark

1. Limit of quantification (LOQ) for each synthetic cannabinoid is 50 ng/mL.
2. Further validation should be performed when interference compounds appear in the samples.

Reference

Freijo Jr., T. D., Harris, S. E. and Kala, S. V. 2014. A rapid quantitative method for the analysis of synthetic cannabinoids by liquid chromatography–tandem mass spectrometry. J. Anal. Toxicol. 38: 466-478.

Reference chromatograms




Figure. The MRM chromatograms of 34 synthetic cannabinoids and 13 internal standards analyzed by LC-MS/MS.




Figure. The MRM chromatograms of 34 synthetic cannabinoids and 13 internal standards analyzed by LC-MS/MS (Continued).




Figure. The MRM chromatograms of 34 synthetic cannabinoids and 13 internal standards analyzed by LC-MS/MS (Continued).




Figure. The MRM chromatograms of 34 synthetic cannabinoids and 13 internal standards analyzed by LC-MS/MS (Continued).

Table. Parameters of selective ions for 34 synthetic cannabinoids and 13 internal standards.

Analyte	Ionization mode	Ion pair	Declustering potential (V)	Collision energy (eV)	Internal Standard
		Precursor ion (<i>m/z</i>) > product ion (<i>m/z</i>)			
1-Adamantyl-(1-pentylindol-3-yl)methanone (AB-001)	ESI ⁺	350 > 135* 350 > 77	154 154	39 111	JWH-210-d ₉
<i>N</i> [(1 <i>S</i>)-1-(Aminocarbonyl)-2-methylpropyl]-1-(cyclohexylmethyl)-1 <i>H</i> -indazole-3-carboxamide (AB-CHMINACA)	ESI ⁺	357 > 241* 357 > 312	59 59	35 23	AB-CHMINACA-d ₄
<i>N</i> [(1 <i>S</i>)-1-(Aminocarbonyl)-2-methylpropyl]-1-[(4-fluorophenyl)methyl]-1 <i>H</i> -indazole-3-carboxamide (AB-FUBINACA)	ESI ⁺	369 > 253* 369 > 324	83 83	32 21	AB-FUBINACA-d ₄
<i>N</i> [(1 <i>S</i>)-1-(Aminocarbonyl)-2-methylpropyl]-1-pentyl-1 <i>H</i> -indazole-3-carboxamide (AB-PINACA)	ESI ⁺	331 > 215* 331 > 286	40 40	36 20	AB-PINACA-d ₉
<i>N</i> -(1-Amino-3,3-dimethyl-1-oxobutan-2-yl)-1-[(4-fluorophenyl)methyl]indazole-3-carboxamide (ADB-FUBINACA)	ESI ⁺	383 > 338* 383 > 253	72 72	21 36	AB-FUBINACA-d ₄
<i>N</i> [[1-[(4-Fluorophenyl)methyl]-1 <i>H</i> -indazol-3-yl]carbonyl]-L-valine, ethyl ester (AEB-FUBINACA)	ESI ⁺	398 > 109* 398 > 253	124 124	59 33	AB-FUBINACA-d ₄
<i>N</i> -(1-Adamantyl)-1-pentylindazole-3-carboxamide (AKB48)	ESI ⁺	366 > 135* 366 > 93	76 76	26 66	AKB48-d ₉
[1-(5-Fluoropentyl)indol-3-yl]-naphthalen-1-ylmethanone (AM-2201)	ESI ⁺	360 > 155* 360 > 127	117 117	37 70	AM-2201-d ₅
<i>N</i> -(1-Adamantyl)-1-(5-bromopentyl)indazole-3-carboxamide (5-Bromo-AKB48)	ESI ⁺	444 > 135* 444 > 79	36 36	30 96	AKB48-d ₉
<i>N</i> -(1-Adamantyl)-1-(5-chloropentyl)indazole-3-carboxamide (5-Chloro-AKB48)	ESI ⁺	400 > 135* 400 > 93	39 39	26 71	AKB48-d ₉
(1-(5-Chloropentyl)-1 <i>H</i> -indol-3-yl)(2,2,3,3-tetramethylcyclopropyl)methanone (5-Chloro-UR-144)	ESI ⁺	346 > 248* 346 > 144	100 100	37 51	UR-144-d ₅

*Quantitative ion

Table. Parameters of selective ions for 34 synthetic cannabinoids and 13 internal standards (Continued).

Analyte	Ionization mode	Ion pair	Declustering potential (V)	Collision energy (eV)	Internal Standard
		Precursor ion (<i>m/z</i>) > product ion (<i>m/z</i>)			
<i>rel</i> -5-(1,1-Dimethylheptyl)-2-[(1 <i>R</i> ,3 <i>S</i>)-3-hydroxycyclohexyl]-phenol (CP47,497)	ESI ⁻	317 > 159* 317 > 245	140 140	71 45	CP47,497- d ₁₁
Methyl (2 <i>S</i>)-2-[[1-(4-fluorobutyl)indazole-3-carbonyl]amino]-3,3-dimethylbutanoate (4F-MDMB-BINACA)	ESI ⁺	364 > 219* 364 > 304	60 60	33 22	AB- CHMINACA-d ₄
Methyl 2-[[1-(5-fluoropentyl)indazole-3-carbonyl]amino]-3,3-dimethylbutanoate (5F-ADB)	ESI ⁺	378 > 233* 378 > 318	97 97	33 23	AM-2201-d ₅
<i>N</i> -(1-Adamantyl)-1-(5-fluoropentyl)-1 <i>H</i> -Indazole-3-Carboxamide (5F-AKB48)	ESI ⁺	384 > 135* 384 > 93	72 72	26 68	AKB48-d ₉
Methyl (2 <i>S</i>)-2-[[1-(5-fluoropentyl)indazole-3-carbonyl]amino]-3-methylbutanoate (5F-AMP)	ESI ⁺	364 > 304* 364 > 233	108 108	22 31	AB- FUBINACA-d ₄
Methyl (2 <i>S</i>)-2-[[1-(5-fluoropentyl)indole-3-carbonyl]amino]-3,3-dimethylbutanoate (5F-MDMB-PICA)	ESI ⁺	377 > 232* 377 > 144	89 89	32 48	XLR-11-d ₅
Quinolin-8-yl 1-[(4-fluorophenyl)methyl]indole-3-carboxylate (FUB-PB-22)	ESI ⁺	397 > 252* 397 > 109	78 78	24 49	AM-2201-d ₅
Naphthalen-1-yl-(1-pentylindol-3-yl)methanone (JWH-018)	ESI ⁺	342 > 155* 342 > 127	110 110	32 72	JWH-018-d ₉
(1-Hexylindol-3-yl)-naphthalen-1-ylmethanone (JWH-019)	ESI ⁺	356 > 127* 356 > 155	78 78	33 64	JWH-210-d ₉
Naphthalen-1-yl-(1-pent-4-enylindol-3-yl)methanone (JWH-022)	ESI ⁺	340 > 155* 340 > 127	142 142	32 60	JWH-250-d ₅
(1-Butylindol-3-yl)-naphthalen-1-ylmethanone (JWH-073)	ESI ⁺	328 > 155* 328 > 127	151 151	32 65	JWH-073-d ₇
(4-Methoxynaphthalen-1-yl)-(1-pentylindol-3-yl)methanone (JWH-081)	ESI ⁺	372 > 185* 372 > 114	102 102	35 103	JWH-081-d ₅
(4-Methylnaphthalen-1-yl)-(1-pentylindol-3-yl)methanone (JWH-122)	ESI ⁺	356 > 169* 356 > 115	57 57	34 86	JWH-210-d ₉

*Quantitative ion

Table. Parameters of selective ions for 34 synthetic cannabinoids and 13 internal standards (Continued).

Analyte	Ionization mode	Ion pair	Declustering potential (V)	Collision energy (eV)	Internal Standard
		Precursor ion (<i>m/z</i>) > product ion (<i>m/z</i>)			
2-(2-Chlorophenyl)-1-(1-pentylindol-3-yl)ethenone (JWH-203)	ESI ⁺	340 > 125*	136	36	JWH-018-d ₉
		340 > 214	136	35	
(4-Ethynaphthalen-1-yl)-(1-pentylindol-3-yl)methanone (JWH-210)	ESI ⁺	370 > 183*	138	33	JWH-210-d ₉
		370 > 214	138	34	
2-(2-Methoxyphenyl)-1-(1-pentylindol-3-yl)ethanone (JWH-250)	ESI ⁺	336 > 121*	129	27	JWH-250-d ₅
		336 > 91	129	65	
Methyl (2 <i>S</i>)-2-[[1-(cyclohexylmethyl)indole-3-carbonyl]amino]-3,3-dimethylbutanoate (MDMB-CHMICA)	ESI ⁺	385 > 240*	80	25	XLR-11-d ₅
		385 > 144	80	55	
Methyl (2 <i>S</i>)-2-[[1-(cyclohexylmethyl)indole-3-carbonyl]amino]-3-methylbutanoate (MMB-CHMICA)	ESI ⁺	371 > 240*	32	23	AB-CHMINACA-d ₄
		371 > 144	32	50	
Methyl 2-[[1-(4-fluorophenyl)methyl]indazole-3-carbonyl]amino]-3-methylbutanoate (MMB-FUBINACA)	ESI ⁺	384 > 253*	106	32	AB-FUBINACA-d ₄
		384 > 109	106	66	
Naphthalen-1-yl 1-(5-fluoropentyl)indole-3-carboxylate (NM-2201)	ESI ⁺	376 > 232*	55	27	JWH-210-d ₉
		376 > 144	55	52	
[1-(5-Fluoropentyl)indazol-3-yl]-naphthalen-1-ylmethanone (THJ-2201)	ESI ⁺	361 > 233*	126	24	JWH-210-d ₉
		361 > 145	126	48	
(1-Pentylindol-3-yl)-(2,2,3,3-tetramethylcyclopropyl)methanone (UR-144)	ESI ⁺	312 > 125*	56	30	UR-144-d ₅
		312 > 214	56	31	
[1-(5-Fluoropentyl)indol-3-yl]-(2,2,3,3-tetramethylcyclopropyl)methanone (XLR-11)	ESI ⁺	330 > 125*	126	30	XLR-11-d ₅
		330 > 232	126	33	
<i>N</i> -[(2 <i>S</i>)-1-Amino-3-methyl-1-oxobutan-2-yl]-1-(cyclohexylmethyl)indazole-3-carboxamide-d ₄ (AB-CHMINACA-d ₄) (I.S.)	ESI ⁺	361 > 245	114	35	-

*Quantitative ion

Table. Parameters of selective ions for 34 synthetic cannabinoids and 13 internal standards (Continued).

Analyte	Ionization mode	Ion pair		Declustering potential (V)	Collision energy (eV)	Internal Standard
		Precursor ion (<i>m/z</i>)	> product ion (<i>m/z</i>)			
<i>N</i> -[(2 <i>S</i>)-1-Amino-3-methyl-1-oxobutan-2-yl]-1-[(4-fluorophenyl)methyl]indazole-3-carboxamide- <i>d</i> ₄ (AB-FUBINACA- <i>d</i> ₄) (I.S.)	ESI ⁺	373	> 328	68	21	-
<i>N</i> -[(1 <i>S</i>)-1-(Aminocarbonyl)-2-methylpropyl]-1-pentyl-1 <i>H</i> -indazole-3-carboxamide- <i>d</i> ₉ (AB-PINACA- <i>d</i> ₉) (I.S.)	ESI ⁺	340	> 224	40	36	-
<i>N</i> -(1-Adamantyl)-1-pentylindazole-3-carboxamide- <i>d</i> ₉ (AKB48- <i>d</i> ₉) (I.S.)	ESI ⁺	375	> 135	69	15	-
[1-(5-Fluoropentyl)indol-3-yl]-naphthalen-1-ylmethanone- <i>d</i> ₅ (AM-2201- <i>d</i> ₅) (I.S.)	ESI ⁺	365	> 127	31	73	-
<i>rel</i> -5-(1,1-Dimethylheptyl)-2-[(1 <i>R</i> ,3 <i>S</i>)-3-hydroxycyclohexyl]-phenol- <i>d</i> ₁₁ (CP47,497- <i>d</i> ₁₁) (I.S.)	ESI ⁻	328	> 256	72	47	-
Naphthalen-1-yl-(1-pentylindol-3-yl)methanone- <i>d</i> ₉ (JWH-018- <i>d</i> ₉) (I.S.)	ESI ⁺	351	> 155	36	33	-
(1-Butylindol-3-yl)-naphthalen-1-ylmethanone- <i>d</i> ₇ (JWH-073- <i>d</i> ₇) (I.S.)	ESI ⁺	335	> 155	76	33	-
(4-Methoxynaphthalen-1-yl)-(1-pentylindol-3-yl)methanone- <i>d</i> ₉ (JWH-081- <i>d</i> ₉) (I.S.)	ESI ⁺	381	> 185	51	35	-
(4-Ethyl-naphthalen-1-yl)-(1-pentylindol-3-yl)methanone- <i>d</i> ₉ (JWH-210- <i>d</i> ₉) (I.S.)	ESI ⁺	379	> 183	36	33	-
2-(2-Methoxyphenyl)-1-(1-pentylindol-3-yl)ethanone- <i>d</i> ₅ (JWH-250- <i>d</i> ₅) (I.S.)	ESI ⁺	341	> 121	46	27	-
(1-Pentylindol-3-yl)-(2,2,3,3-tetramethylcyclopropyl)methanone- <i>d</i> ₅ (UR-144- <i>d</i> ₅) (I.S.)	ESI ⁺	317	> 125	161	31	-
[1-(5-fluoropentyl)indol-3-yl]-(2,2,3,3-tetramethylcyclopropyl)methanone- <i>d</i> ₅ (XLR-11- <i>d</i> ₅) (I.S.)	ESI ⁺	335	> 125	61	31	-

*Quantitative ion