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
Microorganism	Preparation of Test Strain	Growth Promotion	
		Total Aerobic Microbial Count	Total Yeasts and Molds Count
<i>Staphylococcus aureus</i> such as ATCC 6538, NCIMB 9518, CIP 4.83, or NBRC 13276	Soybean-Casein Digest Agar or Soybean-Casein Digest Broth 30°-35° 18-24 hours	Soybean-Casein Digest Agar and Soybean-Casein Digest Broth ≤ 100 cfu 30°-35° ≤ 3 days	
<i>Pseudomonas aeruginosa</i> such as ATCC 9027, NCIMB 8626, CIP 82.118, or NBRC 13275	Soybean-Casein Digest Agar or Soybean-Casein Digest Broth 30°-35° 18-24 hours	Soybean-Casein Digest Agar and Soybean-Casein Digest Broth ≤ 100 cfu 30°-35° ≤ 3 days	
<i>Bacillus subtilis</i> such as ATCC 6633, NCIMB 8054, CIP 52.62, or NBRC 3134	Soybean-Casein Digest Agar or Soybean-Casein Digest Broth 30°-35° 18-24 hours	Soybean-Casein Digest Agar and Soybean-Casein Digest Broth ≤ 100 cfu 30°-35° ≤ 3 days	
<i>Candida albicans</i> such as ATCC 10231, NCPF 3179, IP 48.72, or NBRC 1594	Sabouraud Dextrose Agar or Sabouraud Dextrose Broth 20°-25° 2-3 days	Soybean-Casein Digest Agar ≤ 100 cfu 30°-35° ≤ 5 days	Sabouraud Dextrose Agar ≤ 100 cfu 20°-25° ≤ 5 days
<i>Aspergillus niger</i> such as ATCC 16404, IMI 149007, IP 1431.83, or NBRC 9455	Sabouraud Dextrose Agar or Potato-Dextrose Agar 20°-25° 5-7 days, or until good sporulation is achieved	Soybean-Casein Digest Agar ≤ 100 cfu 30°-35° ≤ 5 days	Sabouraud Dextrose Agar ≤ 100 cfu 20°-25° ≤ 5 days

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1.05458.0500 Cryptic soy agar Casein-peptone soymeal-peptone agar
for microbiology (According harm. EP/USP/JP and
ISO)

Batch VM166358

	Spec. Values	Batch Values
Appearance		
clearness	clear	clear
colour	yellowish brown	yellowish-brown
pH-value (25 °C)	7.1 - 7.5	7.3
Solidification behaviour (2 hrs., 45 °C)	liquid	liquid
Stability test (Colour and hemolysis)	non-hemolytic	non-hemolytic

*Typical composition (g/litre): Peptone from casein 15.0; Peptone from soymeal 5.0; Sodium chloride 5.0;
Agar-agar 15.0.*

Growth promotion test in accordance with the harmonised method of EP, USP and JP.

	Spec. Values	Batch Values
Inoculum		
Staphylococcus aureus ATCC 6538	10 - 100	93
Bacillus subtilis ATCC 6633	10 - 100	16
Escherichia coli ATCC 8739	10 - 100	16
Streptococcus pyogenes ATCC 21059	10 - 100	98
Pseudomonas aeruginosa ATCC 9027	10 - 100	28
Candida albicans ATCC 10231	10 - 100	94
Aspergillus brasiliensis ATCC 16404 (formerly A. niger)	10 - 100	64

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
Product Name	: BOTTLE VOGEL JOHNSON AGAR 500G		
Catalog Number	: 256220	Manufacture Date	: 2008/12/05
Batch Number	: 9127032		
Expiration Date	: 2013/10/31		

01. Dehydrated Medium Appearance: Pink, free-flowing, homogeneous
02. Solubility: 6% solution, soluble in distilled or deionized water on boiling
03. Solution Appearance: Red, very slightly to slightly opalescent, may have a slight white precipitate
04. Cultural Response: Medium was prepared per label instructions using 20 ml/l of Tellurite Solution 1%. Plates were inoculated with the test organisms and incubated at $35 \pm 2^{\circ}\text{C}$ for 18-48 hours.


TEST ORGANISMS	RECOVERY	TELLURITE REDUCTION	MANNITOL FERMENTATION
Escherichia coli ATCC® 25922	markedly to completely inhibited	-	-
Proteus mirabilis ATCC® 25933	partial to complete inhibition	+	-
Staphylococcus aureus ATCC® 25923	good	+	+
Staphylococcus epidermidis ATCC® 12228	fair to good	±	-

KEY: + = positive, black colony for tellurite reduction
yellow medium for mannitol fermentation

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



僅准許被授權人員進入
ADMITTANCE TO AUTHORIZED PERSONNEL ONLY II

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

- 凍乾管開封活化和分裝儲存記錄
 - Name of microorganism
 - BCRC strain no./lot no.
 - Date of resuscitation
 - Media used/lot no.
 - Incubation temperature/time
 - Date of preparation of cryovials
 - Total numbers of cryovials prepared

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

- 標準菌種儲存 (Storage of stock cultures)
 - The “Seed-Lot” technique is recommended for storage of stock cultures.
 - Aliquots of stock culture (the first transfer or passage) are suspended in a cryoprotective medium, transfer to vials, and frozen at -30°C or below, until use.
 - If stored at -70°C , or in lyophilized form, strains may be kept indefinitely.
 - The freezer should be calibrated and mapped.

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
- Those frozen stocks can then be used to inoculated monthly or weekly working cultures.
- Once opened, do not refreeze unused cell suspensions after culturing a working suspension.
- The number of transfers of working cultures should be tracked and not more than five passages.

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
- One passage is defined as the transfer of organisms from a viable culture to a fresh medium with growth of the microorganisms.
- Any form of subculturing is considered to be a transfer/passage.

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
- 冷凍保存管使用記錄
 - Date of use
 - Cryovial no.
 - Use by
 - Subculture medium/Lot no.

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- Confirmation of the purity of the culture and the identity of the culture should be performed prior to its use in quality control testing.
- 菌體懸浮液於TSA Agar上，做劃線分離培養，用無菌L型玻棒均勻塗抹，以檢驗菌種之純度

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