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


Approval standards	ISO 21149, Cosmetics — Microbiology — Enumeration and detection of aerobic mesophilic bacteria	ASEAN microbial limit test for cosmetic products ACM THA 06	China: Standard methods of microbiological examination for cosmetics , General rules 化妆品微生物标准检验方法 细菌总数测定 GB 7918.1-87	U.S. Pharmacopeia - USP 37- <61> Microbial examination of nonsterile products: microbial enumeration tests	COLIPA, Guidelines on Microbial Quality Management, published by the European Cosmetic, Toiletry and Perfumery Association, 1997	E P, Microbiological Examination of non-sterile products, 4th edition, published by the European Pharmacopoeia, 2002	FDA, Bacteriological Analytical Manual, published by the U.S. Food and Drug Administration, 2017
Recommended media list	Bacteria: Soybean-casein digest agar (SCDA) ; Eugon LT 100 agar	Bacteria: Modified letheen agar (MLA) Mould and yeast: Malt extract agar (MEA) ; Potato dextrose agar (PDA)	Bacteria: SCDLP medium (Soybean-casein digest agar ; SCDA with polysorbate 80, lecithin)	Bacteria: Soybean-casein digest agar (SCDA) ; Trypticase soy agar (TSA) Sabouraud dextrose agar	Bacteria: Soybean-casein digest agar (SCDA)/buffered sodium chloride-peptone solution with neutralizing agents. e.g. polysorbate 20 or 80, lecithin, thiosulphate	Bacteria: Fluid thioglycollate medium ; Soybean-casein digest agar (SCDA)	Bacteria: Modified letheen agar (MLA) (M78) Mould and yeast: Malt extract agar (MEA) (M93) Commercial bacterial identification kit (API or equivalent)
Validation test	Validation of the efficiency of neutralizer	Validation of the efficiency of neutralizer	Not mentioned	Promotion tests for media are required	Not mentioned	Promotion and validation tests for media are required	Not mentioned

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Negative control:
dilute broth/ no
contamination

Positive control:
Growth promotion

**Treatment
methods:**
Neutralization
/ filtration/ dilution

TABLE 3. Culture Conditions for Preparation of Standardized Cell Suspensions and for Microbial Recovery Studies in Method Validation Experiment (United States Pharmacopeia Chapter <61>).¹¹


ORGANISM	CULTURE MEDIUM	TEMP (°C)	TIME (CELL SUSPENSION)	TIME (RECOVERY)
<i>S. aureus</i>	Soybean-Casein Digest (broth, agar)	30 to 35	18 to 24 hours	≤3 days
<i>P. aeruginosa</i>	Soybean-Casein Digest (broth, agar)	30 to 35	18 to 24 hours	≤3 days
<i>B. subtilis</i>	Soybean-Casein Digest (broth, agar)	30 to 35	18 to 24 hours	≤3 days
<i>C. albicans</i>	Sabouraud Dextrose (agar, broth)	20 to 25	2 to 3 days	≤5 days ^a
<i>A. brasiliensis</i>	Sabouraud Dextrose Agar or Potato-Agar	20 to 25	5 to 7 days	≤5 days ^a

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