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# Safety Gate: the EU rapid alert system for dangerous non-food products

Notifying country	France
Product category	Cosmetics
Product	Shampoo
Name	Shampooing
Brand	Marcapar barbier
Batch number	BSH1/Y091 BSH2/X051
Type / number of model	250ML - 007670 500ML - 007671
Barcode	3760172621836 3760172621843
Product description	Liquid shampoo, verbena/bergamot fragrance intended for men. This product is supplied in 250 and 500 ml bottle. Product sold online.
Packaging description	Plastic bottle.
Country of origin	France
Risk type	Microbiological
Risk description	The product is microbiologically contaminated with the bacteria <i>Pseudomonas aeruginosa</i> (measured value up to: 10 000 CFU/g). The microorganism may potentially cause illness in some consumers, particularly those with a weakened immune system; and if used on damaged skin, or if it comes into contact with the eyes, the product may cause infection or irritation.

1 Media files

[View all](#)

<https://ec.europa.eu/safety-gate-alerts/screen/webReport/alertDetail/10009816?lang=en>

# Microbiological contamination of cosmetic products – observations from Europe, 2005–2018

Species of bacteria found in microbiologically contaminated cosmetic products registered in the Rapex database in 2005 –2018

Microorganisms identified	Number or reported products	Per cent of all of the reported products
<b>Gram-negative bacteria</b>	<b>62</b>	<b>59.62</b>
<b>Achromobacter</b>	<b>1</b>	<b>0.96</b>
<i>Achromobacter xylosoxidans</i>	1	0.96
<b>Burkholderia</b>	<b>4</b>	<b>3.85</b>
<i>Burkholderia cepacia</i>	3	2.88
<i>Burkholderia species</i>	1	0.96
<b>Citrobacter</b>	<b>1</b>	<b>0.96</b>
<i>Citrobacter freundii</i>	1	0.96
<b>Enterobacter</b>	<b>12</b>	<b>11.54</b>
<i>Enterobacter cloacae</i>	4	3.85
<i>Enterobacter gergoviae</i>	8	7.69
<b>Klebsiella</b>	<b>3</b>	<b>2.88</b>
<i>Klebsiella oxytoca</i>	1	0.96
<i>Klebsiella pneumoniae</i>	2	1.92
<b>Pantoea</b>	<b>2</b>	<b>1.92</b>
<i>Pantoea agglomerans</i>	1	0.96
<i>Pantoea species</i>	1	0.96
<b>Pseudomonas</b>	<b>37</b>	<b>35.58</b>
<i>Pseudomonas aeruginosa</i>	30	28.84
<i>Pseudomonas putida</i>	1	0.96
<i>Pseudomonas species</i>	6	5.77

<b>Rhizobium</b>	<b>1</b>	<b>0.96</b>
<i>Rhizobium radiobacter</i>	1	0.96
<b>Serratia</b>	<b>1</b>	<b>0.96</b>
<i>Serratia marcescens</i>	1	0.96
<b>Gram-positive bacteria</b>	<b>7</b>	<b>6.73</b>
<b>Bacillus</b>	<b>1</b>	<b>0.96</b>
<i>Bacillus firmus</i>	1	0.96
<b>Enterococcus</b>	<b>2</b>	<b>1.92</b>
<i>Enterococcus faecium</i>	1	0.96
<i>Enterococcus species</i>	1	0.96
<b>Staphylococcus</b>	<b>4</b>	<b>3.85</b>
<i>Staphylococcus aureus</i>	4	3.85
<b>Fungi</b>	<b>2</b>	<b>1.92</b>
<b>Candida</b>	<b>2</b>	<b>1.92</b>
<i>Candida albicans</i>	2	1.92
<b>Unspecified microbiological risk</b>	<b>45</b>	<b>43.37</b>
Mould	2	1.92
Unspecified bacteria	36	34.62
Unspecified bacteria and mould	1	0.96
Unspecified bacteria, yeast, and mould	4	3.85
Unspecified microbiological contamination	1	0.96
Yeast	1	0.96

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## Annex D (informative)

### Neutralizers of antimicrobial activity of preservatives and rinsing liquids

Preservative	Chemical compounds able to neutralize preservative's antimicrobial activity	Examples of suitable neutralizers and of rinsing liquids <sup>[6] [7]</sup> (for membrane filtration methods)
Phenolic compounds: parabens, phenoxyethanol, phenylethanol, etc. Anilides	Lecithin, Polysorbate 80, Ethylene oxide condensate of fatty alcohol Non-ionic surfactants	Polysorbate 80, 30 g/l + lecithin, 3 g/l. Ethylene oxide condensate of fatty alcohol, 7 g/l + lecithin, 20 g/l + polysorbate 80, 4 g/l. D/E neutralizing broth <sup>a</sup> . <i>Rinsing liquid: distilled water; tryptone, 1 g/l + NaCl 9 g/l; polysorbate 80, 5 g/l.</i>
Quaternary ammonium compounds Cationic surfactants	Lecithin, saponin, polysorbate 80, sodium dodecyl sulphate Ethylene oxide condensate of fatty alcohol	Polysorbate 80, 30 g/l + sodium dodecyl sulphate, 4 g/l + lecithin, 3 g/l. Polysorbate 80, 30 g/l + saponin, 30 g/l + lecithin, 3 g/l. D/E neutralizing broth <sup>a</sup> . <i>Rinsing liquid : distilled water; tryptone, 1 g/l + NaCl, 9 g/l; polysorbate 80, 5 g/l.</i>

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## H. Microbiological enumeration and growth

### 1. Aerobic plate count (APC). FDA BAM 23

Use spread plate technique to facilitate recognition of different colony types. Decimally dilute the cosmetic preparation (see G, above) in MLB to obtain a complete dilution series from  $10^{-1}$  to  $10^{-3}$ . Use a new sterile pipet to transfer 1.0 ml of the current dilution into 9 ml of fresh MLB to make the next decimal dilution.

### 9.3.1 Dilutions for counting methods ISO 21149 : 2017

Usually, the initial suspension is the first counted dilution. If needed, additional serial dilutions (e.g. 1:10 dilution) may be performed from the initial suspension using the same diluent (according to the expected level of contamination of the product).

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

## 10 Counting of colonies (plate counts and membrane filtration methods)

After incubation, count the colonies:

- in Petri dishes containing 30 colonies to 300 colonies; if less than 30 colonies are counted, see 12.2.3;
- on membranes containing 15 colonies to 150 colonies; if less than 15 colonies are counted, see 12.2.3.

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
In addition, plate 1 ml of  $10^{-1}$  dilution onto one set of either two (0.5 ml each), or three (0.3, 0.3, 0.4 ml each) MLA plates. Regardless of the number of plates used in plating, the total volume should add up to 1 ml. Do the same for a duplicate set. The dilution factor is 10.

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