

Our Complete Product Range from 18G to 29G

Grooved transparent hub

- provides rapid detection of C.S.F. flashback confirming accurate placement

Colour-coded stylet according to ISO and DIN

- ensures rapid identification of the needle gauge
- enables perfect matching: each stylet fits exactly into the respective needle's lumen

Optimum needle wall thickness (ID/OD)

- assures highest possible strength and stability of the needle along its entire length
- allows the maximum liquor flow rate for any given needle gauge
- thin-wall spinal needles available in 26G, 27G and 29G

The complete range

- provides a wide selection of spinal needles for all indications and for every field of application
- takes into account modern working practices

